

# TCM4/40 operator's manual



TCM4/40  
operator's  
manual



# TCM4/40 monitoring systems

## Operator's manual

From software version 3.0



**RADIOMETER** 

Table of contents

1. Introduction

2. What is what

3. Menu structure and setup programs

4. Installation and maintenance

5. Calibration

6. In vivo monitoring

7. Troubleshooting

8. Specifications and ordering information

9. Functional description

Index

Date of issue

### **System performance**

The procedures described in this manual must be observed in order to ensure proper system performance, and to avoid hazards.

Radiometer cannot provide or verify system performance characteristics if the system is not installed, used and maintained in accordance with Radiometer procedures or if accessories not meeting the specifications provided by Radiometer are used.

Radiometer warrants that the data media on which the software included in the system is furnished is free from defects in material and workmanship under normal use for three (3) months from the date of delivery as evidenced by a copy of invoice or receipt.

### **Third-party software and trademarks**

The TCM4/40 monitoring systems comprise the Microsoft® Windows® CE 5.0 program.

By using the system, you accept the terms of the Software License Agreement(s) of the provider(s) of the above software as shown in the End User License Agreement(s) included in this manual. If you cannot accept the terms of the Software License Agreement(s), you should not use the system, but immediately contact your provider for a return of the system and a refund of the purchase price.

Microsoft® and Windows® are trademarks of Microsoft Corporation.

### **Warranties and disclaimer**

Radiometer makes no warranties, express or implied, other than expressly stated.

Any warranties expressly stated in this document are conditional upon the system being installed, used and maintained in accordance with Radiometer procedures, including that only accessories meeting the specifications provided by Radiometer are used.

Radiometer disclaims any liability for system performance if the system is not installed, used and maintained in accordance with Radiometer procedures or if accessories not meeting the specifications provided by Radiometer are used.

Further, Radiometer disclaims any liability for loss of data and direct, consequential or other damages, including loss of profit or loss of business, whether such claim for damages is based upon contract, negligence or tort (including strict liability), and even if Radiometer has knowledge of the possibility of the potential damage or loss.

### **Confidentiality**

The contents of this document shall not be reproduced or communicated to any third party without the prior written consent of Radiometer.

### **Changes**

This document is subject to change without notice and you are urged to contact Radiometer to verify whether the document has been changed.

While every effort is made to ensure the correctness of the information provided in this document as changed from time to time, Radiometer disclaims any liability for errors and omissions.

## End user license agreement for Microsoft Software

- You have acquired a device ("TCM4/40 monitor") that includes software licensed by Radiometer Medical ApS from Microsoft Licensing Inc. or its affiliates ("MS"). Those installed software products of MS origin, as well as associated media, printed materials and "online" or electronic documentation ("SOFTWARE") are protected by international intellectual property laws and treaties. The SOFTWARE is licensed, not sold. All rights reserved.
- IF YOU DO NOT AGREE TO THIS END USER LICENSE AGREEMENT ("EULA"), DO NOT USE THE DEVICE OR COPY THE SOFTWARE. INSTEAD, PROMPTLY CONTACT Radiometer Medical ApS FOR INSTRUCTIONS ON RETURN OF THE UNUSED DEVICE(S) FOR A REFUND. **ANY USE OF THE SOFTWARE, INCLUDING BUT NOT LIMITED TO USE ON THE DEVICE, WILL CONSTITUTE YOUR AGREEMENT TO THIS EULA (OR RATIFICATION OF ANY PREVIOUS CONSENT).**
- **GRANT OF SOFTWARE LICENSE.** This EULA grants you the following license:
  - You may use the SOFTWARE only on the DEVICE.
  - **NOT FAULT TOLERANT.** THE SOFTWARE IS NOT FAULT TOLERANT. Radiometer Medical ApS HAS INDEPENDENTLY DETERMINED HOW TO USE THE SOFTWARE IN THE DEVICE, AND MS HAS RELIED UPON Radiometer Medical ApS TO CONDUCT SUFFICIENT TESTING TO DETERMINE THAT THE SOFTWARE IS SUITABLE FOR SUCH USE.
  - **NO WARRANTIES FOR THE SOFTWARE. THE SOFTWARE is provided "AS IS" and with all faults. THE ENTIRE RISK AS TO SATISFACTORY QUALITY, PERFORMANCE, ACCURACY, AND EFFORT (INCLUDING LACK OF NEGLIGENCE) IS WITH YOU. ALSO, THERE IS NO WARRANTY AGAINST INTERFERENCE WITH YOUR ENJOYMENT OF THE SOFTWARE OR AGAINST INFRINGEMENT.** IF YOU HAVE RECEIVED ANY WARRANTIES REGARDING THE DEVICE OR THE SOFTWARE, THOSE WARRANTIES DO NOT ORIGINATE FROM, AND ARE NOT BINDING ON, MS.
  - **Note on Java Support.** The SOFTWARE may contain support for programs written in Java. Java technology is not fault tolerant and is not designed, manufactured, or intended for use or resale as online control equipment in hazardous environments requiring fail-safe performance, such as in the operation of nuclear facilities, aircraft navigation or communication systems, air traffic control, direct life support machines, or weapons systems, in which the failure of Java technology could lead directly to death, personal injury, or severe physical or environmental damage. Sun Microsystems, Inc. has contractually obligated MS to make this disclaimer.
  - **No Liability for Certain Damages. EXCEPT AS PROHIBITED BY LAW, MS SHALL HAVE NO LIABILITY FOR ANY INDIRECT, SPECIAL, CONSEQUENTIAL OR INCIDENTAL DAMAGES ARISING FROM OR IN CONNECTION WITH THE USE OR PERFORMANCE OF THE SOFTWARE. THIS LIMITATION SHALL APPLY EVEN IF ANY REMEDY FAILS OF ITS ESSENTIAL PURPOSE. IN NO EVENT SHALL MS BE LIABLE FOR ANY AMOUNT IN EXCESS OF U.S. TWO HUNDRED FIFTY DOLLARS (U.S.\$250.00).**
  - **Limitations on Reverse Engineering, Decompilation, and Disassembly.** You may not reverse engineer, decompile, or disassemble the SOFTWARE, except and only to the extent that such activity is expressly permitted by applicable law notwithstanding this limitation.
  - **SOFTWARE TRANSFER ALLOWED BUT WITH RESTRICTIONS.** You may permanently transfer rights under this EULA only as part of a permanent sale or transfer of the Device, and only if the recipient agrees to this EULA. If the SOFTWARE is an upgrade, any transfer must also include all prior versions of the SOFTWARE.
  - **EXPORT RESTRICTIONS.** You acknowledge that SOFTWARE is of US-origin. You are to comply with all applicable international and national laws that apply to the SOFTWARE, including the U.S. Export Administration Regulations, as well as end-user, end-use and country destination restrictions issued by U.S. and other governments. For additional information on exporting the SOFTWARE, see <http://www.microsoft.com/exporting/>.





# Contents

<b>1. Introduction.....</b>	<b>1-1</b>
Names, intended use and limitations.....	1-2
Symbols used on the monitor.....	1-4
Symbols used in the manual.....	1-6
<b>2. What is what.....</b>	<b>2-1</b>
TCM4/40 monitoring systems.....	2-2
Monitor – top and front.....	2-3
Monitor – rear.....	2-4
The screen: general elements.....	2-5
Online tutorials.....	2-8
Touch key glossary.....	2-9
<b>3. Menu structure and setup programs.....</b>	<b>3-1</b>
Menu structure.....	3-2
List of setup programs.....	3-3
<b>View setup.....</b>	<b>3-4</b>
Normal view.....	3-5
Trend table view.....	3-7
Trend curve view.....	3-8
<b>Parameter setup.....</b>	<b>3-9</b>
$p\text{CO}_2$ .....	3-10
$p\text{O}_2$ .....	3-11
SpO <sub>2</sub> /Pulse.....	3-12
Blood gas.....	3-13
<b>Calibration setup.....</b>	<b>3-14</b>
SmartCal.....	3-15
Calibration status.....	3-16
<b>Printer setup.....</b>	<b>3-17</b>
<b>Technical setup.....</b>	<b>3-19</b>
Technical settings.....	3-20
Date/time.....	3-22
Default values.....	3-23
<b>4. Installation and maintenance.....</b>	<b>4-1</b>
Operating requirements.....	4-2
Installing the monitoring system.....	4-3
Shutting down the monitor.....	4-5

Cleaning the monitor .....	4-6
Maintenance of the monitor .....	4-7
<b>Preparation and maintenance of the tc sensors .....</b>	<b>4-9</b>
General information about membraning .....	4-10
Cleaning the sensor head of the E5480 sensor .....	4-11
Membraning the E5480 sensor .....	4-12
Membraning the E5280 and E5260 sensors .....	4-13
Cleaning, disinfection and storage of sensors .....	4-15
<b>Maintenance of the SpO<sub>2</sub> sensors .....</b>	<b>4-16</b>
<b>5. Calibration .....</b>	<b>5-1</b>
General information .....	5-2
Calibration .....	5-3
Checking zero current and sensitivity .....	5-4
<b>6. Patient monitoring .....</b>	<b>6-1</b>
General information .....	6-2
Application and removal of tc sensors .....	6-6
Application of SpO <sub>2</sub> sensors .....	6-8
Patient monitoring (In vivo monitoring) .....	6-11
Patient DMS .....	6-13
Results in Normal view .....	6-16
Results in Trend table view .....	6-18
Results in Trend curve view .....	6-19
Analog output .....	6-21
Continuous data output: standard .....	6-22
Continuous data output: VueLink .....	6-24
Continuous data output: MonLink .....	6-29
Data export: serial .....	6-30
Data export: USB .....	6-33
Alarms .....	6-34
How to print .....	6-38
Blood gas comparison .....	6-44
In vivo calibration .....	6-45
<b>7. Troubleshooting .....</b>	<b>7-1</b>
The TCM4/40 systems .....	7-2
<b>8. Specifications and ordering information .....</b>	<b>8-1</b>
Specifications .....	8-2
Accessories .....	8-14

<b>9. Functional description .....</b>	<b>9-1</b>
<b>tcpCO<sub>2</sub>/tcpO<sub>2</sub> measurement.....</b>	<b>9-2</b>
Construction of sensors .....	9-3
pCO <sub>2</sub> measuring principle .....	9-5
pO <sub>2</sub> measuring principle.....	9-6
Calibration of sensor .....	9-7
Solutions and calibration gases .....	9-10
Traceability certificates .....	9-11
<b>Pulse oximetry measurement .....</b>	<b>9-15</b>
Measuring principle .....	9-16
Calibration of sensor .....	9-17



**Index****Date of issue**



# 1. Introduction

Names, intended use and limitations.....	1-2
Symbols used on the monitor .....	1-4
Symbols used in the manual.....	1-6

## Names, intended use and limitations

<b>Proprietary names</b>	TCM4 monitor and TCM40 monitor.
<b>Common names</b>	Transcutaneous $p\text{CO}_2/p\text{O}_2$ monitoring system (TCM4) and Transcutaneous $p\text{CO}_2/p\text{O}_2/\text{SpO}_2/\text{Pulse}$ monitoring system (TCM40)
<b>TCM4 series monitors</b>	The TCM4 and TCM40 monitoring systems are part of the TCM4 Series monitoring system.
<b>Reference</b>	This operator's manual for the TCM4/40 monitoring systems is intended for use as a reference. It provides detailed operating instructions and answers to relevant questions about your monitoring systems.
<b>All rights reserved</b>	At the time of printing, the manual is in conformity with the systems. All rights are reserved for instruments, circuits, techniques and names appearing in the manual.
<b>Intended use</b>	<p>The TCM4 monitoring system is intended for continuous transcutaneous monitoring of carbon dioxide (<math>\text{tcpCO}_2</math>) and oxygen (<math>\text{tcpO}_2</math>) partial pressures. It is indicated for use on neonates, pediatrics, and adults not under gas anesthesia.</p> <p>The TCM40 monitoring system is intended for continuous transcutaneous monitoring of carbon dioxide (<math>\text{tcpCO}_2</math>) and oxygen (<math>\text{tcpO}_2</math>) partial pressures as well as of oxygen saturation of arterial hemoglobin (<math>\text{SpO}_2</math>) and pulse rate. It is indicated for use on neonates, pediatrics and adults not under gas anesthesia.</p>
<b>Sensors</b>	<p>The sensors are provided non-sterile and are intended for reuse.</p> <p>For information on site locations, see the relevant section about application of sensors in chapter 6 <i>In vivo monitoring</i>.</p>
<b>Environment of use</b>	In hospital/clinical environment.
<b>Operator profile</b>	Only trained health care personnel are permitted to use the monitor.
<b>Limitations</b>	<p>Transcutaneous monitoring is intended only as an adjunct in patient assessment and must be used in conjunction with clinical signs and symptoms.</p> <div style="margin-top: 10px;">  <b>WARNING – Risk of incorrect measurements</b>  <math>\text{tcpCO}_2/\text{tcpO}_2</math> monitoring should not be used on patients in a compromised hemodynamic state as this may cause incorrect measurements. </div> <div style="margin-top: 10px;">  <b>WARNING – Risk of incorrect measurements</b>  The DS100A <math>\text{SpO}_2</math> sensor is contraindicated for use on active patients or for prolonged use. It is not designed for long-term monitoring. Using this sensor for long-term monitoring may result in incorrect measurements. </div>

---

*Continued on next page*

## Names, intended use and limitations, *Continued*

### Limitations (continued)

**WARNING – *Risk of allergic reactions***

The OXIband A/N and P/I SpO<sub>2</sub> sensors are contraindicated for use on patients who exhibit allergic reactions to the pressure-sensitive adhesive on the wraps.

**CAUTION – *US federal law restriction***











Federal law restricts this device to sale by or on the order of a physician.

**NOTICE:** This equipment is not a blood gas device.

### Legal notices

- Instruments should be repaired by authorized service personnel or by Radiometer-certified representatives only.
- Purchase of the TCM40 monitoring system confers no express or implied license under any Nellcor patent to use this instrument with any oximetry sensor that is not manufactured or licensed by Nellcor.




## Symbols used on the monitor

Symbol	Explanation
	CSA approved
	Caution, consult accompanying documents
	Temperature limitation
	Indicates that the product complies with the requirements of the Medical Device Directive 93/42/EEC June 1993. This product is a class IIa device.
	Serial number
	Non-ionizing radiation
	Type BF equipment (body floating)
	Manufacturer
	USB
	<p><b>Waste of Electrical and Electronic Equipment (WEEE)</b></p> <p>The symbol indicates that:</p> <ul style="list-style-type: none"> <li>• Radiometer Medical ApS and its distributors within the European Union (EU) and associated states have taken the necessary steps to comply with the directive 2002/96/EC on waste electrical and electronic equipment (WEEE)</li> <li>• The instrument, when reaching its end of life, must be collected and recycled separately from other waste according to national requirements. Please contact your local Radiometer distributor for instructions.</li> </ul> <p>Environmental implications: WEEE contains materials that are potentially hazardous to the environment and to human health.</p>

*Continued on next page*







## Symbols used on the monitor, *Continued*

Symbol	Explanation
IOIOI	COM gate
	Ethernet interface connection to network. Not for phone connection.
○	Off (Power: disconnection from the mains)
■	On (Power: connection to the mains)
	Monitor on/off
	Fuse

## Symbols used in the manual

This manual contains alerts, which are important and should be read carefully before performing the related procedures. The manual also contains non-safety information.

Symbol	Signal word	Explanation
	WARNING	A warning alerts the reader about a situation which, if not avoided, could result in death or serious injury. It may also describe potential serious adverse reactions and safety hazards.
	CAUTION	A caution alerts the reader about a potentially hazardous situation which, if not avoided, may result in minor or moderate injury to the user or the patient or damage to the equipment or other property. It may also be used to alert against unsafe practices. This includes the special care necessary for the safe and effective use of the device and the care necessary to avoid damage to the device that may occur as a result of use or misuse.
	NOTICE	Addresses practical information that is not related to personal injury ("need-to-know information").
		Manufacturer
		Indicates that the product complies with the requirements of the Medical Device Directive 93/42/EEC June 1993.  This product is a class IIa device.

## 2. What is what

TCM4/40 monitoring systems.....	2-2
Monitor – top and front.....	2-3
Monitor – rear .....	2-4
The screen: general elements .....	2-5
Online tutorials.....	2-8
Touch key glossary .....	2-9

## TCM4/40 monitoring systems

### Introduction

The TCM4 monitoring system includes:

- Base unit
- $\text{tcpCO}_2/\text{tcpO}_2$  module
- Combined  $\text{tcpCO}_2/\text{tcpO}_2$  sensors or single  $\text{tcpCO}_2$  sensor


The TCM40 monitoring system includes:

- Base unit
- $\text{tcpCO}_2/\text{tcpO}_2$  module and  $\text{SpO}_2$  module
- Sensors for  $\text{tcpCO}_2/\text{tcpO}_2$  module: Combined  $\text{tcpCO}_2/\text{tcpO}_2$  sensors or single  $\text{tcpCO}_2$  sensor
- Sensors for  $\text{SpO}_2$  module:  $\text{SpO}_2$  sensors (Nellcor DS100A, Nellcor Oxiband A/N or Nellcor Oxiband P/I)



**NOTICE:** For ordering information, see the section *Accessories* in chapter 8.

## Monitor – top and front

### Front view




### Parts and functions

Part	Function
Handle	For carrying the monitor.
Touch screen	For easy communication between operator and monitor.
Calibration chamber	For storage and calibration of tcpCO <sub>2</sub> /tcpO <sub>2</sub> sensor.
tcpCO <sub>2</sub> /tcpO <sub>2</sub> module	For transcutaneous measurements of tcpCO <sub>2</sub> and tcpO <sub>2</sub> and calibration of the sensor.
Monitor support and holder for sensor cable	For supporting the monitor and holding the sensor cable.
Monitor-elevating support	For supporting the monitor at an elevated angle.
 ON/OFF button	For turning the monitor ON and OFF.
 Battery indicator	For indicating whether the battery is being recharged or not. If the light is on, the battery is being recharged. <b>NOTICE:</b> Only possible if the power switch on the back of the monitor is in the ON position.
SpO <sub>2</sub> module (TCM40 monitor only)	For measurements of SpO <sub>2</sub> and pulse rate.

## Monitor – rear

### Rear view



### Parts and functions

Part	Function
Battery	For allowing monitoring during transport and power failure.
SpO <sub>2</sub> sensor socket (TCM40 monitor only)	For connecting an SpO <sub>2</sub> sensor to the monitor.
Line fuse	1.25 AT. For preventing a short circuit.  <div style="display: flex; align-items: center;"> <div> <b>WARNING – Risk of fire</b>            Replace fuse only as recommended by Radiometer.            Otherwise you risk that the monitor catches fire.         </div> </div>
Power socket	For connecting a power cord to the monitor.
Serial port (RS232)	For connecting the monitor to an external computer.
Analog output	For connecting the monitor to a peripheral device such as a polysomnograph.
USB ports	For connecting the monitor to an external printer and a memory stick.  <b>NOTICE:</b> It is only possible to use one of the ports at a time.
Module release	For releasing the sensor module from the monitor with a release key.
tcpCO <sub>2</sub> /tcpO <sub>2</sub> sensor socket	For connecting a tcpCO <sub>2</sub> /tcpO <sub>2</sub> sensor to the monitor.
Power switch	For turning the power supply to the monitor ON and OFF.  <b>NOTICE:</b> If the power supply is ON, the battery will be recharged when needed.


## The screen: general elements

### Screen types





There are two types of screens on the TCM4/40 monitors: view (i.e. Normal, Trend table and Trend Curve) and menu screens.

### View screen configuration

The view screens are divided into the following fields:




**NOTICE:** SpO<sub>2</sub> and pulse rate are available on the TCM40 monitor only. On the TCM4 monitor, the lower part of the curve display will either show the power curve or be empty, and the Pulse/SpO<sub>2</sub> display will always be empty (see examples on next page).

Part	Shows
Time and date	The real time (24-hour cycle) and date
System/view status	Normal view, Trend table view, Trend curve view, New patient, Alert or Alarm
Bat.	The battery level when the monitor is running on battery. Full  , almost full  , low  and critically low 
ID	Patient identification number. Gives access to Patient DMS.
Curve display	TCM4 monitor: $p\text{CO}_2$ , $p\text{O}_2$ and power TCM40 monitor: $p\text{CO}_2$ , $p\text{O}_2$ and SpO <sub>2</sub> Sensor status (e.g. Calibrating) and gas level (only displayed during calibration and if there is 10 % or less gas left in the gas cylinder). Barometric pressure (only displayed during calibration).
$p\text{CO}_2/p\text{O}_2$ display	$p\text{CO}_2$ , $p\text{O}_2$ , Corr., Power, Temp, SmartHeat, In vivo calibration active, ▲ (alarm is ON) and ✖ (alarm is OFF)
Pulse/SpO <sub>2</sub> display	Pulse, SpO <sub>2</sub> , "HI OFF" (SpO <sub>2</sub> alarm high is disabled), ▲ (alarm is ON) and ✖ (alarm is OFF)


Continued on next page

## The screen: general elements, *Continued*

**Example of  
Normal view  
screen on  
TCM40 monitor**



**Example of  
Normal view  
screen on TCM4  
monitor**



**NOTICE:** In Measuring mode, the *Event* touch key replaces the *Tutorial* key.

**Touch keys in  
Normal view**

Touch key	Function
<i>ID</i>	Gives access to Patient DMS.
<i>Site time xx:xx</i>	Resets the site timer to its preset value (see Parameter setup).
<i>Tutorial</i>	Gives access to instruction videos – when not monitoring.
<i>Event</i>	Marks an event during monitoring.
<i>Calibrate</i>	Starts a calibration of the sensor (and the SmartCal period, if the function is set to ON in the setup).
<i>Print</i>	Gives access to the Printer start/stop time screen.
<i>Setup</i>	Gives access to all the setup menus and submenus.
<i>Alarm silence</i>	Silences/resets the alarm system.


*Continued on next page*



## The screen: general elements, *Continued*


### Menu screen configuration

The menu screens are divided into the following fields:



The menu screens contain the following tabs, which give access to the setup menus and submenus: View, Parameter, Calibration, Printer and Technical (password protected).

### Example of menu screen



**NOTICE:** SpO<sub>2</sub> is available on the TCM40 monitor only.


### Menu screen touch keys

Touch key	Function
<b>OK</b>	Accepts the changes and returns to the main screen.
<b>Apply</b>	Accepts the changes without leaving the menu screen.
<b>Cancel</b>	Returns to the main screen without saving changes.

## Online tutorials







**Tutorials** The tutorials are short video sequences of commonly used procedures. To access the tutorials, press ***Tutorial***.

**NOTICE:** In Measuring mode, the ***Event*** touch key replaces the ***Tutorial*** key.



**NOTICE:** The DS100A sensor tutorial can only be selected on the TCM40 monitor.

### Touch keys



Touch key	Function
	Pauses the selected tutorial.
	Plays the selected tutorial.
 and 	Scroll backward/forward one step in the selected tutorial.
	Returns to the TCM tutorials screen.
	Exits the tutorial and returns to Normal view.

## Touch key glossary

### Touch keys

In the text, touch keys are written in ***bold italic*** throughout the manual.

The glossary of all the touch keys used in the software with their description is given in alphabetical order in the table below:

Touch key	Function
	Adds a check mark to the highlighted patient ID/session number (in Patient DMS).
	Displays detailed information about the highlighted patient ID/session number (in Patient DMS).
<b><i>Alarm silence</i></b>	Silences/resets the alarm system.
<b><i>Apply</i></b>	Accepts the changes made in the settings without leaving the respective menu screen.
<b><i>Blood gas</i></b>	Gives access to the Blood gas setup where blood gas values can be keyed in.
<b><i>Cal. status</i></b>	Gives information about the last calibration.
<b><i>Calibrate</i></b>	Starts a calibration of the electrode.
<b><i>Cancel</i></b>	Returns to the main screen without saving the changes made in the settings.
<b><i>Cursor</i></b>	Adds/removes a cursor in the Trend curve view.
<b><i>Date/time</i></b>	Gives access to the Date/time setup (password protected).
<b><i>Default values</i></b>	Changes all settings to Radiometer default values (password protected).
<b><i>Delete config. file</i></b>	Part of the Service setup (only for service purposes).
<b><i>Enter</i></b>	Registers the entered password.
<b><i>Event</i></b>	Marks an event during monitoring.
<b><i>Export</i></b>	Exports the patient information and measuring data for the patient IDs/session numbers with a check mark (in Patient DMS).
<b><i>ID</i></b>	Gives access to Patient DMS.
<b><i>In vivo calibration</i></b>	Starts an in vivo calibration, i.e. a correction of the measured tcpCO <sub>2</sub> /tcpO <sub>2</sub> values with the keyed-in blood gas values.
<b><i>Normal</i></b>	Gives access to the Normal view setup.

*Continued on next page*

## Touch key glossary, *Continued*











### Touch keys (*continued*)

Touch key	Function
<b><i>OK</i></b>	Accepts the changes made in the settings and returns to the main screen.
<b><i>pCO<sub>2</sub></i></b>	In the Parameter menu, it gives access to the <i>pCO<sub>2</sub></i> parameter setup. In the Trend curve setup, it gives access to the <i>pCO<sub>2</sub></i> curve range screen.
<b><i>pO<sub>2</sub></i></b>	In the Parameter menu, it gives access to the <i>pO<sub>2</sub></i> parameter setup. In the Trend curve setup, it gives access to the <i>pO<sub>2</sub></i> curve range screen.
<b><i>Power</i></b>	Gives access to the Power range screen.
<b><i>Print</i></b>	Gives access to the Printer start/stop time screen.
<b><i>Service menu</i></b>	Gives access to the Service setup (password protected).
<b><i>Setup</i></b>	Gives access to all the Setup menus.
<b><i>Site time</i></b>	Resets the site timer to its preset value (see Parameter setup).
<b><i>SmartCal</i></b>	Keeps the electrode ready for use (i.e. calibrated for max. 12 hours).
<b><i>SpO<sub>2</sub>/Pulse</i></b>	In the Parameter menu, it gives access to the SpO <sub>2</sub> /Pulse parameter setup. In the Trend curve setup, it gives access to the SpO <sub>2</sub> /Pulse curve range screen.
<b><i>System info</i></b>	Part of the Service setup (only for service purposes).
<b><i>Tech. settings</i></b>	Gives access to the Technical settings (password protected).
<b><i>Test</i></b>	Part of the Service setup (only for service purposes).
<b><i>Touch screen calibration</i></b>	Part of the Service setup (only for service purposes).
<b><i>Trend curve</i></b>	Gives access to the Trend curve setup.
<b><i>Trend table</i></b>	Gives access to the Trend table setup.
<b><i>Tutorial</i></b>	Gives access to instruction videos.







*Continued on next page*

## Touch key glossary, *Continued*

### Arrow touch keys

Touch key	Function
 or 	Changes the settings of the highlighted option.
 or 	Scrolls the displayed screen/parameters upward or downward.
 or 	Scrolls the displayed screen/parameters forward or backward.
	Scrolls quickly to the most recent result.
	Deletes one character at a time.
 or 	Moves one character at a time to the left or the right.

### Tutorial touch keys

Touch key	Function
	Pauses the selected tutorial.
	Plays the selected tutorial.
 and 	Scroll backward/forward one step in the selected tutorial.
	Returns to the TCM Tutorials screen.
	Exits the tutorial and returns to Normal view.




### 3. Menu structure and setup programs

Menu structure .....	3-2
List of setup programs.....	3-3
<b>View setup .....</b>	<b>3-4</b>
Normal view .....	3-5
Trend table view.....	3-7
Trend curve view .....	3-8
<b>Parameter setup .....</b>	<b>3-9</b>
$p\text{CO}_2$ .....	3-10
$p\text{O}_2$ .....	3-11
SpO <sub>2</sub> /Pulse .....	3-12
Blood gas.....	3-13
<b>Calibration setup .....</b>	<b>3-14</b>
SmartCal.....	3-15
Calibration status.....	3-16
<b>Printer setup .....</b>	<b>3-17</b>
<b>Technical setup .....</b>	<b>3-19</b>
Technical settings.....	3-20
Date/time .....	3-22
Default values .....	3-23

## Menu structure

**Menu structure** The following diagram illustrates the menu structure of the TCM4/40 monitors.



\* Available on the TCM40 monitor only.

\*\* Available on the TCM4 monitor only.

**NOTICE:** If no touch key is activated for 30 seconds, the main screen (Normal view) is displayed.



## List of setup programs

### Accessing the Setup menus

Press the **Setup** touch key to get access to the Setup menus:

- View
- Parameter
- Calibration
- Printer
- Technical

Detailed information about the five main Setup menus is given in the following sections.

## View setup

Normal view .....	3-5
Trend table view .....	3-7
Trend curve view .....	3-8

## Normal view

### Introduction


In the Normal view setup, it is possible to select which parameters will be shown in Normal view, to set the time span of the graphical display and to activate the screen saver.

**NOTICE:** If no touch key is activated for 30 seconds, the main screen (Normal view) is displayed.


### Accessing the program

To access the Normal view setup, press **Setup** → **Normal**.

Normal view setup on TCM40 monitor:



Normal view setup on TCM4 monitor:



### Selecting Normal view options

Step	Action
1.	Use the <b>Up</b> and <b>Down</b> arrows to select the curves that are to be shown in the curve area in Normal view.
<b>NOTICES:</b>	
<ul style="list-style-type: none"> <li>Parameter values will always be shown.</li> <li>The number of parameters to choose from depends on the sensor type installed (combined <math>p\text{CO}_2/p\text{O}_2</math> or <math>p\text{CO}_2</math> only), and whether the <math>p\text{O}_2</math> parameter is set to ON or OFF in Technical setup.</li> </ul>	

*Continued on next page*

## Normal view, *Continued*

Selecting Normal view options ( <i>continued</i> )	Step	Action
	2.	Use the <i>Up</i> and <i>Down</i> arrows to select the time span.
		<b>NOTICE:</b> The time span selection only applies to the $p\text{CO}_2$ , $p\text{O}_2$ and Power curves.
	3.	Select the curve range options of the individual parameters by pressing each of the parameters and, in the appearing screens, selecting the high and low values.
	4.	Press <b>OK</b> to accept the changes and return to the Normal view setup screen, or press <b>Cancel</b> to return to the Normal view setup screen without saving changes.
	5.	If required, select (✓) the screen saver.
		<b>NOTICES:</b> <ul style="list-style-type: none"> <li>• To activate the screen saver, it is necessary to press <b>OK</b> before leaving the Normal view setup.</li> <li>• The option is only visible if it has been activated in Technical settings.</li> <li>• The screen saver is primarily meant for sleep labs, to reduce the backlight from the display.</li> <li>• Touching the screen deactivates the screen saver; and to reactivate it, the option must be selected in Normal view setup.</li> </ul>
	6.	Press:
		• <b>OK</b> to accept the changes and return to the main screen
		• <b>Apply</b> to accept the changes without leaving the Normal view setup
		• <b>Cancel</b> to return to the main screen without saving changes

## Trend table view

**Introduction** In the Trend table setup, it is possible to select the time interval between each record in the Trend table.

**Accessing the program** To access the Trend table setup, press *Setup* → *Trend table*.

Selecting Trend table options	Step	Action
	1.	Use the <i>Up</i> and <i>Down</i> arrows to select the time interval.
	2.	Press: <ul style="list-style-type: none"><li>• <i>OK</i> to accept the changes and return to the main screen</li><li>• <i>Apply</i> to accept the changes without leaving the Trend table setup</li><li>• <i>Cancel</i> to return to the main screen without saving changes</li></ul>


## Trend curve view

**Introduction** In the Trend curve setup, it is possible to set the time span of the graphic display and the ranges for  $p\text{CO}_2$ ,  $p\text{O}_2$ , Power and  $\text{SpO}_2$ /Pulse, and to select which parameter(s) (maximum two) should be displayed on the Trend curve.

**NOTICE:**  $\text{SpO}_2$  and pulse are available on the TCM40 monitor only.

### Accessing the program

To access the Trend curve setup, press **Setup** → **Trend curve**.



### Selecting Trend curve options

- | Step | Action   |
|------|--|
| 1.   | Use the <b>Up</b> and <b>Down</b> arrows to select the time span.  |
| 2.   | Select maximum two parameters to be displayed as trend curves.<br><b>NOTICE:</b> If two parameters are selected, one of these must be deselected in order to be able to select a new parameter.  |
| 3.   | Select the curve range options of the individual parameters by pressing each of the parameters and, in the appearing screens, selecting the high and low values, using the arrow keys.   |
| 4.   | Press <b>OK</b> to accept the changes and return to the Trend curve setup screen, or press <b>Cancel</b> to return to the Trend curve setup screen without saving changes.   |
| 5.   | In the Trend curve setup screen, press: <ul style="list-style-type: none"> <li>• <b>OK</b> to accept the changes and return to the main screen</li> <li>• <b>Apply</b> to accept the changes without leaving the Trend curve setup</li> <li>• <b>Cancel</b> to return to the main screen without saving changes</li> </ul> |

## Parameter setup

**WARNING – *Risk of incorrect monitoring***

Make sure to select the alarm limits carefully. Setting alarm limits to extreme values can render the alarm system useless.

$p\text{CO}_2$ .....	3-10
$p\text{O}_2$ .....	3-11
$\text{SpO}_2/\text{Pulse}$ .....	3-12
Blood gas.....	3-13

## **pCO<sub>2</sub>**

### Accessing the program

To access the pCO<sub>2</sub> setup, press **Setup** → **Parameter** → **pCO<sub>2</sub>**.

### pCO<sub>2</sub> settings

Settings	Options
pCO <sub>2</sub> alarm	ON/OFF <b>NOTICE:</b> The alarm is triggered if the parameter value exceeds or is equal to the alarm limit (high or low). It consists of a visual part (the parameter value and an alarm text will be flashing) and an acoustic part (a discontinuous tone).
pCO <sub>2</sub> alarm high	6-200 mmHg (in steps of 1); 0.8-26.7 kPa (in steps of 0.1)
pCO <sub>2</sub> alarm low	5-99 mmHg (in steps of 1); 0.7-9.9 kPa (in steps of 0.1)
Alarm sound level <sup>1</sup>	1-5 (1 is minimum)
Sensor temp. <sup>2</sup>	37.0-45.0 °C (in steps of 0.5) <b>NOTICE:</b> If the set sensor temperature is changed during measurement, a new calibration is required.
SmartHeat <sup>2</sup>	ON/OFF <b>NOTICE:</b> If SmartHeat is set to ON, it adds +1 °C (max. temp. 45 °C) to the set sensor temperature for 5 minutes after the sensor has been removed from the calibration chamber.
Site time <sup>2</sup>	OFF; ½-12 hours (in steps of ½ hour)
Site time heat <sup>2</sup>	ON/OFF <b>NOTICE:</b> If Site time heat is set to OFF, the sensor heat is switched off when the site timer reaches zero; if set to ON, the heat continues.

<sup>1</sup> The setting of this option is common to pCO<sub>2</sub>, pO<sub>2</sub> and SpO<sub>2</sub>.

<sup>2</sup> The setting of this option is common to pCO<sub>2</sub> and pO<sub>2</sub>.

### Selecting pCO<sub>2</sub> settings

Step	Action
1.	Select the relevant option with the ▼ touch key.
2.	Choose the settings of that option with the <b>Up</b> or <b>Down</b> arrow key.
3.	Follow steps 1-2 for the remaining options.
4.	Press: <ul style="list-style-type: none"> <li>• <b>OK</b> to accept the changes and return to the main screen</li> <li>• <b>Apply</b> to accept the changes without leaving the pCO<sub>2</sub> setup</li> <li>• <b>Cancel</b> to return to the main screen without saving changes</li> </ul>



## ***pO<sub>2</sub>***

### **Accessing the program**

To access the *pO<sub>2</sub>* setup, press **Setup** → **Parameter** → ***pO<sub>2</sub>***.

### ***pO<sub>2</sub>* settings**

<b>Settings</b>	<b>Options</b>
<i>pO<sub>2</sub></i> alarm	ON/OFF  <b>NOTICE:</b> The alarm is triggered if the parameter value exceeds or is equal to the alarm limit (high or low). It consists of a visual part (the parameter value and an alarm text will be flashing) and an acoustic part (a discontinuous tone).
<i>pO<sub>2</sub></i> alarm high	1-800 mmHg (in steps of 1); 0.1-99.9 kPa (in steps of 0.1)
<i>pO<sub>2</sub></i> alarm low	0-99 mmHg (in steps of 1); 0.0-9.9 kPa (in steps of 0.1)
Alarm sound level <sup>1</sup>	1-5 (1 is minimum)
Sensor temp. <sup>2</sup>	37.0-45.0 °C (in steps of 0.5)  <b>NOTICE:</b> If the set sensor temperature is changed during measurement, a new calibration is required.
SmartHeat <sup>2</sup>	ON/OFF  <b>NOTICE:</b> If SmartHeat is set to ON, it adds +1 °C (max. temp. 45 °C) to the set sensor temperature for 5 minutes after the sensor has been removed from the calibration chamber.
Site time <sup>2</sup>	OFF; ½-12 hours (in steps of ½ hour)
Site time heat <sup>2</sup>	ON/OFF  <b>NOTICE:</b> If Site time heat is set to OFF, the sensor heat is switched off when the site timer reaches zero; if set to ON, the heat continues.

<sup>1</sup> The setting of this option is common to *pCO<sub>2</sub>*, *pO<sub>2</sub>* and *SpO<sub>2</sub>*.

<sup>2</sup> The settings of these options are common to *pCO<sub>2</sub>* and *pO<sub>2</sub>*.

### **Selecting *pO<sub>2</sub>* settings**


<b>Step</b>	<b>Action</b>
1.	Select the relevant option with the ▼ touch key.
2.	Choose the settings of that option with the <b>Up</b> or <b>Down</b> arrow key.
3.	Follow steps 1-2 for the remaining options.
4.	Press: <ul style="list-style-type: none"> <li>• <b>OK</b> to accept the changes and return to the main screen</li> <li>• <b>Apply</b> to accept the changes without leaving the <i>pO<sub>2</sub></i> setup</li> <li>• <b>Cancel</b> to return to the main screen without saving changes</li> </ul>

## SpO<sub>2</sub>/Pulse

### Accessing the program

To access the SpO<sub>2</sub>/Pulse setup, press **Setup** → **Parameter** → **SpO<sub>2</sub>/Pulse**.

### SpO<sub>2</sub>/Pulse settings

Settings	Options
SpO <sub>2</sub> alarm	ON/OFF  <b>NOTICE:</b> The alarm is triggered if the parameter value exceeds or is equal to the alarm limit (high or low). It consists of a visual part (the parameter value and an alarm text will be flashing) and an acoustic part (a discontinuous tone).
SpO <sub>2</sub> alarm high	OFF/86-100 % (in steps of 1)   <b>WARNING – Risk of patient injury</b> Make sure to select the upper alarm limit for oxygen saturation carefully and in accord with accepted clinical standards. High oxygen levels may predispose a premature infant to develop retinopathy.  <b>NOTICE:</b> SpO <sub>2</sub> alarm high can be disabled ("HI OFF" will be displayed beneath the alarm symbol) while SpO <sub>2</sub> alarm low is kept active.
SpO <sub>2</sub> alarm low	85-99 % (in steps of 1)
SatSeconds	OFF/10-100 (in steps of 10)
Alarm sound level <sup>1</sup>	1-5 (in steps of 1)
Pulse alarm	ON/OFF
Pulse alarm high	35-240 bpm (in steps of 5)
Pulse alarm low	30-235 bpm (in steps of 5)

<sup>1</sup> The setting of this option is common to pCO<sub>2</sub>, pO<sub>2</sub> and SpO<sub>2</sub>.

### Selecting SpO<sub>2</sub>/Pulse settings

Step	Action
1.	Select the relevant option with the ▼ touch key.
2.	Choose the settings of that option with the <b>Up</b> or <b>Down</b> arrow key.
3.	Follow steps 1-2 for the remaining options.
4.	Press: <ul style="list-style-type: none"> <li>• <b>OK</b> to accept the changes and return to the main screen</li> <li>• <b>Apply</b> to accept the changes without leaving the SpO<sub>2</sub>/Pulse setup</li> <li>• <b>Cancel</b> to return to the main screen without saving changes</li> </ul>

## Blood gas

### Introduction

In the Blood gas setup it is possible to:

- key in the blood gas values of a patient to compare these with transcutaneous measurements from the same patient. The blood gas values are displayed as blood drops in all views, as well as on all printed reports.
- calibrate the transcutaneous measurements against the blood gas values (in vivo calibration), and all results will then be displayed as red stars.

### Accessing the program

To access the Blood gas setup, press **Setup** → **Parameter** → **Blood gas**.

**NOTICE:** The system must be monitoring.

**NOTICE:** The ***In vivo calibration*** touch key is only visible if the option has been activated in Technical settings (see *Technical settings* later in this chapter).

### Procedures

For the procedures on how to key in a blood gas value and how to perform an in vivo calibration, see *Blood gas comparison* and *In vivo calibration* in chapter 6.

## Calibration setup

SmartCal.....	3-15
Calibration status.....	3-16

## SmartCal

**Purpose** The SmartCal function makes sure that the monitor is always ready for monitoring by calibrating the sensor when needed.

**Accessing the program** To access the SmartCal setup, press *Setup* → *Calibration* → *SmartCal*.

### SmartCal settings

Settings	Options	Function
SmartCal	ON/OFF	Enables/disables the SmartCal function. <b>NOTICES:</b> <ul style="list-style-type: none"> <li>When enabled, the monitor will calibrate automatically. The calibration interval may vary from 5 minutes to 2 hours; and the maximum measuring time will be reduced accordingly.</li> <li>Although the option is enabled, it is necessary to press <i>Calibrate</i> to start a SmartCal period.</li> </ul>
SmartCal duration	Forever/ 1-12 hours (in steps of 1)	Defines the duration of the SmartCal period. <b>NOTICE:</b> If set to forever, the monitor will calibrate whenever the sensor is placed in the calibration chamber; if set to 1-12 hours, the monitor will calibrate within the selected period, and after a SmartCal period, it will be necessary to press <i>Calibrate</i> to start a new SmartCal period.
Elapsed time	0:00-12:00 (hours:minutes)	Shows how much time has elapsed of the active SmartCal period.

### Selecting SmartCal settings

Step	Action
1.	Select the relevant option with the ▼ touch key.
2.	Choose the settings of that option with the <i>Up</i> or <i>Down</i> arrow key.
3.	Follow steps 1-2 for the remaining options.
4.	Press: <ul style="list-style-type: none"> <li><i>OK</i> to accept the changes and return to the main screen</li> <li><i>Apply</i> to accept the changes without leaving the SmartCal setup</li> <li><i>Cancel</i> to return to the main screen without saving changes</li> </ul>

## Calibration status

**Introduction**      The Calibration status screen shows the status of the last calibration.

**Accessing the program**      To access the Calibration status screen, press *Setup* → *Calibration* → *Cal. Status*.

### Calibration status information

Status info	Unit	Description
Last cal.	hour:minutes	Shows the time of the last calibration.
Set temp.	°C	Shows the sensor temperature during the last calibration.
Barometer	mmHg/kPa	Shows the barometric pressure during the last calibration.
Cal. value O <sub>2</sub>	mmHg/kPa	Shows the O <sub>2</sub> calibration value.
Cal. value CO <sub>2</sub>	mmHg/kPa	Shows the CO <sub>2</sub> calibration value.
Gas level	%	Shows how much calibration gas is left in the cylinder.  <b>NOTICE:</b> The level will be displayed as "High" until there is less than 10 % left in the gas cylinder, and then as a percentage.

Press *OK* or *Cancel* to return to the main screen.

## Printer setup

### Introduction

In the Printer setup, it is possible to:

- set the time interval between two values on the table printout
- set the time span of the measurements that are to be shown on the printout
- choose whether to print out the report as a table, a curve or both
- adjust the set curve ranges
- connect a color or a black-and-white (B/W) printer to the monitor

### Accessing the program

To access the Printer setup, press **Setup** → **Printer**.

**NOTICE:** SpO<sub>2</sub> and pulse are available on the TCM40 monitor only.

### Selecting printer settings

Step	Action
1.	Select the time interval with the <b>Up</b> and <b>Down</b> arrow keys.
2.	Select the time span with the <b>Up</b> and <b>Down</b> arrow keys.  <b>NOTICE:</b> The selected time span will influence the start time of the Printer start/stop time screen, as the interval between the printer start time and the printer stop time corresponds to the time span; i.e. if the time span is set to 1 hour, the interval between the printer start time and the printer stop time will also be 1 hour (see also <i>How to print</i> in chapter 6: <i>In vivo monitoring</i> ).
3.	To change the curve range for pCO <sub>2</sub> , pO <sub>2</sub> , Power or SpO <sub>2</sub> /Pulse, press the relevant parameter touch key and select the high and low values, using the <b>Up</b> and <b>Down</b> arrows.

*Continued on next page*

## Printer setup, *Continued*

Selecting printer settings ( <i>continued</i> )	Step	Action
	4.	Press <b>OK</b> to accept the changes and return to the Printer setup, or press <b>Cancel</b> to return to the Printer setup without saving the changes.
	5.	Select at least one report type: Table and/or Curve.
	6.	Select printer type: color or B/W (i.e. black and white).
		<b>NOTICES:</b>
		<ul style="list-style-type: none"><li>• It is only possible to connect a local printer to the monitor. Printing over the network is not supported.</li><li>• Radiometer recommends that you use an HP printer with PCL3 protocol.</li></ul>
	7.	Press: <ul style="list-style-type: none"><li>• <b>OK</b> to accept the changes and return to the main screen</li><li>• <b>Apply</b> to accept the changes without leaving the Printer setup</li><li>• <b>Cancel</b> to return to the main screen without saving changes</li></ul>



## Technical setup

Technical settings.....	3-20
Date/time .....	3-22
Default values .....	3-23

## Technical settings

### Accessing the program

To access the Technical settings, press **Setup** → **Technical** (enter password and press **Enter**) → **Tech. Settings**.

### Technical settings

Settings	Options	Function
Meta. corr. factor	0-15 mmHg (in steps of 1 mmHg); 0-1 kPa (in steps of 0.1 kPa)	Defines the metabolic correction factor. To indicate that CO <sub>2</sub> values have been corrected, "Corr." is displayed together with the pCO <sub>2</sub> value in Normal view.
Severinghaus corr.	ON/OFF	If set to ON, all CO <sub>2</sub> values are corrected with the Severinghaus correction factor, and "Corr." is displayed together with the pCO <sub>2</sub> value in Normal view.
In vivo calibration	ON/OFF	Enables/disables access to the <b>In vivo calibration</b> touch key in the Blood gas setup.
Cal. gas mix pO <sub>2</sub>	0.0-100.0 % (in steps of 0.1)	A calibration constant.
Cal. gas mix pCO <sub>2</sub>	0.0-10.0 % (in steps of 0.1)	A calibration constant.
Unit pCO <sub>2</sub> /pO <sub>2</sub>	mmHg/kPa	Defines the pCO <sub>2</sub> /pO <sub>2</sub> unit.
pO <sub>2</sub> parameter	ON/OFF	Defines whether to display the pO <sub>2</sub> parameter or not.  <b>NOTICE:</b> This option only has an effect when a combined tcpCO <sub>2</sub> /tcpO <sub>2</sub> sensor is connected to the monitor.
Alarm mode	Latching/ non-latching	Defines whether the alarm is latching (the monitor remains in alarm status even though the alarm condition ceases to exist) or non-latching (the monitor resets itself as soon as the alarm condition ceases to exist).
Continuous data output	OFF, Standard, VueLink, MonLink or Raw data	Gives four possibilities for data output: standard, VueLink, MonLink and raw data.  See chapter 6 for detailed information.  <b>NOTICE:</b> Raw data is for service purposes only. For more information, see the <i>TCM4 Series service manual</i> .

*Continued on next page*

## Technical settings, *Continued*

### Technical settings (*continued*)

Option	Range (default)	Function
Data export	OFF, USB or Serial	Enables access to the <b>Export</b> touch key (in Patient DMS), which is used to export a dump of the trend data to a memory stick or an external PC.  <b>NOTICE:</b> "Serial" is only available if "Continuous data output" is set to OFF.
Data export interval	2, 10, 30 or 60 seconds	Defines the interval between the export of data.
Display brightness	10-100 % (in steps of 10)	Defines the backlight brightness of the display.  <b>NOTICE:</b> Full backlight reduces the lifetime of the display.
Screen saver	ON/OFF	Makes a screen saver option available in Normal view setup.
$pO_2$ analog range	0-200 or 0-800 mmHg; 0.0-26.7 or 0.0-99.9 kPa	Defines the $pO_2$ analog output range.
$pCO_2$ analog range	0-100 or 0-200 mmHg; 0.0-13.3 or 0.0-26.7 kPa	Defines the $pCO_2$ analog output range.

### Selecting technical settings

Step	Action
1.	Select the relevant option with the ▼ touch key.
2.	Choose the settings of that option with the <b>Up</b> or <b>Down</b> arrow key.
3.	Follow steps 1-2 for the remaining options.
4.	Press: <ul style="list-style-type: none"> <li>• <b>OK</b> to accept the changes and return to the main screen</li> <li>• <b>Apply</b> to accept the changes without leaving the Technical settings</li> <li>• <b>Cancel</b> to return to the main screen without saving changes</li> </ul>

## Date/time

### Accessing the program

To access the Date/time settings, press **Setup** → **Technical** (enter password and press **Enter**) → **Date/time**.



#### **WARNING – Risk of data loss**

If the Date/time settings are changed backward in time, only the measurements that have been performed *prior* to the new date/time will be kept in the memory. Other data will be deleted.

**NOTICE:** Date and time cannot be set during measurement. Trying to do so will result in an alert.

### Changing date/time settings

The "Date/time settings" input field reflects the current settings.

- | Step | Action   |
|------|--|
| 1.   | Select date format: <ul style="list-style-type: none"> <li>• MM-DD-YYYY (month-day-year)</li> <li>• DD-MM-YYYY (day-month-year)</li> </ul>   |
| 2.   | Move one character at a time in the "Date/time settings" input field with the << and >> touch keys. Enter the new settings with the numeric keypad.  |
| 3.   | Press: <ul style="list-style-type: none"> <li>• <b>OK</b> to accept the changes and return to the main screen</li> <li>• <b>Apply</b> to accept the changes without leaving the Date/time settings</li> <li>• <b>Cancel</b> to return to the main screen without saving changes</li> </ul> |

## Default values

**Introduction** The Default values function changes all parameter and monitor settings to factory defaults. The default values are listed below.

**NOTICE:** When changing all settings to factory defaults, all kPa settings will change to mmHg.

### Changing settings to factory defaults

Step	Action
1.	Press <b>Setup</b> → <b>Technical</b> (enter password and press <b>Enter</b> ) → <b>Default values</b> .
2.	A dialog box with the text "This will return the monitor into default setup and current setup will be lost" appears.  Press <b>OK</b> to change all parameter and monitor settings to factory defaults or <b>Cancel</b> to exit without changing settings.

### List of default settings

Menus	Settings	Default
Normal view	Curve selection	TCM40 monitor: $p\text{CO}_2/p\text{O}_2/\text{SpO}_2$
		TCM4 monitor: $p\text{CO}_2/p\text{O}_2/\text{Power}$
	Time span	1 hour
Table view	Time interval	15 minutes
Curve view	Time span	8 hours
	Curve selection	TCM40 monitor: $p\text{CO}_2$ and $\text{SpO}_2$
		TCM4 monitor: $p\text{CO}_2$ and $p\text{O}_2$
Curve range	$p\text{CO}_2$ high	60 mmHg/8.0 kPa
	$p\text{CO}_2$ low	0 mmHg/0.0 kPa
	$p\text{O}_2$ high	200 mmHg/26.7 kPa
	$p\text{O}_2$ low	0 mmHg/0.0 kPa
	Power high	400 mW
	Power low	0 mW
	$\text{SpO}_2$ high	100 %
	$\text{SpO}_2$ low	90 %
	Pulse high	200 bpm
	Pulse low	0 bpm

*Continued on next page*

## Default values, *Continued*

### List of default settings (*continued*)

Menus	Settings	Default
$p\text{CO}_2$	$p\text{CO}_2$ alarm	ON
	$p\text{CO}_2$ alarm high	50 mmHg/6.7 kPa
	$p\text{CO}_2$ alarm low	30 mmHg/4.0 kPa
	Alarm sound level	2
	Sensor temp.	43.0 °C
	SmartHeat	OFF
	Site time	4 hours
	Site time heat	OFF
$p\text{O}_2$	$p\text{O}_2$ alarm	ON
	$p\text{O}_2$ alarm high	95 mmHg/12.7 kPa
	$p\text{O}_2$ alarm low	60 mmHg/8.0 kPa
	Alarm sound level	2
	Sensor temp.	43.0 °C
	SmartHeat	OFF
	Site time	4 hours
	Site time heat	OFF
SpO <sub>2</sub> /Pulse	SpO <sub>2</sub> alarm	ON
	SpO <sub>2</sub> alarm high	OFF
	SpO <sub>2</sub> alarm low	85 %
	SatSeconds	OFF
	Alarm sound level	2
	Pulse alarm	ON
	Pulse alarm high	170 bpm
	Pulse alarm low	40 bpm
Calibration	SmartCal	OFF
	SmartCal duration	4 hours

*Continued on next page*

## Default values, *Continued*

### List of default settings (*continued*)

Menus	Settings	Default
Printer	Time interval	15 minutes
	Time span	1 hour
	Report type	Table and curve
	Printer type	Color
Technical	Meta. corr. factor	7 mmHg/1.0 kPa
	Severinghaus corr.	ON
	In vivo calibration	OFF
	Cal. gas mix $pO_2$	20.9 %
	Cal. gas mix $pCO_2$	7.5 %
	Unit $pCO_2/pO_2$	mmHg
	$pO_2$ parameter	ON
	Alarm mode	Non-latching
	Continuous data output	OFF
	Data export	OFF
	Data export interval	10 seconds
	Display brightness	70 %
	Screen saver	OFF
	$pCO_2$ analog range	0-200 mmHg/0.0-26.7 kPa
	$pO_2$ analog range	0-800 mmHg/0.0-99.9 kPa
Date/time	Date format	DD-MM-YYYY





## 4. Installation and maintenance

Operating requirements .....	4-2
Installing the monitoring system.....	4-3
Shutting down the monitor.....	4-5
Cleaning the monitor.....	4-6
Maintenance of the monitor .....	4-7
<b>Preparation and maintenance of the tc sensors.....</b>	<b>4-9</b>
General information about remembraning .....	4-10
Cleaning the sensor head of the E5480 sensor.....	4-11
Membraning the E5480 sensor.....	4-12
Membraning the E5280 and E5260 sensors.....	4-13
Cleaning, disinfection and storage of sensors .....	4-15
<b>Maintenance of the SpO<sub>2</sub> sensors.....</b>	<b>4-16</b>

## Operating requirements

**WARNING – Risk of incorrect measurements**

Do not use the monitor adjacent to or stacked with other equipment as these can cause electromagnetic interference and thereby result in incorrect measurements. If stacking or use adjacent to other equipment is necessary, the monitor should be observed to verify normal operation before used on patients. See the section *EMC approvals and compliance* in chapter 8.

**WARNING – Risk of incorrect measurements**

When installing, operating or servicing the monitor, special consideration should be given to the information regarding the electromagnetic precautions for this equipment given in the section *EMC approvals and compliance* in chapter 8. Otherwise the monitor may be affected by electromagnetic interference, causing incorrect measurements.

**WARNING – Risk of fire**

Do not place the monitor in an enriched oxygen atmosphere or inside a hyperbaric chamber as it may cause a fire hazard.

**WARNING – Risk of explosion**

Do not use the monitor in the presence of flammable anesthetics or gases as it may cause an explosion.

**WARNING – Risk of fire**

Do not cover the ventilator as this may cause it to seize up.

**WARNING – Risk of electrical shock**

Do not use the monitor if it is damaged as this may result in electrical shock.

**Environmental requirements**

The following environmental requirements must be observed at all times:


- ambient temperature: 5-40 °C
- relative humidity: 20-80 %.

**NOTICE:** Do not operate the monitor at ambient temperatures below 5 °C or above 40 °C and relative humidity below 20 % or above 80 %. Operating the monitor outside these limits may affect the readings of the device.

**Ventilation requirements**

The monitor must be located in a well-ventilated dust-free atmosphere.

## Installing the monitoring system

Procedure	Step	Action
	1.	<p><b>NOTICE:</b> While installing the modules, the monitor must be turned off.</p> <p>Install the relevant module(s) (i.e. tc and SpO<sub>2</sub>) by pressing each module completely into the base unit until a click is heard.</p> <p><b>NOTICE:</b> If using an E5480 sensor, the gasket in the calibration chamber must be code no. 837-488; and if using an E5280 or E5260 sensor, the gasket in the calibration chamber must be code no. 837-159 (see package insert for instructions).</p>
	2.	<p>For TCM40 monitor only: Plug the SpO<sub>2</sub> sensor into the monitor.</p> <p><b>NOTICE:</b> It is not necessary for the monitor to be turned off while the sensor is being connected.</p>
	3.	<p>Install the calibration gas cylinder and the battery according to the procedures described later in this chapter under <i>Maintenance of the monitor</i>.</p>
	4.	<p>Connect the monitor power cord to</p> <ul style="list-style-type: none"> <li>• the power socket at the rear of the monitor and</li> <li>• an appropriate power supply</li> </ul>
	5.	<p>Connect the system to external equipment, if required.</p> <p> <b>WARNING – Risk of personal injury</b></p> <p>Before connecting other equipment to the TCM monitor, the manufacturer of the equipment or a qualified engineer must be consulted to ensure that the equipment is compatible and that the safety of the patient, the operator or the environment will not be impaired. The resulting combined system must comply with EN 60601-1-1.</p>
	6.	<p>Turn on the monitor by pressing the power switch to the ON position at the rear of the monitor and then pressing the ON/OFF button on the front of the monitor.</p> <p><b>NOTICE:</b> Every time the monitor is turned on, a beep sounds, indicating that it has been checked that the sound of the alarm system is working. If the sound is not working, an error message is shown.</p>
	7.	<p>Check that the date and time in the display correspond with the actual date and time. Otherwise correct them in Technical setup.</p>
	8.	<p>Membrane the tcpCO<sub>2</sub>/tcpO<sub>2</sub> sensor as described later in this chapter.</p>
	9.	<p>Connect the tcpCO<sub>2</sub>/tcpO<sub>2</sub> sensor plug to the sensor socket at the rear of the tc module, and place the sensor in the calibration chamber at the front.</p>

*Continued on next page*


## Installing the monitoring system, *Continued*


Procedure ( <i>continued</i> )	Step	Action
	10.	Check that the tc module is functioning: The message "Calibration required" is displayed on the screen. Leave it until step 12.
	11.	Change View or Setup settings, if required, by pressing <b>Setup</b> . See chapter 3: <i>Menu structure and setup programs</i> .
		<b>NOTICE:</b> The monitor is delivered with default settings (see these in the section <i>Default values</i> in chapter 3).
	12.	Calibrate the sensor as described in chapter 5: <i>Calibration</i> .

## Shutting down the monitor

### Shutting down the monitor

Step	Action
------	--------

- |    |  |
|----|--|
| 1. | Press the  button on the front of the monitor.<br>The following system message will be displayed: |
|----|--|



- |    |   |
|----|---|
| 2. | Press <b>OK</b> to shut down the monitor – the message "Saving data. Please wait." will be displayed – or <b>Cancel</b> to return to the main screen without shutting down the monitor. |
|----|---|

## Cleaning the monitor

### Cleaning the exterior

When cleaning the monitor:

- Shut down the monitor by following the procedure described earlier in this chapter
- Use a cloth that is lightly dampened with soapy water
- Do not use abrasive cleansers or pads: the finish may become damaged
- Do not use aggressive detergents. Extensive use may cause the plastic to become brittle and cracks may occur.

### Cleaning the touch screen

A dry or lightly dampened soft, lint-free cloth may be used to clean the monitor's touch screen. Simply wipe the screen gently to remove fingerprints and/or dirt. To avoid streaking, an approved screen cleaner is recommended.

### Disinfection of outer surfaces

Disinfection of the monitor exterior and touch screen is performed when appropriate. The disinfection frequency depends on local requirements and the use of the monitor.

**NOTICE:** Follow legal requirements and local rules for safe work practices with chemicals.

The following disinfectants, dissolved in water, may be used to disinfect the monitor exterior and touch screen:

- 70 % isopropyl alcohol
- 70 % ethanol
- 4 % Diversol BX solution

To disinfect the monitor exterior and touch screen, wipe the surfaces using one of the above disinfectants on a paper towel or tissue.

### NOTICE

Do not spray, pour or spill any liquid on the monitor or any of the accessories, connectors, switches or openings in the chassis.

## Maintenance of the monitor

### Battery

The following battery type must be installed on the monitor: Standard 12 V 2 Ah lead-acid battery (code no. 431-018, available from Radiometer)

Contact your local battery supplier for available type.



**CAUTION – *Risk of patient not being monitored***

For data safety reasons, a battery must always be connected to the system.



**CAUTION – *Risk of patient not being monitored***

Replace battery only with the types recommended by Radiometer.



**CAUTION – *Handling of biohazardous waste***

Dispose of the battery according to local procedures to avoid personal injury or pollution of the environment.

The monitor can function for approx. one hour on battery supply (depending on battery type). When the monitors run on battery, it is indicated on the display, and the battery level indicator shows the battery level. If the battery level is low, an alert sound is heard and the battery level indicator is flashing. If the battery level is critically low, an alert message will furthermore be displayed.



**CAUTION – *Risk of patient not being monitored***

Make sure the battery level never becomes critically low, as this will prevent data from being saved onto the disk.

To avoid a complete discharge of the battery, reconnect the monitor to the mains as soon as possible to recharge the battery. Recharging the battery takes approximately 8 hours.


### Changing the battery

Step	Action
1.	Unscrew the battery cover at the rear of the monitor.
2.	Detach the battery from the connector.
3.	Attach a new battery to the connector.
4.	Reattach the battery cover.

*Continued on next page*

## Maintenance of the monitor, *Continued*

### Changing the calibration gas

Step	Action
1.	<p>Unscrew the old calibration gas cylinder.</p> <p> <b>WARNING – Risk of explosion</b>            Calibration gas cylinder: Contents under pressure. Do not puncture. Do not use or store near heat or open flame. Exposure to temperatures above 54 °C (for CAL2) and 50 °C (for CAL1) may cause contents to vent or cause bursting. Never discard container into fire or incinerator as it may cause an explosion.</p>
2.	<p>Screw the new calibration gas cylinder clockwise as far as possible into the socket.</p> <p><b>NOTICE:</b> Excessive force may damage threads and result in a leakage and thereby an increase in consumption.</p>

### Changing the gasket in the calibration chamber

Step	Action
1.	Remove the gasket from the calibration chamber.
2.	Clean the calibration chamber with a soft cloth moistened with skin antiseptic, e.g. 70 % alcohol.
3.	Mount a new gasket into the calibration chamber.
	<b>NOTICE:</b> Be sure to align the gasket correctly.



## Preparation and maintenance of the tc sensors

General information about membraning .....	4-10
Cleaning the sensor head of the E5480 sensor .....	4-11
Membraning the E5480 sensor.....	4-12
Membraning the E5280 and E5260 sensors.....	4-13
Cleaning, disinfection and storage of sensors .....	4-15

## General information about membraning

To obtain reliable measurements, remembrane the sensor every week.



**WARNING – *Risk of incorrect measurements (sensor failure)***




After remembraning the sensor, check that the membrane has been clicked completely and centrally onto the sensor housing. Otherwise it may cause the sensor to fail.

When a sensor has been remembraned, connect the sensor plug to the sensor socket at the rear of the monitor. Check that the system shows "Calibration required" and calibrate the sensor twice as described in chapter 5: *Calibration*.

## Cleaning the sensor head of the E5480 sensor

### Cleaning the sensor head







**NOTICE:** Once a month, the sensor head must be cleaned to remove old electrolyte before the sensor is remembraned.

Step		Action
1.		Slide the sensor into the end of the membraning tool.
2.		Pull the sensor upwards to remove the membrane.
3.		Clean the sensor head with cleaning paper and remembrane the sensor as described on next page.

## Membraning the E5480 sensor







### Membraning the sensor

To membrane the sensor, use a membraning kit (code no. 905-805).

Step	Action
1.	 <p>Apply two drops of electrolyte solution to the membraning tool.</p> <p><b>NOTICE:</b> Ensure that there are no air bubbles in the electrolyte solution. If air bubbles are present, wait a few seconds and check again.</p>
2.	 <p>Place the sensor (without the protection cap) in the sensor slot.</p>
3.	 <p>a) To remove the old membrane, grip the membraning tool firmly at both ends.</p>  <p>b) Pull in the direction of the arrows until only one arrow is visible in the sensor slot.</p>
4.	 <p>To click on the new membrane, pull forcefully in the direction of the arrows until the tool locks and no arrows are visible in the sensor slot.</p>
5.	 <p>Remove the sensor and wipe off the surplus electrolyte solution with cleaning paper.</p>



## Membraning the E5280 and E5260 sensors

**Membraning the sensor** To membrane the sensor, use a membraning kit (code no. 904-892).

Step		Action
1.		Remove the protection cap. Then remove the old O-rings by sliding the O-ring remover under the O-ring, just above the arrow on the sensor house.
2.		Turn the O-ring remover clockwise to release the two O-rings.
3.		Peel off the old membranes.
4.		Clean the sensor surface: Absorb the old electrolyte solution with the cleaning paper.
5.		Rub the sensor measuring surface carefully two or three times to remove the thin layer of silver that has precipitated on the sensor.
6.		Apply two drops of the $\text{tcpCO}_2/\text{tcpO}_2$ electrolyte solution on the surface of the sensor. <b>NOTICE:</b> Check that the electrolyte solution covers the entire surface without air bubbles.

*Continued on next page*

## Membraning the E5280 and E5260 sensors, *Continued*

Membraning the sensor ( <i>continued</i> )	Step	Action
	7.	 <ul style="list-style-type: none"><li>• Place the membrane unit on a hard and stable surface.</li><li>• Turn the sensor slowly so that the measuring surface faces downwards.</li><li>• Insert the sensor head into the top of the green membrane unit.</li></ul>
	8.	 <ul style="list-style-type: none"><li>• Press the sensor firmly into the unit until a click is heard.</li><li>• Remove the sensor from the unit and wipe off the surplus electrolyte solution with the cleaning paper.</li></ul>

## Cleaning, disinfection and storage of sensors

### Cleaning the sensors

Wipe the following parts gently with a soft cloth moistened with skin antiseptic, e.g. 70 % alcohol:

- the sensor head
- the cable

**NOTICE:** Constant use of hand lotion containing isopropanol/propylalcohol and alcohol prior to handling the sensor may damage the cable. To avoid transferring lotion to the cable, dry hands prior to handling the sensor.

### Disinfection



#### **WARNING – Risk of infection**

Due to the nature and seriousness of diseases such as the Human Immunodeficiency Virus (HIV) (the causative agent of AIDS) and Hepatitis B, always regard equipment and accessories that can come into contact with human or animal tissues or fluids (particularly blood) as contaminated and potentially hazardous to avoid the risk of infection.

Disinfection of the sensor is carried out with a membraned sensor and normally prior to calibration.

Disinfection of the sensor and the cable can be carried out by immersing the sensor and the cable in a water-based disinfection solution, such as a germicidal water-based solution (e.g. one based on n-alkyl dimethyl benzyl ammonium chloride with isopropanol).



#### **WARNING – Risk of incorrect measurements (sensor failure)**

Do not immerse the sensor plug in disinfection solution. Immersing the sensor plug in disinfection solution will cause the sensor to fail.



#### **WARNING – Risk of incorrect measurements (sensor failure)**

Do not heat sterilize as the sensor cannot tolerate temperatures exceeding 70 °C. Heat sterilization of the sensor will cause it to fail.

As the number of disinfection solutions is increasing and may vary from country to country and from hospital to hospital, it is not possible to come up with a complete list of water-based disinfection solutions that may be used.

However, Radiometer proposes that one of the following water-based disinfection solutions is used:

- MadaCide-FD (MADA Inc.)
- Control III (Maril Products Inc.)
- Hi-Tor Plus (Ecolab Inc.)
- 3.2 % aqueous glutaraldehyde solution (Surgikos - Johnson & Johnson Medical)

Solutions solely based on isopropanol/propylalcohol should be avoided as frequent use of isopropanol/propylalcohol may damage the electrode cable.

**NOTICE:** To establish the correct disinfection procedure for your hospital, it is important that you consult the hygiene committee at your hospital.

### Storage

You may keep the sensors stored in the calibration chambers. For longer periods of storage, keep the sensors mounted with a protective cap.

**NOTICE:** Add two drops of electrolyte solution to the cap.

## Maintenance of the SpO<sub>2</sub> sensors

### Cleaning the SpO<sub>2</sub> sensors

The SpO<sub>2</sub> sensors may be surface-cleaned with a solution such as 70 % isopropyl alcohol. If low-level disinfection is required, use a 1:10 bleach solution. Do not use undiluted bleach (5-5.25 % sodium hypochlorite) or any cleaning solution other than those recommended here because permanent damage to the sensor could occur.



**WARNING – Risk of incorrect measurements (sensor failure)**

Do not expose connector pins to cleaning solution as this may damage the SpO<sub>2</sub> sensor.



**WARNING – Risk of incorrect measurements (sensor failure)**

Do not sterilize the SpO<sub>2</sub> sensor by irradiation, steam or ethylene oxide as this may damage the sensor.

There are two recommended cleaning/disinfection methods:

- The wipe method (for all three SpO<sub>2</sub> sensors)
- The soak method (for Oxiband A/N and P/I sensors only)

### Wipe method procedure

Step	Action
1.	Saturate a clean, dry gauze pad with the cleaning solution. Wipe all surfaces of the sensor and cable with this gauze pad.
2.	Saturate another clean, dry gauze pad with sterile or distilled water. Wipe all surfaces of the sensor and cable with this gauze pad.
3.	Dry the sensor and cable by wiping all surfaces with a clean, dry gauze pad.


### Soak method procedure

For Oxiband A/N and P/I sensors only.



**WARNING – Risk of incorrect measurements (sensor failure)**

Do not immerse or wet the DS100A sensor as this may damage the sensor.

Step	Action
1.	Place the sensor in the cleaning solution, such that the sensor head(s) and desired length of cable are completely immersed.
	 <b>WARNING – Risk of incorrect measurements (sensor failure)</b> Do not immerse the connector end of an SpO <sub>2</sub> sensor cable as this may damage the sensor.
2.	Dislodge air bubbles by gently shaking the sensor and cable.
3.	Soak the sensor and the cable for 10 minutes.
4.	Remove from cleaning solution.
5.	Place the sensor and the cable in room-temperature sterile or distilled water for 10 minutes.
6.	Remove from the water.
7.	Dry the sensor and cable by wiping all surfaces with a clean, dry gauze pad.



## 5. Calibration

General information .....	5-2
Calibration.....	5-3
Checking zero current and sensitivity .....	5-4

## General information

To ensure accurate and safe performance of the sensor, it must be calibrated according to the recommendations below.

### Calibration material

To calibrate the sensor, use the CAL1 standard calibration gas mixture (7.5 % CO<sub>2</sub>, 20.9 % O<sub>2</sub>, balance N<sub>2</sub>).

### SmartCal

The SmartCal function makes sure that the monitor is always ready by calibrating when needed for a period of time that is defined by the operator (1-12 hours or forever).

### Calibration frequency

SmartCal	Explanation
ON	<p>If SmartCal is ON and the duration is set to forever or 1-12 hours, a calibration will be performed automatically every 4 hours (i.e. whenever the sensor is placed in the calibration chamber and within the selected SmartCal period).</p> <p><b>NOTICE:</b> If the duration is set to 1-12 hours, it will be necessary to press <b>Calibrate</b> to start a new SmartCal period after the first period ended.</p>
OFF	<p>If SmartCal is OFF, Radiometer recommends performing a manual calibration:</p> <ul style="list-style-type: none"> <li>• every time the sensor has been remembraned</li> <li>• prior to each monitoring period</li> <li>• when changing measuring sites</li> <li>• every 4 hours</li> </ul> <p><b>NOTICE:</b> If the sensor was not in use for an extended time period, it is recommended to select SmartCal ON and to allow at least 4 hours for stabilization of the sensor in the calibration chamber.</p>

### Recommendation

Check the barometer of the monitor against a known calibrated barometer periodically. See ranges in chapter 8: *Specifications and ordering information*.

## Calibration

**Calibrating the sensor manually** If the SmartCal function is OFF, the sensor must be calibrated manually using the procedure below.



**WARNING – Risk of incorrect measurements**

When performing a gas calibration, make sure the calibration gas mixture is set to 7.5 % for  $p\text{CO}_2$  and 20.9 % for  $p\text{O}_2$  as incorrect calibration values may cause incorrect measurements.

See Technical settings in chapter 3: *Menu structure and setup programs*.

Step	Action
1.	Make sure the sensor is placed in the calibration chamber.
2.	Press <b>Calibrate</b> .  A bar appears on the screen, showing the progress of the calibration.  Furthermore, the screen contains information on barometric pressure and, if there is 10 % or less gas left in the gas cylinder, a bar will show the gas level.  <b>NOTICES:</b> <ul style="list-style-type: none"> <li>• If SmartCal has been activated in the setup, pressing the <b>Calibrate</b> touch key will start the SmartCal period.</li> <li>• During calibration it is not possible to reactivate the <b>Calibrate</b> touch key, which will be grayed out.</li> </ul>
3.	When the calibration is complete, the Ready screen appears.



**NOTICES**

- If the sensor has been remembraned or if it has not been used for 24 hours or more, it must be calibrated twice. Either calibrate the sensor, leave it in the calibration chamber for 30 minutes and then calibrate it again; or enable the SmartCal function, press **Calibrate** and then leave the sensor in the calibration chamber for 30 minutes.
- If SmartCal is OFF and the sensor is not removed from the calibration chamber within 30 minutes after Ready, the heat to the sensor will be switched off and a new calibration will be required.

## Checking zero current and sensitivity

**Recommendations** The zero current of  $pO_2$  and the sensitivity of  $pCO_2$  should be checked if the sensor performance appears to deteriorate.

**Items required** The following items are required:

- CAL2 standard calibration gas (10 %  $CO_2$  with  $N_2$  as balance)
- Gas adapter for CAL2 gas

### Checking $pO_2$ zero current and $pCO_2$ sensitivity

Step	Action
1.	Before calibrating the sensor, set the metabolic correction factor to zero and the Severinghaus correction factor to "OFF" in the Technical setup.
2.	Calibrate the sensor.
3.	Attach the adapter to the CAL2 gas cylinder and place the sensor in the calibration chamber of the gas adapter.
4.	The $pCO_2$ reading is displayed on the screen. It should be within 73-79 mmHg (9.7-10.5 kPa) within 10 minutes after the sensor has been placed in the calibration chamber of the gas adapter.
5.	Read the $tcpO_2$ value on the screen: <ul style="list-style-type: none"> <li>• If under 5 mmHg (0.7 kPa), the sensor is in good condition.</li> <li>• If equal to or above 5 mmHg (0.7 kPa), the sensor must be remembraned and recalibrated. Then repeat steps 1-3 of this procedure.</li> </ul>
	<b>NOTICE:</b> If the $tcpO_2$ reading is still equal to or above 5 mmHg (0.7 kPa), the sensor is defective.
6.	After zero current and sensitivity have been checked, unscrew the CAL2 gas cylinder from the gas adapter.
7.	Remember to reset the metabolic correction factor to "7" and the Severinghaus correction factor to "ON" in the Technical setup.

**NOTICE:** Lack of sensitivity may be due to insufficient removal of used electrolyte solution during the membraning procedure. In such cases, it is recommended to remembrane the sensor according to the procedure described in chapter 4: *Installation and maintenance*.

## 6. Patient monitoring

General information .....	6-2
Application and removal of tc sensors .....	6-6
Application of SpO <sub>2</sub> sensors .....	6-8
Patient monitoring (In vivo monitoring) .....	6-11
Patient DMS .....	6-13
Results in Normal view .....	6-16
Results in Trend table view .....	6-18
Results in Trend curve view .....	6-19
Analog output .....	6-21
Continuous data output: standard .....	6-22
Continuous data output: VueLink .....	6-24
Continuous data output: MonLink .....	6-29
Data export: serial .....	6-30
Data export: USB .....	6-33
Alarms .....	6-34
How to print .....	6-38
Blood gas comparison .....	6-44
In vivo calibration .....	6-45

## General information

### Measuring sites *tcpCO<sub>2</sub>/tcpO<sub>2</sub>* measurements:

Clinical studies have shown the abdomen and chest to be the best measuring sites for both neonates and adults.

#### Saturation measurements:

- Nellcor Durasensor DS100A sensor:
  - Patient (> 40 kg): The preferred site is the index finger, or alternatively a smaller finger, but *not* the thumb.
- Nellcor Oxiband A/N sensor:
  - Adult (> 40 kg): The preferred site is around an index finger, with the cable positioned along the top of the finger. Alternatively, use a thumb or another finger, with the cable positioned along the palm; or around a great toe, with the cable positioned along the sole of the foot.
  - Neonate (< 3 kg): The preferred site is the foot, below the toes, with the cable positioned along the sole of the foot. Alternatively, place the sensor around the palm of a hand, below the fingers, with the cable positioned along the palm.
- Nellcor Oxiband P/I sensor:
  - Pediatric (15-40 kg): The preferred site is around an index finger, with the cable positioned along the top of the finger. Alternative sites are around the thumb or another finger, with the cable positioned along the top of the finger; or a great toe, with the cable positioned along the sole of the foot.
  - Infant (3-15 kg): The preferred site is around a great toe, with the cable positioned along the sole of the foot.

**NOTICE:** The monitor must only be connected to sensors for one patient at a time.

### *tcpCO<sub>2</sub>/tcpO<sub>2</sub>* sensor temperature

For neonates, a sensor temperature between 42 and 44 °C is recommended.

For adults, a sensor temperature between 43 and 45 °C is recommended.



#### **WARNING – Risk of burns**

Do not allow the *tcpCO<sub>2</sub>/tcpO<sub>2</sub>* sensor temperature to exceed 43 °C for neonates and 44 °C for adults if sensors are attached to skin for more than four hours as this may otherwise cause burns.

### General alerts



#### **WARNING – Risk of incorrect measurements**

Do not use the sensors during MRI scanning. Conducted current may cause burns. Also, the sensors may affect the MRI image, and the MRI unit may affect the accuracy of oximetry measurements.



#### **WARNING – Risk of strangulation**

As with all medical equipment, carefully route and affix patient cabling using the cable clip to reduce the possibility of patient entanglement or strangulation.



#### **WARNING – Risk of skin damage**

To avoid the risk of skin damage, make sure to set SmartHeat to OFF before applying the sensor to a neonate.

---

*Continued on next page*

## General information, *Continued*

### General alerts (*continued*)

**WARNING – Risk of incorrect measurements**

Remove the sensors from the patient immediately if the system or patient is exposed to a defibrillator, electrocautery or other high-frequency electrical signals, as these may affect the device and may cause injury to the patient.

**WARNING – Risk of fire**

Do not place the monitor in an enriched oxygen atmosphere or inside a hyperbaric chamber as it may cause a fire hazard.

**WARNING – Risk of personal injury**

Make sure to select the upper alarm limit for oxygen saturation carefully and in accord with accepted clinical standards. High oxygen levels may predispose a premature infant to develop retinopathy.

**WARNING – Risk of skin damage**

Long-term hyperthermia may blister skin. When producing local hyperemia by means of hyperthermia, a certain risk of applying temperatures harmful to the skin is always present, although the risk is limited due to the control system of the instrument. Always pay attention to the use of hyperthermia for special patients – e.g. patients in shock, patients with low blood pressure, and patients with vascular constrictions.

**WARNING – Risk of bruises**

When applying a fixation ring to a patient, make sure to place it so that the patient does not lie on top of it, as this may cause the fixation ring to leave bruises on the patient.

**WARNING – Risk of incorrect measurements**

Always keep protection caps on the tcpCO<sub>2</sub>/tcpO<sub>2</sub> sensors (except when applied to skin, placed in the calibration chamber or during handling and maintenance). Exposing sensor membranes to light (for example light from incubators) may cause elevated pCO<sub>2</sub> values.

**WARNING – Risk of incorrect measurements**

Do not use the monitor adjacent to or stacked with other equipment as these can cause electromagnetic interference and thereby result in incorrect measurements. If stacking or use adjacent to other equipment is necessary, the monitor should be observed to verify normal operation before used on patients. See the section *EMC approvals and compliance* in chapter 8.

**WARNING – Risk of incorrect measurements**

When installing, operating or servicing the monitor, special consideration should be given to the information regarding the electromagnetic precautions for this equipment given in the section *EMC approvals and compliance* in chapter 8. Otherwise the monitor may be affected by electromagnetic interference, causing incorrect measurements.

**WARNING – Risk of incorrect measurements**

The use of sensors, cables and accessories other than those specified may result in increased emission and/or decreased immunity and inaccurate readings of the monitor.

---

*Continued on next page*

## General information, *Continued*

### General alerts (continued)



**WARNING – Risk of incorrect measurements**

tcpCO<sub>2</sub>/tcpO<sub>2</sub> monitoring should not be used on patients in a compromised hemodynamic state as this may cause incorrect measurements.



**WARNING – Risk of incorrect measurements**

Using the SpO<sub>2</sub> sensors in the presence of high ambient light may result in inaccurate measurements. In such cases, cover the sensor site with an opaque material.



**WARNING – Risk of incorrect measurements**

Remove any intravascular dyes or externally applied coloring such as nail polish, dye or pigmented cream as they may lead to inaccurate pulse oximetry measurements.



**WARNING – Risk of incorrect measurements**

Pulse oximetry readings and pulse signal can be affected by certain ambient environmental conditions, sensor application errors and certain patient conditions. See the appropriate sections of the manual for specific safety information.



**WARNING – Risk of incorrect measurements**

Always select the measuring site carefully to avoid selecting a site with low perfusion or low signal quality, which can cause incorrect measurements.



**WARNING – Risk of patient not being monitored**

Do not use a NIBP cuff or other constricting devices on the same appendage as the sensor. A NIBP cuff will interrupt the patient's circulatory blood flow and result in no pulse found or loss of pulse.

### TCM40 performance considerations

Certain patient conditions can affect the measurements of the monitor and cause the loss of the pulse signal.

Inaccurate measurements can be caused by:

- prolonged patient movement
- venous pulsations
- intravascular dyes, such as indocyanine green or methylene blue
- defibrillation

**Dysfunctional hemoglobins** such as carboxyhemoglobin, methemoglobin, and sulfhemoglobin are unable to carry oxygen. SpO<sub>2</sub> readings may appear normal; however, a patient may be hypoxic because less hemoglobin is available to carry oxygen. Further assessment beyond pulse oximetry is recommended.

**Anemia** causes decreased arterial oxygen content. Although SpO<sub>2</sub> readings may appear normal, an anemic patient may be hypoxic. Correcting anemia can improve arterial oxygen content. The monitor may fail to provide an SpO<sub>2</sub> if hemoglobin levels fall below 5 g/dL.

**Saturation:** The monitor displays saturation levels between 1 and 100 %.

**Pulse rates:** The monitor displays pulse rates between 20 and 300 beats per minute. Detected pulse rates outside the range of 20 to 300 beats per minute are displayed as the closest value within the range.

*Continued on next page*



## General information, *Continued*

### SpO<sub>2</sub> sensor performance considerations

Inaccurate measurements may be caused by:

- incorrect application of the sensor
- placement of the sensor on an extremity with a blood pressure cuff, arterial catheter, or intravascular line
- ambient light
- prolonged patient movement

Loss-of-pulse signal can occur for the following reasons:

- The sensor is applied too tightly
- A blood pressure cuff is inflated on the same extremity as the one with the sensor attached
- There is arterial occlusion proximal to the sensor

Use only Radiometer-recommended sensors and sensor cables.

Select an appropriate sensor, apply it as directed and observe all warnings and cautions presented in the directions for use accompanying the sensor. Clean and remove any substances such as nail polish from the application site. Periodically check to ensure that the sensor remains properly positioned on the patient.



#### **WARNING – Risk of skin damage**

Inspect the sensor site as directed in the sensor directions for use. Incorrect application or inappropriate duration of use of an SpO<sub>2</sub> sensor can cause skin damage.



#### **WARNING – Risk of incorrect measurements**

Make sure the sensor is applied correctly. Incorrect application of the SpO<sub>2</sub> sensor can cause incorrect measurements.

High ambient light sources such as surgical lights (especially those with a xenon light source), bilirubin lamps, fluorescent lights, infrared heating lamps, and direct sunlight can interfere with the performance of an SpO<sub>2</sub> sensor. To prevent interference from ambient light, ensure that the sensor is properly applied and cover the sensor site with opaque material.

If patient movement presents a problem, try one or more of the following remedies to correct the problem:

- Verify that the sensor is properly and securely applied
- Move the sensor to a less active site
- Use an adhesive sensor that tolerates some patient motion
- Use a new sensor with fresh adhesive backing

## Application and removal of tc sensors

### Required materials for application of sensors

The following items, included in the fixation kits, are required to apply the tc sensors to a patient:

- Fixation ring



#### **WARNING – Risk of infection and inaccurate results**




Reuse of single-use devices may lead to infection of patients and inaccurate results.

- Contact liquid

### Prior to sensor application

Step	Action
1.	Calibrate the sensor as described in chapter 5: <i>Calibration</i> .
2.	Clean the selected measuring site with alcohol or other skin-preparation solution.
3.	Dry the site well with a gauze pad.




### Application of E5480 sensor

Step	Action
1.	Take a fixation ring and remove the protective film.
2.	 <p>Apply the fixation ring to clean and dry skin by pressing the center of the fixation ring onto the measuring site with a finger and then running a finger around the rim circumference to ensure a good seal.</p>
3.	 <p>Apply 2 drops of contact liquid in the center of the ring.</p>
4.	 <p>Place the sensor (without the protection cap) in the fixation ring with the sensor cord pointing in the opposite direction of the tab on the ring. Then turn the sensor a quarter of a turn clockwise to fasten it in the fixation ring.</p>

*Continued on next page*

## Application and removal of tc sensors, *Continued*

### Application of E5280 and E5260 sensors

Step	Action	
1.	Take a fixation ring and remove the protective film.	
2.		Apply the fixation ring to the measuring site by pressing the center of the fixation ring onto the measuring site with a finger and then running a finger around the rim circumference.  <b>NOTICE:</b> Press firmly to prevent leaks.
3.		Fill the center of the fixation ring with 3-5 drops of the contact liquid.
4.		Remove the sensor from the calibration chamber, align the arrow on the sensor with one of the marks on the fixation ring and turn the sensor a quarter of a turn clockwise to fasten it in the fixation ring.

### Removal of sensor

Step	Action
1.	Remove the sensor from the fixation ring by turning it anticlockwise.
2.	Remove the fixation ring by lifting the tab.
3.	Clean the sensor surface carefully with an alcohol swab.
4.	Then place the sensor in the calibration chamber.

## Application of SpO<sub>2</sub> sensors

### Required materials for SpO<sub>2</sub> monitoring

Saturation monitoring is only available on the TCM40 monitor.

The following items are required to apply an SpO<sub>2</sub> sensor to a patient:

- Nellcor SpO<sub>2</sub> sensor
- Adhesive wrap model ADH-A/N or FOAM-A/N (for Oxiband A/N and Oxiband P/I sensors)
- Nellcor DOC-10 pulse oximetry cable



#### **WARNING – Risk of incorrect measurements**

Use only Nellcor oximetry brand SpO<sub>2</sub> sensors and sensor cables. And before use, carefully read the sensor directions for use, including all warnings, cautions and instructions. Incorrect use can cause incorrect measurements.



#### **WARNING – Risk of incorrect measurements**

Use only the DOC-10 pulse oximetry cable and do not increase the length of the sensor by connecting an extra sensor cable. Use of another sensor cable or more than one cable will have an adverse effect on performance.



#### **WARNING – Risk of skin damage**

Use only Nellcor oximetry brand wraps designed for use with the SpO<sub>2</sub> sensor. Do not use tape. Use of additional tape or various other types of wraps can cause skin damage.

### NOTICES:





- Radiometer provides a starter kit containing sensor (Nellcor DS100A, Nellcor Oxiband A/N or Nellcor Oxiband P/I), wraps and cable. New sensors (or other Nellcor sensors than the one in the starter kit), wraps and cables must be ordered from your local Nellcor agent.
- When selecting an SpO<sub>2</sub> sensor, consider the patient's weight and activity level, the adequacy of perfusion and the available sensor sites, the need for sterility and the anticipated duration of monitoring.
- When selecting a sensor site, priority should be given to an extremity free of an arterial catheter, blood pressure cuff or intravascular infusion line.
- Before applying the sensor, clean and remove any substances, such as nail polish, from the application site.

---




*Continued on next page*

## Application of SpO<sub>2</sub> sensors, *Continued*

### Application of DS100A sensor

Step		Action
1.		Place the patient's index finger over the sensor window of the DS100A sensor with the fingertip against the stop.
2.		If the fingernail is long, the nail tip will extend over the finger stop.
3.		Spread open the rear tabs of the sensor to provide even force over the length of the pads. Check the position of the sensor. If an index finger cannot be positioned correctly, or is not available, a smaller finger can be used, or use another sensor.
4.		<b>NOTICE:</b> Do not use the DS100A on a thumb or toe or across a child's hand or foot.  The sensor should be oriented in such a way that the cable is positioned along the top of the hand.

### Application of Oxiband A/N and P/I sensors

Step		Action
1.		Place the wrap on a flat surface, adhesive (sticky) side up, with the tab on the left.
2.		Remove the small piece of paper backing and enough of the large piece of backing to expose both alignment holes.
3.		Position the sensor on the wrap so the alignment bumps protrude through the holes. Press down firmly in the center of the sensor to ensure that the wrap adheres tightly.
4.		Remove the remainder of the paper backing from the wrap. The sensor and wrap are now ready to be applied to the patient. Select an appropriate site as described previously in this chapter under <i>General information</i> .

*Continued on next page*

## Application of SpO<sub>2</sub> sensors, *Continued*

Application of Oxiband A/N and P/I sensors ( <i>continued</i> )	Step	Action
	5.	Position the sensor so that the notches are centered on the side of the measuring site.  The cable end must be positioned along the appropriate surface, as described previously in this chapter under <i>General information</i> .
	6.	Wrap the sensor around the measuring site so the optical components (and alignment bumps) oppose each other across the site. Press the sensor to ensure that it adheres snugly.
	7.	Wrap the remaining portion of the sensor wrap around the site, loosely enough to ensure good circulation.

**Removal of wraps** The disposable wraps used with the Oxiband A/N and P/I sensors should be removed by peeling them away from the sensor.



**WARNING – Risk of sensor damage**

Do not use excessive force when removing the disposable wrap from the sensor, as it may damage the sensor.

## Patient monitoring (In vivo monitoring)

Patient monitoring	Step	Action
	1.	Apply the sensor(s) as described earlier in this chapter and wait for a stable reading.
		<b>NOTICES:</b>
		<ul style="list-style-type: none"> <li>• The physiological stabilization time of a patient is 10-20 minutes for the tcpO<sub>2</sub> reading. Longer time may indicate incorrect sensor application or a poorly selected measuring site.</li> <li>• The physiological stabilization time of a patient is 5-10 minutes for the tcpCO<sub>2</sub> reading. Longer time may indicate incorrect sensor application or a poorly selected measuring site.</li> <li>• Excessive motion may compromise performance. In such cases, try to keep the patient still or change the sensor site to one with less motion.</li> <li>• If SmartHeat is set to ON, it adds +1 °C (max. temp. 45 °C) to the set sensor temperature for 5 minutes after the sensor has been removed from the calibration chamber.</li> <li>• If the SpO<sub>2</sub> sensor does not track the pulse reliably, it may be incorrectly positioned or the sensor site may be too thick, thin or deeply pigmented, or otherwise deeply colored (for example as a result of externally applied coloring such as nail polish, dye or pigmented cream) to permit appropriate light transmission. If any of these situations occurs, reposition the sensor or choose an alternate sensor for use on a different site.</li> <li>• The visual indication of the plethysmograph curve is not proportional to the pulse volume.</li> </ul>
	2.	Press <b>Site time</b> to reset the site timer to the value selected in Parameter setup, if required. The site timer applies to tcpCO <sub>2</sub> /tcpO <sub>2</sub> only.
		<b>NOTICES:</b>
		<ul style="list-style-type: none"> <li>• The site timer will count down to zero at 1-minute intervals, and when it reaches zero, the message "Site time end" will be displayed on the screen.</li> <li>• If Site time heat is set to OFF in the setup, the sensor heat is switched off when the site timer reaches zero and the monitor stops monitoring; if set to ON, the heat continues.</li> </ul>
	3.	To mark an event, press <b>Event</b> . The text "Set event mark number x?" will be displayed.
		Press <b>OK</b> to mark the event, or press <b>Cancel</b> to return to Normal view without marking the event.
		<b>NOTICE:</b> It is possible to mark up to 99 events.

*Continued on next page*

## Patient monitoring (In vivo monitoring), *Continued*

### Patient monitoring (*continued*)

Step	Action
4.	To link the session (i.e. the measurement in progress) to a specific patient ID, press <b>ID</b> in the top right corner of the view screens and make the necessary changes in the Patient DMS. For more information, see <i>Patient DMS</i> later in this chapter.
5.	View the readings and adjust the settings, if necessary. See chapter 3: <i>Menu structure and setup programs</i> .



#### **WARNING – Risk of skin damage**

Sensors must be moved to a new site at least every four hours. Because individual skin condition affects the ability of the skin to tolerate sensor placement, it may be necessary to change the sensor site more frequently with some patients. If skin integrity changes, move the sensor to another site.



#### **WARNING – Risk of incorrect measurements**

Always keep protection caps on the  $\text{tcpCO}_2/\text{tcpO}_2$  sensors (except when applied to skin, placed in the calibration chamber or during handling and maintenance). Exposing sensor membranes to light (for example light from incubators) may cause elevated  $p\text{CO}_2$  values.

**NOTICE:** After removing the sensor from the patient, make sure to wipe it clean from contact liquid before placing it in the calibration chamber.



## Patient DMS

The Patient DMS (Data Management System) manages all session/patient ID data, and it helps avoiding the risk of patient data mix-up.

From the Patient DMS screen it is possible to:

- change the automatically generated session number to a unique patient ID
- view data from one or more sessions with the same ID in the Trend table or Trend curve view
- print one or more sessions
- export one or more sessions
- delete sessions
- see detailed information about a session


### Session

A session is a collection of data starting when the sensor is removed from the calibration chamber and ending when the sensor is placed back in the calibration chamber.




Each session gets a unique number, which can be linked to a specific patient ID at any time. It is thereby possible to link several sessions to the same patient ID.

### Accessing Patient DMS

Press **ID** in the top right corner of the view screens to enter the Patient DMS.





### Touch keys

Touch key	Function
 and 	Scroll between the patient IDs/session numbers in the list.
	Adds a check mark to the highlighted patient ID/session no. <b>NOTICE:</b> Only patient IDs/session numbers with a check mark can be viewed, deleted, printed and exported.

*Continued on next page*

## Patient DMS, *Continued*

### Touch keys (*continued*)

Touch key	Function
	Shows detailed information about the highlighted patient ID/session number. 
<b>Change ID</b>	Enables changing the highlighted patient ID/session number with the on-screen alphanumeric keyboard (see next page). Maximum 64 characters. <b>NOTICE:</b> When changing a patient ID/session number, make sure to change it to something unique for each patient. If the same ID is used for two patients, the DMS will not be able to tell them apart.
<b>Delete</b>	Deletes the patient IDs/session numbers with a check mark. <b>NOTICES:</b> <ul style="list-style-type: none"> <li>• A dialog warns users that the action cannot be undone.</li> <li>• The session in progress cannot be deleted.</li> </ul>
<b>Delete all</b>	Deletes all patient IDs/session numbers in the list. <b>NOTICES:</b> <ul style="list-style-type: none"> <li>• A dialog warns users that the action cannot be undone.</li> <li>• The session in progress cannot be deleted.</li> </ul>
<b>Print</b>	Prints a report (one or more report types) for the patient IDs/session numbers with a check mark.
<b>Export</b>	Exports the patient information and measuring data for the patient IDs/session numbers with a check mark to an external PC or a memory stick (see the sections <i>Data export: serial</i> and <i>Data export: USB</i> later in this chapter). <b>NOTICE:</b> The <b>Export</b> touch key is only active if "Data export" has been selected in Technical settings.
<b>OK</b>	Loads sessions with a check mark into the Trend table and Trend curve views. <b>NOTICES:</b> <ul style="list-style-type: none"> <li>• Only sessions with the same patient ID can be loaded into the same view.</li> <li>• If changes are made to the sessions that have been loaded into a view, all the sessions will have to be loaded into the view again.</li> </ul>
<b>Cancel</b>	Returns to the view screen.

*Continued on next page*

## Patient DMS, *Continued*

### On-screen alphanumeric keyboard


The following screen shows the basic keyboard layout:



Touch Key	Function
	Escape. To return to the Patient DMS screen without making changes.
	To switch between general and language-specific character sets. Placed in the upper right corner of the screen. <b>NOTICE:</b> Not applicable in English.
	Backspace. To delete one character at a time from right to left.
	Shift. To shift between different keyboards in the same character set.
or	To scroll left/right in the text edit field.
	Space. To add a space.
	Enter. To accept the changes made in the text edit field and return to the Patient DMS screen.

## Results in Normal view

### Example of Normal view screen during monitoring




### NOTICES:

- The SpO<sub>2</sub> plethysmograph displays data at a speed of 25 mm/s in Normal view.
- The visual indication of the plethysmograph curve is not proportional to the pulse volume.
- The pCO<sub>2</sub>, pO<sub>2</sub>, SpO<sub>2</sub> and Pulse parameter values are updated every 2 seconds.
- SpO<sub>2</sub> and pulse are available on the TCM40 monitor only.
- On the TCM4 monitor, the lower part of the curve display will either show the power curve or be empty, and the Pulse/SpO<sub>2</sub> display will always be empty.


### Curve display

In Normal view, it is possible to change the curve ranges setup and the time span setup during measurement by pressing the respective areas in the curve display.



### pCO<sub>2</sub>/pO<sub>2</sub> display

In Normal view, it is possible to change the parameter setup during measurement by pressing the respective areas in the pCO<sub>2</sub>/pO<sub>2</sub> display.



*Continued on next page*

## Results in Normal view, *continued*

### *pCO<sub>2</sub>/pO<sub>2</sub>* display (*continued*)

The *pCO<sub>2</sub>/pO<sub>2</sub>* display also contains the following information:

Part	Information
Corr.	Metabolic or Severinghaus corrections are enabled.
Power	Current power delivered to the sensor.
Temp.	Current sensor temperature. <b>NOTICE:</b> If the sensor temperature is flashing, the sensor temperature differs from the set sensor temperature with $\pm 0.6$ °C.
SmartHeat	SmartHeat has been enabled in the setup.
In vivo calibration active	Current <i>pCO<sub>2</sub>/pO<sub>2</sub></i> values are in vivo calibrated.

### Pulse/SpO<sub>2</sub> display

In Normal view it is possible to change the Pulse/SpO<sub>2</sub> setup during measurement by pressing the Pulse/SpO<sub>2</sub> display.



The Pulse/SpO<sub>2</sub> display contains the following information:

Part	Information
Pulse	Pulse rate in beats per minute (bpm).
SpO <sub>2</sub>	Saturation level of oxygenated hemoglobin in %.
HI OFF	SpO <sub>2</sub> alarm high limit is disabled.

### Markings in Normal view

The following markings can be seen in Normal view:

Marking	Information
●	Blood gas values have been entered.
*	An in vivo calibration has been performed.
1-99	An event mark number has been added.





## Results in Trend table view

### Viewing results in Trend table view




#### Step Action

1. Press **Setup** → **View** → **Trend table** → **OK**.

The following is an example of a trend table on a TCM40 monitor:

14 04 22/01/2009	Trend table view					ID 22012009-0014
		pCO <sub>2</sub> mmHg	pO <sub>2</sub> mmHg	Power mW	SpO <sub>2</sub> %	Pulse bpm
Events	Current:	36	156	260	97	64
	14:00:00 22/01	36	157	260	---	---
	13:45:00 22/01	36	157	258	---	---
	13:30:00 22/01	36	157	262	---	---
	13:15:00 22/01	36	157	262	---	---
	13:00:00 22/01	36	157	263	---	---
	12:45:00 22/01	36	157	259	---	---
	12:30:00 22/01	36	156	258	---	---
	12:15:00 22/01	36	157	261	---	---
	12:00:00 22/01	36	157	263	---	---
	11:45:00 22/01	36	157	266	---	---
	11:30:00 22/01	36	157	266	---	---
	11:15:00 22/01	36	157	268	---	---
<div>     <span>Setup</span> <span>Alarm silence</span> </div>						

**NOTICE:** SpO<sub>2</sub> and pulse are available on the TCM40 monitor only.

2. Press  or  to move up or down one line in the list of results, or press  to see the most recent results.
3. To print the results, press **Print**. See also *How to print* later in this chapter.

### Markings in Trend table view

The following markings can be seen in the Trend table:

Marking	Indication
●	Blood gas values have been entered.
*	An in vivo calibration has been performed.
-	No value is available.
1-99	An event mark number has been added.


## Results in Trend curve view

### Viewing results in trend curve view


#### Step Action

1. Press **Setup** → **View** → **Trend curve** → **OK**.

The following is an example of a trend curve on a TCM40 monitor:



2. Press ◀ or ▶ to see more results to the left or the right.
3. Press **Cursor** to add a cursor to the screen. This will also change the function of the **Left** and **Right** arrow keys: they now move the cursor instead of time. If the cursor is moved all the way to one side, the time axis moves, which makes it possible to go back and forth in time.




The cursor values ( $p\text{CO}_2$ ,  $p\text{O}_2$ , Power,  $\text{SpO}_2$ , Pulse, Time and Date) are shown in the lower right part of the screen, whereas current values are shown in the upper right part of the screen.

Press **Cursor** again to remove the cursor from the screen.

*Continued on next page*

## Results in Trend curve view, *Continued*

### Viewing results in trend curve view (*continued*)

Step	Action
4.	To change the high or low values of the parameter ranges, press  .
5.	To print the results, press <b><i>Print</i></b> . See also <i>How to print</i> later in this chapter.

### Markings in Trend curve view

The following markings can be seen in the Trend curve view:

Marking	Indication
●	Blood gas values have been entered.
*	An in vivo calibration has been performed.
-	No value is available.
1-99	An event mark number has been added.



## Analog output

**Introduction** The monitor is equipped with an output for direct connection to an external chart recorder or polysomnograph.

**Setting up analog output** Follow the steps below to connect the TCM4/40 monitor to an external chart recorder or polysomnograph:

Step	Action																											
1.	Connect the TCM4xx ETX Analog Adapter (code no. 636-650) to the analog port.  <b>NOTICE:</b> Do not use other adapters than 636-650.																											
2.	Connect the wires as follows:																											
	<table border="1"> <thead> <tr> <th>Analog output</th> <th>Wire color</th> <th>Range</th> </tr> </thead> <tbody> <tr> <td>Ground</td> <td>Blue</td> <td>N/A</td> </tr> <tr> <td>Alarm</td> <td>Pink</td> <td>Activated: 1000 mV, <math>\pm 10</math> mV Not activated: 0 mV, <math>\pm 10</math> mV</td> </tr> <tr> <td>Heat</td> <td>Brown</td> <td>1 mV/mW @ 10-1000 mW, <math>\pm 10</math> mV (<math>\pm 10</math> mW) &lt; 10 mW <math>\approx</math> 0 mV</td> </tr> <tr> <td>Temp</td> <td>Green</td> <td>20 mV/<math>^{\circ}</math>C @ 10-50 <math>^{\circ}</math>C, <math>\pm 10</math> mV (<math>\pm 0.50</math> <math>^{\circ}</math>C)</td> </tr> <tr> <td>tcpCO<sub>2</sub></td> <td>Grey</td> <td>10 mV/mmHg @ 0-100 mmHg, <math>\pm 10</math> mV (<math>\pm 1</math> mmHg) 5 mV/mmHg @ 0-200 mmHg, <math>\pm 10</math> mV (<math>\pm 2</math> mmHg)</td> </tr> <tr> <td>tcpO<sub>2</sub></td> <td>Yellow</td> <td>5 mV/mmHg @ 0-200 mmHg, <math>\pm 10</math> mV (<math>\pm 2</math> mmHg) 1 mV/mmHg @ 0-800 mmHg, <math>\pm 10</math> mV (<math>\pm 10</math> mmHg)</td> </tr> <tr> <td>SpO<sub>2</sub></td> <td>White</td> <td>10 mV/% SpO<sub>2</sub> @ 0-100 % (<math>\pm 1</math> % SpO<sub>2</sub>)</td> </tr> <tr> <td>Pulse</td> <td>Red</td> <td>4 mV/bpm @ 20-250 bpm <math>\pm 10</math> mV (<math>\approx \pm 3</math> bpm)</td> </tr> </tbody> </table>	Analog output	Wire color	Range	Ground	Blue	N/A	Alarm	Pink	Activated: 1000 mV, $\pm 10$ mV Not activated: 0 mV, $\pm 10$ mV	Heat	Brown	1 mV/mW @ 10-1000 mW, $\pm 10$ mV ( $\pm 10$ mW) < 10 mW $\approx$ 0 mV	Temp	Green	20 mV/ $^{\circ}$ C @ 10-50 $^{\circ}$ C, $\pm 10$ mV ( $\pm 0.50$ $^{\circ}$ C)	tcpCO <sub>2</sub>	Grey	10 mV/mmHg @ 0-100 mmHg, $\pm 10$ mV ( $\pm 1$ mmHg) 5 mV/mmHg @ 0-200 mmHg, $\pm 10$ mV ( $\pm 2$ mmHg)	tcpO <sub>2</sub>	Yellow	5 mV/mmHg @ 0-200 mmHg, $\pm 10$ mV ( $\pm 2$ mmHg) 1 mV/mmHg @ 0-800 mmHg, $\pm 10$ mV ( $\pm 10$ mmHg)	SpO <sub>2</sub>	White	10 mV/% SpO <sub>2</sub> @ 0-100 % ( $\pm 1$ % SpO <sub>2</sub> )	Pulse	Red	4 mV/bpm @ 20-250 bpm $\pm 10$ mV ( $\approx \pm 3$ bpm)
Analog output	Wire color	Range																										
Ground	Blue	N/A																										
Alarm	Pink	Activated: 1000 mV, $\pm 10$ mV Not activated: 0 mV, $\pm 10$ mV																										
Heat	Brown	1 mV/mW @ 10-1000 mW, $\pm 10$ mV ( $\pm 10$ mW) < 10 mW $\approx$ 0 mV																										
Temp	Green	20 mV/ $^{\circ}$ C @ 10-50 $^{\circ}$ C, $\pm 10$ mV ( $\pm 0.50$ $^{\circ}$ C)																										
tcpCO <sub>2</sub>	Grey	10 mV/mmHg @ 0-100 mmHg, $\pm 10$ mV ( $\pm 1$ mmHg) 5 mV/mmHg @ 0-200 mmHg, $\pm 10$ mV ( $\pm 2$ mmHg)																										
tcpO <sub>2</sub>	Yellow	5 mV/mmHg @ 0-200 mmHg, $\pm 10$ mV ( $\pm 2$ mmHg) 1 mV/mmHg @ 0-800 mmHg, $\pm 10$ mV ( $\pm 10$ mmHg)																										
SpO <sub>2</sub>	White	10 mV/% SpO <sub>2</sub> @ 0-100 % ( $\pm 1$ % SpO <sub>2</sub> )																										
Pulse	Red	4 mV/bpm @ 20-250 bpm $\pm 10$ mV ( $\approx \pm 3$ bpm)																										
3.	On the monitor, press <b>Setup</b> → <b>Technical</b> . Enter the password and press <b>Enter</b> and <b>Tech. settings</b> .																											
4.	Use the <b>Down</b> arrow to select the pCO <sub>2</sub> analog range and choose the 0-100 or 0-200 range.																											
5.	Use the <b>Down</b> arrow to select the pO <sub>2</sub> analog range and choose the 0-200 or 0-800 range.																											
6.	Press <b>OK</b> to accept the changes and return to the main screen.  <b>NOTICE:</b> To test the analog output, calibrate the monitor. When the calibration is complete, the analog output will correspond to the tcpCO <sub>2</sub> and tcpO <sub>2</sub> readouts on the display.																											

1. Connect the TCM4xx ETX Analog Adapter (code no. 636-650) to the analog port.

**NOTICE:** Do not use other adapters than 636-650.

2. Connect the wires as follows:

Analog output	Wire color	Range
Ground	Blue	N/A
Alarm	Pink	Activated: 1000 mV, $\pm 10$ mV Not activated: 0 mV, $\pm 10$ mV
Heat	Brown	1 mV/mW @ 10-1000 mW, $\pm 10$ mV ( $\pm 10$ mW) < 10 mW $\approx$ 0 mV
Temp	Green	20 mV/ $^{\circ}$ C @ 10-50 $^{\circ}$ C, $\pm 10$ mV ( $\pm 0.50$ $^{\circ}$ C)
tcpCO <sub>2</sub>	Grey	10 mV/mmHg @ 0-100 mmHg, $\pm 10$ mV ( $\pm 1$ mmHg) 5 mV/mmHg @ 0-200 mmHg, $\pm 10$ mV ( $\pm 2$ mmHg)
tcpO <sub>2</sub>	Yellow	5 mV/mmHg @ 0-200 mmHg, $\pm 10$ mV ( $\pm 2$ mmHg) 1 mV/mmHg @ 0-800 mmHg, $\pm 10$ mV ( $\pm 10$ mmHg)
SpO <sub>2</sub>	White	10 mV/% SpO <sub>2</sub> @ 0-100 % ( $\pm 1$ % SpO <sub>2</sub> )
Pulse	Red	4 mV/bpm @ 20-250 bpm $\pm 10$ mV ( $\approx \pm 3$ bpm)

3. On the monitor, press **Setup** → **Technical**. Enter the password and press **Enter** and **Tech. settings**.
4. Use the **Down** arrow to select the pCO<sub>2</sub> analog range and choose the 0-100 or 0-200 range.
5. Use the **Down** arrow to select the pO<sub>2</sub> analog range and choose the 0-200 or 0-800 range.
6. Press **OK** to accept the changes and return to the main screen.

**NOTICE:** To test the analog output, calibrate the monitor. When the calibration is complete, the analog output will correspond to the tcpCO<sub>2</sub> and tcpO<sub>2</sub> readouts on the display.

## Continuous data output: standard

**Introduction** The monitor is equipped with a continuous data output for direct connection to an external PC.

### Connecting to external PC

- | Step | Action   |
|------|--|
| 1.   | Connect the TCM4xx ETX Serial Adapter (code no. 636-649) to the serial port on the monitor and to the PC.<br><br>The RS232 output from the monitor is transmitted continuously every 2 seconds in ASCII code in the following format:<br>[Timestamp]; [O <sub>2</sub> ]; [CO <sub>2</sub> ]; [Heater power]; [Temperature]; [SpO <sub>2</sub> ]; [Pulse] |
| 2.   | On the monitor, press <b>Setup</b> → <b>Technical</b> , enter the password and press <b>Enter</b> and <b>Tech. settings</b> .  |
| 3.   | Use the <b>Up</b> and <b>Down</b> arrows to select "Continuous data output" and choose the option "Standard".  |
| 4.   | Press <b>OK</b> to accept the changes and return to the main screen.   |
| 5.   | On the PC, select <b>Start</b> → <b>Programs</b> → <b>Accessories</b> → <b>Communications</b> → <b>HyperTerminal</b> and open the EXE file.  |
| 6.   | Type in a name for the connection, e.g. TcData, and then click <b>OK</b> .   |
| 7.   |  |




Connect using COM1 and click **OK**.

*Continued on next page*

## Continuous data output: standard, *Continued*


Connecting to  
external PC  
(*continued*)

- | Step | Action   |
|------|--|
| 8.   | Select the following port settings and click <b>OK</b> . |



9. Collect the data in HyperTerminal by selecting **Transfer** → **Capture text**. Then name the file, e.g. CAPTURE.txt, and click **Start**.
10. Stop the collection of data to the file in HyperTerminal by selecting **Transfer** → **Capture text** → **Stop**.

The following shows a screen from HyperTerminal:



### NOTICES:

- Data is easily imported into other programs, e.g. Microsoft Excel (see the procedure under *Data export: serial*).
- It is possible to have analog and serial output at the same time. The analog wires are short-circuit proofed.

## Continuous data output: VueLink

**Introduction** With a VueLink interface module type B, parameter values and attached alarm limit values can be exported from a TCM monitor to a Philips Patient Monitoring System (PPMS). The presentation of the data corresponds to Normal view.

**NOTICE:** Only the parameter values are displayed on the PPMS. The alarm limit values can only be used by a data management system connected to the PPMS.

**Required items**

- TCM4/40 monitors (software version 2.10 or newer)
- PPMS monitor (connected to one or more module racks)
- VueLink interface module type B (M1032A #A05)
- Connecting cable (Philips code no. M1032-61654) and TCM4xx ETX VueLink Adapter (Radiometer code no. 636-651) for connecting the VueLink interface module to the TCM monitor


**NOTICE:** In the following, the term PPMS includes the VueLink interface module.

**Possible PPMS monitors** The following patient monitors are supported:

- CMS (software version C or newer)
- IntelliVue MP40/50/60/70/90 (all software versions)
- Agilent V24/V26 (all software versions)

**NOTICE:** Older HP systems are often called CMS instead of PPMS.

**Connecting to PPMS monitor**



Step	Action
1.	Insert a VueLink interface module in the PPMS monitor module rack.
2.	Connect the VueLink interface module to the TCM monitor with the connecting cable and the VueLink interface cable.
3.	Turn on both monitors.
4.	On the TCM monitor, press <b>Setup</b> → <b>Technical</b> (enter password) → <b>Tech. settings</b> , set "Continuous data output" to "VueLink" and press <b>Apply</b> or <b>OK</b> .

**NOTICE:** Changing settings for Unit will cause a temporary disconnection of communication of up to 60 seconds. The connection will be reestablished automatically.

*Continued on next page*


## Continuous data output: VueLink, *Continued*

### Verification of connection between TCM and PPMS monitors

- | Step | Action   |
|------|--|
| 1.   | <p>On the VueLink interface module, press the <b>VueLink</b> key.</p> <p>The "VueLink-B 1" tab appears on the PPMS monitor screen.</p> <p><b>NOTICES:</b></p> <ul style="list-style-type: none"> <li>• If the tab name is "AUXILIARY PLUS 1", connection to the module has already been established.</li> <li>• It is possible to insert several VueLink interface modules in the PPMS monitor, and they will then be called "VueLink-B 2", "VueLink-B 3", etc.</li> </ul> |
| 2.   | <p>Wait for the text "... Connect device or switch to new device" to disappear.</p> <p><b>NOTICE:</b> If the text has not disappeared after 30 seconds, check the cable connection. If the cable connection is okay, see chapter 7: <i>Troubleshooting</i> for further operator actions.</p>   |
| 3.   | <p>Select the <b>Setup VueLink</b> softkey.</p> <p>The "AUXILIARY PLUS 1" tab containing the parameter values from the TCM monitor appears.</p>  |

### Example of task window on PPMS monitor

The task window "AUXILIARY PLUS 1" is displayed automatically when connection between the TCM and the PPMS monitors has been established.



*Continued on next page*

## Continuous data output: VueLink, *Continued*

### Example of task window on PPMS monitor (continued)

#### NOTICES:

- Transmission of data from the TCM monitor to the PPMS monitor may be delayed with up to four seconds.
- Printout of data from the TCM monitor may cut off the connection between the TCM and the PPMS monitors while printing. The connection will be re-established automatically.
- For examples of other PPMS monitors, please refer to the user instructions for the relevant PPMS monitor.

The task window contains the following information:

Part	Shows
Device type	The type of connected module, e.g. AUXILIARY PLUS 1 (for the VueLink-B 1 module), AUXILIARY PLUS 2 (for the VueLink-B 2 module), etc.
Device name	The name of the connected TCM monitor, e.g. TCM4.
Alarm status indication	Whether alarms on the TCM monitor have been activated.
Parameter	The parameters exported from the TCM monitor.
Asterisk	That data for the following parameter is accessible on the main screen and in external databases.
Unit	The parameter units. <b>NOTICE:</b> Units are only shown in the task window.
Time	At which time the PPMS monitor started receiving data from the TCM monitor.

### Main screen on PPMS monitor

The main screen is displayed by pressing the **Main Screen** key on the PPMS monitor.

The parameters are displayed at the bottom right-hand corner of the PPMS monitor screen:

Parameter on PPMS (TCM monitor)	Color on PPMS (TCM monitor)	Unit on PPMS and TCM monitors
tcpCO <sub>2</sub> (tcpCO <sub>2</sub> )	Green (Blue)	mmHg/kPa
tcpO <sub>2</sub> (tcpO <sub>2</sub> )	Green (Red)	mmHg/kPa
SpO <sub>2</sub> (SpO <sub>2</sub> )	Green (Magenta)	%
PULSE (Pulse)	Green (Brown)	bmp

**NOTICE:** It is not possible to change the colors on the PPMS monitor.

*Continued on next page*

## Continuous data output: VueLink, *Continued*

### Configuration of PPMS monitor

Some features, keys or softkeys may vary from one monitor system to another, but the principles are basically the same. If you encounter differences between your system and what we describe in this manual, please see the technical documentation for the particular PPMS monitor.

The following configuration procedure is based on an Agilent V24C monitor.

Follow these steps to define how you want data presented on your PPMS monitor:

Step	Action
1.	Select the following softkeys: <b>Module Setup</b> → <b>TCM</b> (should have appeared when the TCM monitor was connected) → <b>Select Data</b> .
2.	Mark the module number ( <b>Num</b> ) you wish to change with the up/down arrow keys.
3.	Select the <b>Select Signal</b> softkey.
4.	Select the <b>Select Signal</b> softkey again or use the left/right arrow keys to turn the signal ON (by selecting a parameter) or OFF.
5.	Select the <b>Next Channel</b> softkey.
6.	Repeat steps 2-5 for other modules you wish to change.
7.	Press the <b>Main Screen</b> key to finish and to display the selected signals.

### Alarms

VueLink defines two types of alarms: red and yellow; but a TCM monitor only gives yellow alarms.

A **yellow** alarm indicates a situation where a response of the medical staff is necessary, but which is of less critical importance than a red alarm.

An alarm is displayed as a text on a yellow background at the top center of the PPMS monitor screen, e.g. "TC ALARM", which means that an alarm limit has been exceeded. Details about the specific alarm condition can be seen on the TCM monitor.

### Alerts

Alerts are caused by errors related to the TCM monitor or its accessories.

An alert is displayed as a text on a green background at the top left-hand corner of the PPMS monitor screen.

The following alerts may be displayed on the PPMS monitor screen:

Alert	Interpretation
SEE TC MONITOR	Details about the specific alert condition can be seen on the TCM monitor.
TC BAT LOW	Necessary to connect to mains.

*Continued on next page*

## Continuous data output: VueLink, *Continued*

### Alerts (*continued*)

Depending on the alert, parameter values may be displayed in the following ways:

Parameter reading	Interpretation
? (parameter value)	Data may be wrong
-?-	Provided data is wrong
(Blank)	No data can be provided

**NOTICE:** Alarm/alert delays between the PPMS and TCM monitors are less than two seconds.

### Troubleshooting

Symptom	Cause	Recommended action(s)
It is not possible to establish contact between the TCM and PPMS monitors	Error related to physical connection	<ul style="list-style-type: none"> <li>• Check that both monitors are ON</li> <li>• Check the cable connection</li> <li>• If error remains, contact authorized service personnel</li> </ul>
	Incorrect installation of TCM monitor	<ul style="list-style-type: none"> <li>• Check that the TCM monitor is set up for VueLink transmission</li> <li>• If error remains, contact authorized service personnel</li> </ul>
	Incorrect installation of PPMS monitor	<ul style="list-style-type: none"> <li>• Reset the PPMS to factory defaults, or</li> <li>• Select <b>Module Setup</b> and set "AUXILIARY PLUS 1" or "VueLink-B" to "ON"</li> <li>• If error remains, contact authorized service personnel</li> </ul>



## Continuous data output: MonLink

**Introduction** Enables TCM transmission of real-time data to external equipment. Any interaction is initiated by host request that includes:

Request	Reply
Status	System/parameter status
Measuring values	Measuring values and system/parameter status
Parameter ranges	Analog ranges, selected alarm ranges and parameter alarm status

**Available measured parameters**

TCM configuration	Parameter				
	pCO <sub>2</sub>	pO <sub>2</sub>	SpO <sub>2</sub>	PR	Power
TCM4 monitor	x	x			x
TCM40 monitor	x	x	x	x	x

Detailed alarm and alert notifications are not communicated but have to be revealed on the TCM monitor display.

All data values are updated every 2 seconds. How status and measured values are presented is decided by and is under the responsibility of the vendor of the external equipment.

**NOTICE:** For information on alarm delay and source/identification on external equipment, please refer to the user instructions for the relevant external equipment.

**Detailed communication protocol**

For the detailed communication protocol, please see the document *TCM Communication Protocol Specifications* (code no. 994-038) from Radiometer.

## Data export: serial

**Introduction** With the data export option, a dump of the trend data can be exported to an external PC and presented in spreadsheet format.

**Exporting data to PC** In this example, HyperTerminal version 690170 with Excel version 97 SR-2 is used.


Step	Action
------	--------

1. Connect the TCM4xx ETX Serial Adapter (code no. 636-649) to the serial port on the TCM monitor and to the PC.
2. On the PC, select **Start** → **Programs** → **Accessories** → **Communications** → **HyperTerminal** and open the EXE file.
3. Type in a name for the connection, e.g. TcData, and click **OK**.



Connect using COM1 and click **OK**.

5. Select the following port settings and click **OK**.



*Continued on next page*

## Data export: serial, *Continued*

### Exporting data to PC (continued)

- | Step | Action   |
|------|--|
| 6.   | Collect the data in HyperTerminal by selecting <b>Transfer</b> → <b>Capture text</b> .<br><br>Then name the file, e.g. CAPTURE.txt, and click <b>Start</b> .   |
| 7.   | On the TCM monitor, press <b>Setup</b> → <b>Technical</b> , enter password and press <b>Enter</b> and <b>Tech. settings</b> .  |
| 8.   | Set "Continuous data output" to "OFF", "Data export" to "Serial" and "Data export interval" to 2, 10, 30 or 60 seconds. Then press <b>OK</b> .<br><br>The data is exported to HyperTerminal and can be seen on the PC. |
| 9.   | On the PC, stop the collection of data to the file in HyperTerminal by selecting <b>Transfer</b> → <b>Capture text</b> → <b>Stop</b> .   |

The following shows an extract of a screen from HyperTerminal:

```
TCM4/40 Data Export
15:25:25 30/11/07 - 15:26:56 30/11/07


Patient ID: 30112007-0015

Last calibration: 15:24:53 30/11/07 - O2 : 157 mmHg - CO2 : 56 mmHg

Time;Event;O2 mmHg;CO2 mmHg;Temp °C;Power mW;SpO2 %;Pulse bmp;
15:25:30 30/11/07;;70 ;38 ;42.0;227 ;98;71;
15:25:40 30/11/07;;71 ;38 ;42.0;227 ;97;72;
15:25:50 30/11/07;;71 ;38 ;42.0;228 ;98;71;
15:26:00 30/11/07;;71 ;38 ;42.0;228 ;98;72;
15:26:10 30/11/07;;70 ;38 ;42.0;227 ;98;71;
15:26:20 30/11/07;;70 ;38 ;42.0;228 ;98;72;
15:26:30 30/11/07;;71 ;38 ;42.0;228 ;97;71;
15:26:40 30/11/07;;70 ;38 ;42.0;227 ;98;72;
15:26:50 30/11/07;;70 ;38 ;42.0;227 ;98;71;
```

### Importing data files into Microsoft Excel

- | Step | Action   |
|------|--|
| 1.   | Open Excel and select <b>Data</b> → <b>Get External Data</b> → <b>Import Text File</b> . |
| 2.   | Select the file CAPTURE.txt and click <b>Import</b> .                                    |
| 3.   | Select the following settings and click <b>Next</b> .                                    |




*Continued on next page*


## Data export: serial, *Continued*

Importing data  
files into  
Microsoft Excel  
(*continued*)

- | Step | Action  |
|------|---|
| 4.   | Select the following settings and click <i>Finish</i> . |



- |    |   |
|----|---|
| 5. | Select the following settings and click <i>OK</i> . |
|----|---|



- |    |                            |
|----|----------------------------|
| 6. | Data is imported to Excel. |
|----|----------------------------|

The following shows an extract of a screen from Excel:

	A	B	C	D	E	F	G	H	I
1	TCM4/40 Data Export								
2	15:25:25 30/11/07 - 15:26:56 30/11/07								
3	Patient ID: 30112007-0015								
4	Last calibration: 15:24:53 30/11/07 - O2 : 157 mmHg - CO2 : 56 mmHg								
5	Time	Event	O2 mmHg	CO2 mmHg	Temp °C	Power mW	SpO2 %	Pulse bmp	
6	15:25:30 30/11/07		70	38	42.0	227	98	71	
7	15:25:40 30/11/07		71	38	42.0	227	97	72	
8	15:25:50 30/11/07		71	38	42.0	228	98	71	
9	15:26:00 30/11/07		71	38	42.0	228	98	72	
10	15:26:10 30/11/07		70	38	42.0	227	98	71	
11	15:26:20 30/11/07		70	38	42.0	228	98	72	
12	15:26:30 30/11/07		71	38	42.0	228	97	71	
13	15:26:40 30/11/07		70	38	42.0	227	98	72	
14	15:26:50 30/11/07		70	38	42.0	227	98	71	

- |    |                             |
|----|-----------------------------|
| 7. | Save (Log Excel worksheet). |
|----|-----------------------------|

## Data export: USB

**Introduction** With this option, patient information and a dump of the trend data can be exported to a memory stick through one of the USB ports.

Exporting data to a memory stick	Step	Action
	1.	Connect a memory stick to one of the USB ports on the TCM monitor.
	2.	Press <b>Setup</b> → <b>Technical</b> (enter password and press <b>Enter</b> ) → <b>Tech. Settings</b> .
	3.	Set "Data export" to "USB" and "Data export interval" to 2, 10, 30 or 60 seconds. Then press <b>OK</b> .
	4.	Press <b>ID</b> to access Patient DMS.
	5.	Select the patient ID(s)/session(s) to export and press <b>Export</b> . The selected data will be exported to the USB port.

### NOTICES:

- If a memory stick is not connected to the monitor, a dialog box will be displayed, asking the user to insert a USB storage device. The user can either cancel the data export by pressing **Cancel** or insert a memory stick and press **OK** to initiate the data export.
- While data is being exported, a status dialog will be displayed. It will not be possible to operate the monitor until the export of data has been completed. Removing the memory stick before the data export has been completed will result in incomplete data export.

## Alarms


**Purpose** The alarm system informs the user about physiological and technical errors, or it gives text messages about recommended actions.

**NOTICE:** Every time the monitor is turned on, it performs a test of the alarm system to check that the sound is working.

**Definitions** The alarm system contains three different "alarms":

- Alarms = physiological alarms: one or more of the parameter values have exceeded or are equal to the high or low alarm limits
- Alerts = technical alarms: inform the user if e.g. an error has been detected during calibration
- Messages = pure text, e.g. "Ready"


**Alarms** Alarms are enabled/disabled in Parameter setup.

**NOTICE:** If the  symbol is displayed on the screen, all visual and acoustic alarm indications are disabled. Alerts are not affected.

It is possible to choose between two modes (see Technical settings in chapter 3: *Menu structure and setup programs*):

- Latching: The monitor remains in alarm status even though the alarm condition ceases to exist. So even if all parameter values are back within the alarm limits, the user must reset the alarm.
- Non-latching: The monitor resets itself as soon as the alarm condition ceases to exist.

An alarm is made up of visual indications (the parameter value and an alarm text will be flashing) and an acoustic indication (a discontinuous tone).




*Continued on next page*

## Alarms, *Continued*

### Alerts

An alert cannot be turned off but disappears automatically when the alert condition ceases to exist.

An alert is made up of visual indications ("Alert" displayed in the headline field and a text explaining the problem in the curve view field of the affected parameter) and an acoustic indication (a discontinuous tone).



### Alarm silence

It is possible to silence an alarm or an alert for two minutes by pressing the **Alarm silence** touch key. A progress bar will then be displayed on the **Alarm silence** touch key, showing the progress of the alarm silence period. However, if an alarm/alert is activated for another parameter during the silence period, this will discontinue the silence period and activate the alarm/alert sounds.

**NOTICE:** When affixing the sensor to the patient, i.e. when the monitor goes from Ready to Measuring mode, the alarm silence period starts automatically and lasts for 10 minutes.

### Alarm reset


When the alarm silence period times out, the alarm system is reset.

The alarm system can also be reset by pressing the **Alarm silence** touch key twice (or once, if the alarm silence is already activated).

If the alarm/alert condition is still present when the alarm system is reset, the visual indications will remain for both types of alarms, whereas the acoustic indication will be reactivated for an alarm condition only.

### Alarm tone

An alarm tone consists of bursts of 10 short pulses. Between two bursts there is a 3-second pause. The alarm tone continues until the alarm condition ceases to exist or Alarm silence is activated.



*Continued on next page*

## Alarms, *Continued*

### Alert tone

An alert tone consists of two pulses followed by a 5-second pause. The alert tone continues until the alert condition ceases to exist or Alarm silence is activated.



### End-of-calibration tone

An End-of-calibration tone consists of two pulses.



### TCM40 SatSeconds

With traditional alarm management, upper and lower alarm limits are set for monitoring oxygen saturation. During monitoring, as soon as an alarm limit is violated by as little as one percentage point, an audible alarm immediately sounds. When the SpO<sub>2</sub> level fluctuates near an alarm limit, the alarm sounds each time the limit is violated. Such frequent alarms can be distracting.

Therefore, the TCM40 SpO<sub>2</sub> module utilizes SatSeconds alarm management. With the SatSeconds technique, upper and lower alarm limits are set in the same way as traditional alarm management. However, the clinician also sets a SatSeconds limit that allows monitoring of SpO<sub>2</sub> below the selected lower alarm limit and above the selected upper alarm limit for a period of time before an audible alarm sounds.

The SatSeconds limit controls the time that the SpO<sub>2</sub> level may violate the alarm limits before an audible alarm sounds. SatSeconds is enabled/disabled in SpO<sub>2</sub> parameter setup.

The method of calculation is as follows:

The number of percentage points that the SpO<sub>2</sub> falls outside the alarm limit is multiplied by the number of seconds that the SpO<sub>2</sub> level remains outside that limit. This can be stated as an equation:

$$\text{points} \times \text{seconds} = \text{SatSeconds}$$

where:

- points = SpO<sub>2</sub> percentage points outside of the limit
- seconds = number of seconds that SpO<sub>2</sub> remains at that point outside of the limit

The alarm response time, assuming a SatSeconds limit set at 60 and a lower alarm limit set at 90, is described and illustrated below.

In this example, the SpO<sub>2</sub> level drops to 88 (2 points) and remains there for a period of 2 seconds (2 points  $\times$  2 seconds = 4). The SpO<sub>2</sub> then drops to 86 for 3 seconds and then to 84 for 8 seconds. The resulting SatSeconds are:

SpO <sub>2</sub>	Seconds	SatSeconds
2 $\times$	2 =	4
4 $\times$	3 =	12
6 $\times$	8 =	48
<b>Total SatSeconds =</b>		<b>64</b>


*Continued on next page*



## Alarms, *Continued*

### TCM40 SatSeconds (*continued*)

After approximately 13 seconds, the SatSeconds alarm will sound, because 60 SatSeconds will have been exceeded. See arrow (↑) in the figure below.



Saturation levels may fluctuate rather than remain steady for a period of several seconds. Often, SpO<sub>2</sub> levels may fluctuate above and below the alarm limit, re-entering the non-alarm range several times.

During such fluctuation, the TCM40 monitor integrates the number of SpO<sub>2</sub> points, both positive and negative, until either the SatSeconds limit (SatSeconds setting) is reached or the SpO<sub>2</sub> level returns to within a normal range and remains there.

### SatSeconds "safety net"

The SatSeconds "safety net" is for patients whose saturation levels frequently exceed the high or low alarm limit but do not stay outside the limit long enough for the SatSeconds setting to be reached. When three or more limit violations occur in 60 seconds, an alarm will sound even if the SatSeconds setting has not been reached.

## How to print

### Connecting a printer

Step	Action
1.	Shut down the monitor by following the procedure described in chapter 4: <i>Installation and maintenance</i> .
2.	Connect an HP printer with PCL3 protocol to one of the USB ports on the rear of the monitor.
3.	Switch on the monitor.

### Printing reports

Step	Action
1.	If required, go into Printer setup to change the printer settings (see chapter 3: <i>Menu structure and setup programs</i> ).
2.	Press <b>Print</b> and set the printer start and stop times. <b>NOTICES:</b> <ul style="list-style-type: none"> <li>• The monitor must be connected to a printer with PCL3 protocol.</li> <li>• The memory function holds a total of 48 hours of monitoring data.</li> <li>• The printed report will be performed in accordance with the printer settings selected in the printer setup.</li> <li>• The time span selected in the printer setup will influence the start time of the Printer start/stop time screen, as the interval between the printer start and stop times corresponds to the time span; i.e. if the time span is set to 1 hour, the interval between the printer start and stop times will also be 1 hour (See <i>Printer setup</i> in chapter 3).</li> <li>• The scroll interval depends on the time span selected in the printer setup. If the time span is set to less than 1 hour, the scroll interval will correspond to the exact time span (e.g. 5, 15 or 30 min). If the time span is set to 1 hour or more, the scroll interval will be 1 hour.</li> <li>• The stop time is set to current time and then adjusted to match the time of the latest data record. The stop time can therefore differ from current time with up to 2 hours, depending on the time interval selected in the printer setup.</li> </ul>
3.	Press <b>OK</b> to make a printout, or press <b>Cancel</b> to exit the dialog box without making a printout. <b>NOTICES:</b> <ul style="list-style-type: none"> <li>• The different kinds of event marks are only shown on the printout if they are actually present.</li> <li>• Even if both report types are selected in Printer setup, there will only be one "Comments" page.</li> </ul>

*Continued on next page*

## How to print, *Continued*

### Example of table report


#### RADIOMETER TCM4 SERIES




### Transcutaneous measurement (tcpCO<sub>2</sub>/tcpO<sub>2</sub>/SpO<sub>2</sub>/Pulse)

Measuring unit:	tcpCO <sub>2</sub> /tcpO <sub>2</sub>	mmHg	Sensor temperature:	42.0 °C
	Power	mW	Last calibration	01:27:07 02/03/2004
	SpO <sub>2</sub>	%	Last calibration value:	pCO <sub>2</sub> 80 mmHg
	Pulse	bpm		pO <sub>2</sub> 160 mmHg

Facility name: \_\_\_\_\_

Patient name/ID: \_\_\_\_\_

 Blood gas event mark
 \* In vivo calibration
 - No value available

Time/date	Event	tcpCO <sub>2</sub>	tcpO <sub>2</sub>	Power	SpO <sub>2</sub>	Pulse
02:23:00 02/03/04		34	158	3	97	71
02:24:00 02/03/04		34	158	7	97	71
02:25:00 02/03/04	1	34	158	11	97	71
02:26:00 02/03/04		34	158	-	97	71
02:27:00 02/03/04		-	158	16	96	70
02:28:00 02/03/04	2	34	158	19	96	70
02:29:00 02/03/04		34	158	23	97	72
02:29:13 02/03/04		-	-	-	97	72
02:30:00 02/03/04		34	158	26	97	72
02:31:00 02/03/04	3	34	158	30	97	72
02:32:00 02/03/04		34	158	31	97	72
02:33:00 02/03/04		34	158	33	97	72
02:34:00 02/03/04		34	158	37	97	72
02:35:00 02/03/04		34	158	-	97	72
02:36:00 02/03/04		-	158	40	97	72
02:37:00 02/03/04		34	158	43	97	72
02:38:00 02/03/04		34	158	46	97	72
02:39:00 02/03/04		34	158	51	97	72
02:40:00 02/03/04		34	158	52	97	72
02:41:00 02/03/04		34	158	52	97	72

Printed: 02:53:07 02/03/2004

Page 1

*Continued on next page*

**Example of  
table report  
(continued)**

Transcutaneous measurement (tcpCO<sub>2</sub>/tcpO<sub>2</sub>/SpO<sub>2</sub>/Pulse)

Patient name/ID: \_\_\_\_\_

[illegible]

Physician's signature: \_\_\_\_\_ Date: \_\_\_\_\_

## How to print, *Continued*

### Example of curve report

#### RADIOMETER TCM4 SERIES

#### Transcutaneous measurement (tcpCO<sub>2</sub>/tcpO<sub>2</sub>/SpO<sub>2</sub>/Pulse)



Measuring unit:	tcpCO <sub>2</sub> /tcpO <sub>2</sub>	mmHg	Sensor temperature:	42.0 °C
	Power	mW	Last calibration	01:27:07 02/03/2004
	SpO <sub>2</sub>	%	Last calibration value:	pCO <sub>2</sub> 80 mmHg
	Pulse	bpm		pO <sub>2</sub> 160 mmHg

Facility name: \_\_\_\_\_

Patient name/ID: \_\_\_\_\_

● Blood gas event mark

\* In vivo calibration




Printed: 02:53:07 02/03/2004

Page 1

*Continued on next page*

## How to print, *Continued*

### Example of curve report (*continued*)



Printed: 02:53:07 02/03/2004

Page 2

*Continued on next page*

**Example of  
curve report  
(continued)**

Transcutaneous measurement (tcpCO<sub>2</sub>/tcpO<sub>2</sub>/SpO<sub>2</sub>/Pulse)

---

Physician's signature: \_\_\_\_\_ Date: \_\_\_\_\_

## Blood gas comparison

**Introduction** It is possible to key in the  $p\text{CO}_2$  and  $p\text{O}_2$  blood gas values of a patient to compare these with transcutaneous measurements from the same patient. The blood gas values are displayed as blood drops in all views as well as on all printed reports.


### Keying in blood gas value

- | Step | Action  |
|------|---|
| 1.   | Make sure that the sensor is placed on the patient and that the $\text{tcpCO}_2$ and $\text{tcpO}_2$ values are stable.   |
| 2.   | Before drawing the blood sample, press <b>Setup</b> → <b>Parameter</b> → <b>Blood gas</b> → <b>Sample time</b> . The monitor will store the tc values at that time and the $p\text{CO}_2$ and $p\text{O}_2$ input fields and the numeric keypad are now accessible. |
| 3.   | Draw the blood sample and perform a blood gas analysis.   |
| 4.   | When the blood gas results are ready, enter the Blood gas setup, press the $p\text{CO}_2$ input field and key in the value with the numeric keypad.   |

#### NOTICES:

- The **←** key deletes one character at a time from the right.
  - If the  $p\text{CO}_2$  or  $p\text{O}_2$  input field is dimmed, the parameter is not available, either because of the sensor type or because  $p\text{O}_2$  is disabled in Technical setup.
5. Press the  $p\text{O}_2$  input field and key in the value with the numeric keypad, and then press **OK**.

The blood gas values will be displayed as blood drops in Normal view as in the following example:





## In vivo calibration

### Introduction

It is possible to perform in vivo calibrations, i.e. alterations in the transcutaneous readings, using the results from an arterial blood gas analysis. In vivo calibrations can either be performed on  $p\text{CO}_2$ ,  $p\text{O}_2$  or both parameters.

When the blood gas values have been entered on the TCM4/40 monitor, the new calibration lines are calculated.

The following formulas are used:


$$\text{tcpCO}_2(\text{corr}) = \text{tcpCO}_2 + [p\text{CO}_2(\text{a}) - \text{tcpCO}_2(\text{t})]$$

$$\text{tcpO}_2(\text{corr}) = \text{tcpO}_2 \times [p\text{O}_2(\text{a}) / \text{tcpO}_2(\text{t})]$$


where

- $\text{tcpCO}_2$  is either the measured  $\text{CO}_2$  value or the measured  $\text{CO}_2$  value with correction (if corrections are selected in Technical settings)
- $p\text{CO}_2(\text{a})$  and  $p\text{O}_2(\text{a})$  are the measured arterial values
- $\text{tcpCO}_2(\text{t})$  and  $\text{tcpO}_2(\text{t})$  are the transcutaneous values measured at time t, which is the approximate time when the arterial sample was taken and the **Sample time** touch key was pressed.

On the  $p\text{CO}_2$  part, the corrected (CORR) value line will have the correct ( $45^\circ$ ) slope. The in vivo calibration, therefore, gives an offset for all values – i.e. it moves the line – just as changing the metabolic factor in Technical settings would do.



On the  $p\text{O}_2$  part of the sensor, the in vivo calibration changes the sensitivity (the slope) of the  $p\text{O}_2$  signal, just as with normal calibrations.




*Continued on next page*

## In vivo calibration, *Continued*

### Introduction (*continued*)

#### NOTICES:

- In vivo calibration should not be performed on hemodynamically unstable patients, because the arterial blood gas value and the tc value may fluctuate considerably.
- In vivo calibration of a tc value will bring the displayed value closer to the arterial sample reading. But the correlation between the tc and the arterial values will remain unchanged.
- In vivo calibration does not make the TCM monitor read arterial values. Nor does it eliminate the blood flow dependence on the tc value or any of the other determinants. It gives the same information as the original tc value – only on another level (see the graph below).



### Performing an in vivo calibration

Step	Action
1.	Make sure that access to the <b><i>In vivo calibration</i></b> touch key has been enabled in Technical settings, that the sensor is placed on the patient and that the tc values are stable.
2.	Before drawing the blood sample, press <b>Setup</b> → <b>Parameter</b> → <b>Blood gas</b> → <b>Sample time</b> . The monitor will register the blood sampling time and the pCO <sub>2</sub> and pO <sub>2</sub> input fields and the numeric keypad are now accessible.
3.	Draw the blood sample and perform a blood gas analysis.
4.	When the blood gas result is ready, enter the Blood gas setup, press the pCO <sub>2</sub> input field and key in the value with the numeric keypad.
5.	Then press the pO <sub>2</sub> input field and key in the pO <sub>2</sub> value.
6.	Press <b><i>In vivo calibration</i></b> and <b>OK</b> .


*Continued on next page*

## In vivo calibration, *Continued*

### Performing an in vivo calibration (*continued*)

Step	Action
------	--------

	The curves will now be marked with *, indicating that the measured values are in vivo calibrated and no longer show the original transcutaneous values.
--	---



### NOTICES:

- The in vivo calibration procedure can be repeated during the monitoring period if you wish to alter the values again according to the results of new blood samples.
- To return to the original tc values, simply put the sensor back in the calibration chamber.



## 7. Troubleshooting

The TCM4/40 systems .....	7-2
---------------------------	-----

## The TCM4/40 systems

### Introduction

The monitor contains Radiometer-developed software that has been developed, tested and released according to our certified Quality Assurance System in order to minimize hazards arising from the software. Furthermore, the status of the system is continuously monitored during operation. Should a problem or error occur, it is automatically recorded and presented to the operator on the screen.

This chapter describes possible errors, their causes and the recommended operator actions.

**NOTICE:** The monitor cover should be removed only by authorized service personnel. There are no user-serviceable parts inside.

### Error symptoms

Symptom	Cause	Recommended action(s)
Monitor not operational	<ul style="list-style-type: none"> <li>• Monitor is defective</li> <li>• Battery is exhausted</li> <li>• Monitor is not connected to mains</li> </ul>	<ul style="list-style-type: none"> <li>• Check the connections</li> <li>• Contact authorized service personnel</li> </ul>
Measured values drift when sensor is measuring on cal gas (in vitro)	<ul style="list-style-type: none"> <li>• Insufficient removal of old electrolyte during membranizing procedure</li> <li>• Air bubbles in electrolyte</li> <li>• The sensor membrane has not been changed regularly</li> <li>• Defective sensor membrane (holes or scratches)</li> <li>• Defective O-rings (one or both O-rings are broken)</li> </ul>	<ul style="list-style-type: none"> <li>• Remembrane the sensor</li> <li>• Check sensor zero current and sensitivity</li> </ul>
Measured values drift constantly without any clinical cause when sensor is mounted on patient (in vivo)	<ul style="list-style-type: none"> <li>• Improper attachment of the sensor</li> <li>• Improper attachment of the fixation ring</li> <li>• An inappropriate measurement site has been selected</li> <li>• Inadequate vasodilation</li> </ul>	<ul style="list-style-type: none"> <li>• Detach the sensor and reapply correctly</li> <li>• Detach the fixation ring and reapply correctly</li> <li>• Check sensor zero current and sensitivity</li> </ul>
	<ul style="list-style-type: none"> <li>• Hole in the membrane</li> <li>• Air under the membrane</li> </ul>	Remembrane the sensor

*Continued on next page*

## The TCM4/40 systems, *Continued*

### Error symptoms (*continued*)

Symptom	Cause	Recommended action(s)
Measured values not stable or out of range 20 minutes after application	<ul style="list-style-type: none"> <li>• Patient status unstable</li> <li>• Inadequate vasodilation</li> <li>• Improper attachment of the sensor</li> </ul>	<ul style="list-style-type: none"> <li>• Evaluate patient status</li> <li>• Remount the sensor</li> </ul>
	<ul style="list-style-type: none"> <li>• Sensor exposed to high ambient light</li> <li>• Hole in the membrane</li> <li>• Air under the membrane</li> </ul>	Remembrane the sensor
Measured values change suddenly without any clinical cause	<ul style="list-style-type: none"> <li>• Inaccurate result due to patient movements</li> <li>• Air leakage under fixation ring</li> <li>• Self-adhesive ring does not stick to the skin</li> </ul>	Recalibrate and reapply sensor to the skin, possibly selecting a new measuring site
	<ul style="list-style-type: none"> <li>• Hole in the membrane</li> <li>• Air under the membrane</li> </ul>	Remembrane the sensor
Noise appears on the tension readout	Interference from nearby equipment	Increase the distance from the interfering equipment

### Monitor error message

Error message	Cause	Recommended action
Battery level low	Battery level is low	Connect the monitor to the mains

### tc error messages – in alphabetical order

Error message	Cause	Recommended action(s)
Calibration error. Barometer error during calibration.	Barometer value invalid when collected at start of calibration	<ul style="list-style-type: none"> <li>• Start a new calibration</li> <li>• If error remains, contact authorized service personnel</li> </ul>
Calibration error. Gas flow out of range.	<ul style="list-style-type: none"> <li>• Calibration chamber is blocked</li> <li>• Calibration unit is defective</li> </ul>	<ul style="list-style-type: none"> <li>• Check chamber sealing</li> <li>• Contact authorized service personnel</li> </ul>
Calibration error. Gas level low.	Less than 10 calibrations left	Change gas bottle

*Continued on next page*

## The TCM4/40 systems, *Continued*

tc error  
messages – in  
alphabetical  
order  
(*continued*)

Error message	Cause	Recommended action(s)
Calibration error. Leak error in calibration chamber.	There is a leak in the calibration chamber	<ul style="list-style-type: none"> <li>• Check that a gasket is present and positioned correctly in the chamber</li> <li>• Check that the sensor is positioned correctly in the chamber</li> </ul>
Calibration error. No sensor connected to the chamber.	No sensor detected in the calibration chamber at start of or during calibration	Place a tcpCO <sub>2</sub> /tcpO <sub>2</sub> sensor in the chamber
Calibration error. Sensor drift during calibration.	The sensor has been rejected during the drift check	<ul style="list-style-type: none"> <li>• Start a new calibration</li> <li>• If error remains, remembrane the sensor, and start a new calibration.</li> <li>• If error remains, the sensor is defective. Replace with a new one.</li> </ul>
Calibration error. Sensor sensitivity error during calibration.	The sensor sensitivity is outside the specifications	<ul style="list-style-type: none"> <li>• If the sensor has just been remembraned, start a new calibration</li> <li>• Otherwise remembrane the sensor and start a new calibration</li> <li>• Sensor is defective. Replace with a new one.</li> </ul>
Calibration error. Sensor temperature error during calibration. This may be due to an error in the sensor module.	The selected sensor temperature value cannot be reached	<ul style="list-style-type: none"> <li>• Calibrate the sensor</li> <li>• If error remains, the sensor is defective. Replace with a new one.</li> </ul>
Calibration required	<ul style="list-style-type: none"> <li>• The last calibration has been aborted</li> <li>• 12 hours have passed since the last successful calibration</li> <li>• SmartCal has timed out</li> </ul>	Perform a calibration

*Continued on next page*



## The TCM4/40 systems, *Continued*

**tc error  
messages – in  
alphabetical  
order  
(continued)**

Error message	Cause	Recommended action(s)
Communication error. Sensor module disconnected. Call service.	For unknown reasons, the communication between the monitor and the tc module has not been satisfactory	<ul style="list-style-type: none"> <li>• If error remains, try with another module</li> <li>• If error remains on the new module, contact authorized service personnel</li> </ul>
General error. Call service.	The tc module or sensor is not functioning	<ul style="list-style-type: none"> <li>• Try with another module or sensor</li> <li>• If error remains, contact authorized service personnel</li> </ul>
Heater failure	<ul style="list-style-type: none"> <li>• The measured sensor power is outside the measuring range</li> </ul>	<ul style="list-style-type: none"> <li>• Unplug the sensor from the sensor socket at the rear of the module, and then reconnect it. Calibrate the sensor.</li> <li>• If error remains, the sensor is defective. Replace with a new one.</li> </ul>
	The selected sensor temperature cannot be reached	Check if the patient or the surroundings are very cold
Invalid sensor connected	No tcpCO <sub>2</sub> /tcpO <sub>2</sub> sensor detected	Change the sensor
No sensor connected	<ul style="list-style-type: none"> <li>• No sensor detected in the calibration chamber</li> <li>• The sensor in the calibration chamber is not functioning</li> </ul>	Connect a functional tcpCO <sub>2</sub> /tcpO <sub>2</sub> sensor
pCO <sub>2</sub> out of range	<ul style="list-style-type: none"> <li>• The measured pCO<sub>2</sub> value is outside the measuring range</li> <li>• The pCO<sub>2</sub> sensor is low on battery power</li> </ul>	<ul style="list-style-type: none"> <li>• Change the membrane</li> <li>• Change the sensor</li> </ul>
pO <sub>2</sub> out of range	The measured pO <sub>2</sub> value is outside the measuring range	<ul style="list-style-type: none"> <li>• Change the membrane</li> <li>• Change the sensor</li> </ul>
Site time end	The site time counter has reached zero	<ul style="list-style-type: none"> <li>• Reposition the sensor on the patient</li> <li>• Recalibrate the sensor</li> </ul>

*Continued on next page*

## The TCM4/40 systems, *Continued*

### tc error messages – in alphabetical order (*continued*)

Error message	Cause	Recommended action(s)
Temperature failure	<ul style="list-style-type: none"> <li>Measured temperature is outside the measuring range</li> <li>The two sensor thermistors give different temperatures</li> <li>The sensor temperature is too high</li> </ul>	<ul style="list-style-type: none"> <li>Calibrate the sensor</li> <li>If error remains, the sensor is defective. Replace with a new one.</li> </ul>
Temperature indication is flashing (No text is shown)	The selected sensor temperature has not yet been reached	Wait until it is reached

### SpO<sub>2</sub> error messages – in alphabetical order

Error message	Cause	Recommended action(s)
Communication error. SpO <sub>2</sub> module disconnected. Call service.	For unknown reasons, the communication between the monitor and the SpO <sub>2</sub> module has not been satisfactory	<ul style="list-style-type: none"> <li>If error remains, try with another module</li> <li>If error remains on the new module, contact authorized service personnel</li> </ul>
General error. Call service.	SpO <sub>2</sub> module or sensor is not functioning	<ul style="list-style-type: none"> <li>Try with another module or sensor</li> <li>If error remains, contact authorized service personnel</li> </ul>
No sensor connected	<ul style="list-style-type: none"> <li>No SpO<sub>2</sub> sensor connected to the monitor</li> <li>Invalid sensor connected to the monitor</li> </ul>	Connect an SpO <sub>2</sub> sensor to the monitor

### Dialog box messages – in alphabetical order

A number of messages (system messages and alerts) appear in dialog boxes.

System message	Cause
Battery level critically low	Less than 5 minutes running time on battery
Battery level is critically low. Setup changes cannot be saved.	The user has made changes in the setup while the battery level was critically low

*Continued on next page*

## The TCM4/40 systems, *Continued*

**Dialog box  
messages – in  
alphabetical  
order  
(continued)**

System message	Cause
Battery not connected. For data safety reasons, a proper battery must always be connected.	No battery is connected to the monitor
Battery not connected. Setup changes cannot be saved.	The user has made changes in the setup while the battery was not connected
Blood gas setup is only available in measuring mode. You cannot enter this menu.	The user has tried to enter the Blood gas setup while the monitor is not monitoring
Data error. Shut down system immediately.	RAM failure on the PC unit
Data export completed	Export of data to the USB storage device has been completed
Data export failed	Export of data to the USB storage device has failed. E.g. due to missing USB storage device.
Data from more than one patient cannot be selected into a view	The user has tried to load sessions from different patient IDs into a view
Date and time cannot be set during measurement	The user has tried to enter the Date/time setup during measurement
Exporting data. Please wait.	Data is being exported to the USB storage device
Incorrect password. Try again.	The user has tried to enter the Technical setup with an incorrect password
Incorrect time. Try again.	The user has entered an incorrect time
Insert USB storage device	Export to USB port was selected without a memory stick being connected to the USB port
Monitor temperature too high. Please shut down system immediately.	The temperature in the CPU is too high

*Continued on next page*

## The TCM4/40 systems, *Continued*

### Dialog box messages – in alphabetical order (continued)

System message	Cause
No sessions selected. Select one or more sessions.	The user has tried to view session information (in Patient DMS) without selecting a session
One or more sessions will be deleted from the system. This action cannot be undone.	The user has selected one or more sessions in Patient DMS and pressed <b>Delete</b> or <b>Delete all</b>
$pO_2/pCO_2$ not enabled. You cannot enter this menu.	<ul style="list-style-type: none"> <li>The <math>pO_2</math> parameter is disabled in Technical setup, and the user is trying to enter the <math>pO_2</math> parameter setup</li> <li>A <math>tcpO_2</math> sensor is connected, and the user is trying to enter the <math>pCO_2</math> parameter setup</li> </ul>
Saving data. Please wait.	The monitor is saving data before shutting down
The monitor is measuring. Current session cannot be deleted.	The user has tried to delete the session in progress
The monitor is measuring. No other sessions can be loaded into the view.	The user has tried to select one or more sessions – other than the one in progress – into a view
The sensor must be recalibrated if the temperature is changed. Do you want to change the temperature?	The user has tried to change the sensor temperature during Ready or Measuring state

### Pure text messages – in alphabetical order

Text message	Cause
Calibrating	The sensor is being calibrated
Ready	The sensor has been calibrated and is ready for use

## 8. Specifications and ordering information

Specifications .....	8-2
Accessories.....	8-14

## Specifications

### TCM4/40 monitor specifications

Item	Description
Measured parameters	tcpCO <sub>2</sub> , tcpO <sub>2</sub> , SpO <sub>2</sub> , pulse rate and sensor heating power
Display update period	SpO <sub>2</sub> plethysmograph: 25 mm/sec Numerical values: 2 sec pCO <sub>2</sub> /pO <sub>2</sub> graphs: 2 sec
Display range	tcpCO <sub>2</sub> : 0-200 mmHg or 0.0-26.7 kPa tcpO <sub>2</sub> : 0-800 mmHg or 0.0-99.9 kPa SpO <sub>2</sub> : 0-100 % Pulse rate: 0 and 20-300 beats per minute (bpm)
Measuring range	tcpCO <sub>2</sub> : 5-100 mmHg or 0.7-13.3 kPa tcpO <sub>2</sub> : 0-800 mmHg or 0.0-99.9 kPa SpO <sub>2</sub> : 70-100 % Pulse rate: 20-250 bpm
Sensor heating power range	10-650 mW ± 3 % of reading (< 10 mW ~ 0 mW)
Temperature settings in °C	37.0-45.0 °C in steps of 0.5 °C
Barometer	Built-in: 375-825 mmHg or 50-110 kPa Accuracy: ±5 mmHg or 0.67 kPa
Calibration	Calibration gas (7.5 % CO <sub>2</sub> , 20.9 % O <sub>2</sub> , balance N <sub>2</sub> )
Ambient relative humidity	20-80 %
Ambient temperature	Monitor: 5-40 °C tc sensors: 15-40 °C <b>NOTICES:</b> <ul style="list-style-type: none"> <li>The ambient temperature must always be at least 3 °C lower than the set sensor temperature.</li> <li>If SmartCal is selected, the maximum ambient temperature is 34 °C.</li> </ul> SpO <sub>2</sub> sensors: 5-40 °C
Transport and storage conditions	The monitor and modules can be transported and stored at -20 to +60 °C and < 95 % RH.
Data storage	Up to 48 hours of monitoring data in 2-second data intervals

*Continued on next page*

## Specifications, *Continued*

### TCM4/40 monitor specifications (continued)

Item	Description		
Computer	Screen: 6½" color touch TFT, full VGA (640 × 480) CPU: AMD ETX LX800, 500 MHz (Pentium Class) Software platform: Windows CE 5.0		
Interface connection	Serial output	EIA232, (RS232)	
	Printer output	USB 2.0 (compliant with USB 1.1)	
	Analog output	0-1 V	
Power supply	100-240 V	50-60 Hz	
Battery	Rechargeable Pb battery	Under normal conditions, the battery can operate for 1 hour before recharging is needed	
Dimensions of monitor	Height:	16 cm	6.3 in
	Width:	30.8 cm	12.1 in
	Depth:	23 cm	8.7 in
	Weight:	4 kg	8.8 lbs
Dimensions of tc module	Height:	10.7 cm	4.2 in
	Width:	14.5 cm	5.7 in
	Depth:	14.8 cm	5.8 in
	Weight:	0.6 kg	1.3 lbs
Dimensions of SpO <sub>2</sub> module	Height:	3.5 cm	1.4 in
	Width:	14.5 cm	5.7 in
	Depth:	14.8 cm	5.8 in
	Weight:	0.21 kg	0.5 lb
Alarm sound pressure	At highest alarm sound level	Alarm tone: Alert tone: End-of-calibration tone:	83 dBA 73 dBA 64 dBA
	At lowest alarm sound level	Alarm tone: Alert tone: End-of-calibration tone:	65 dBA 58 dBA 51 dBA

*Continued on next page*

## Specifications, *Continued*

**EMC approvals and compliance** The monitor is intended for use in the electromagnetic environment specified in the tables below. The customer or user of the monitor should assure that it is used in such an environment. The monitor complies with IEC 60601-1-2.

### Guidance and manufacturer's declaration – electromagnetic emissions

Emissions test	Compliance	Electromagnetic environment guidance
RF emissions CISPR 11	Group 1	The monitor's RF emissions are very low and they are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class A	The monitor is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	N/A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	N/A	

### Guidance and manufacturer's declaration – electromagnetic immunity

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	±2 kV for power supply lines ±1 kV for input/output lines	Mains power quality should be that of a typical commercial and/or hospital environment
Surge IEC 61000-4-5	±1 kV differential mode ±2 kV common mode	±1 kV differential mode ±2 kV common mode	Mains power quality should be that of a typical commercial and/or hospital environment

*Continued on next page*



## Specifications, *Continued*


### EMC approvals and compliance (*continued*)

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment guidance
Voltage dips, short interruptions, and voltage variations on power supply input lines IEC 61000-4-11	$< 5 \% U_T$ ( $> 95 \%$ drop in $U_T$ ) for 0.5 cycle  $40 \% U_T$ (60 % drop in $U_T$ ) for 5 cycles  $70 \% U_T$ (30 % drop in $U_T$ ) for 25 cycles  $< 5 \% U_T$ ( $> 95 \%$ drop in $U_T$ ) for 5 seconds	N/A, as the monitor has battery backup	None
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	None
<b>NOTICE:</b> $U_T$ is the AC mains voltage prior to application of the test level.			

*Continued on next page*

## Specifications, *Continued*

### EMC approvals and compliance (continued)

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment guidance
Conducted RF IEC 61000-4-6	3 V <sub>rms</sub> 150 kHz to 80 MHz	3 V <sub>rms</sub>	<p>Portable and mobile RF communications equipment should be used no closer to any part of the monitor, including cables, than the recommended separation distance calculated from the equation appropriate for the frequency of the transmitter.</p> <p><b>Recommended separation distance</b></p> $d = 1.2\sqrt{P} \quad 150 \text{ kHz to } 80 \text{ MHz}$ $d = 1.2\sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$ $d = 2.3\sqrt{P} \quad 800 \text{ MHz to } 2.5 \text{ GHz}$ <p>where <math>P</math> is the output power rating of the transmitter in watts (W) according to the transmitter manufacturer and <math>d</math> is the recommended separation distance in meters (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey<sup>a</sup>, should be less than the compliance level in each frequency range<sup>b</sup>.</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 
<p><sup>a</sup> Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast, and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the monitor is used exceeds the applicable RF compliance level above, the monitor should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the monitor.</p> <p><sup>b</sup> Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.</p> <p><b>NOTICES:</b></p> <ul style="list-style-type: none"> <li>• At 80 MHz and 800 MHz, the higher frequency range applies.</li> <li>• These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.</li> </ul>			

*Continued on next page*

## Specifications, *Continued*

### EMC approvals and compliance (continued)

#### Recommended separation distances between portable and mobile RF communications equipment and the TCM monitor

The monitor is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or user of the monitor can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the monitor as recommended below, according to the maximum output power of the communications equipment.

Rated output power of transmitter Watts	Separation distance according to frequency of transmitter in meters		
	150 kHz to 80 MHz $d = 1.2\sqrt{P}$	80 MHz to 800 MHz $d = 1.2\sqrt{P}$	800 MHz to 2.5 GHz $d = 2.3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23
<p>For transmitters rated at an output power not listed above, the recommended separation distance <math>d</math> in meters can be estimated using the equation applicable to the frequency of the transmitter, where <math>P</math> is the maximum output power rating of the transmitter in the corresponding column.</p> <p><b>NOTICES:</b></p> <ul style="list-style-type: none"> <li>At 80 MHz and 800 MHz, the higher frequency range applies.</li> <li>These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.</li> </ul>			

### Other approvals and compliance

Item	Description
Patient safety	<ul style="list-style-type: none"> <li>The instrument complies with IEC 60601-1 and IEC 60601-2-23.</li> <li>The following test house has approved the instrument: CSA in Canada according to CAN/CSA-C22.2 No. 601.1-M90, 601.1S1-94, 601.1B-98, 601.2.23-02 and UL std. No. 60601-1.</li> </ul>

*Continued on next page*

## Specifications, *Continued*

### Other approvals and compliance (*continued*)

Item	Description
Compliance	<p>The monitor complies with:</p> <ul style="list-style-type: none"> <li>• IEC-60601-1, Medical Electrical Equipment – Part 1: General requirements for safety</li> <li>• IEC-60601-1-2, Medical Electrical Equipment – Part 1-2: General requirements for safety. Collateral standard: Electromagnetic compatibility - requirements and tests</li> <li>• IEC-60601-1-4, General requirements for safety and design of programmable electrical medical systems</li> <li>• IEC-60601-2-23, Medical Electrical Equipment – Part 2-23: Particular requirements for the safety, including essential performance, of transcutaneous partial pressure monitoring equipment</li> <li>• IEC-60601-2-49, Medical Electrical Equipment – Part 2-49: Particular requirements for the safety of multifunction patient monitoring equipment.</li> <li>• IEC-60601-3-1, Medical Electrical Equipment – Part 3-1: Essential performance requirements for transcutaneous oxygen and carbon dioxide partial pressure monitoring equipment</li> <li>• ISO 9919, Medical Electrical Equipment. Particular requirements for the basic safety and essential performance of pulse oximeter equipment for medical use</li> <li>• IEC-60601-1-8, Medical Electrical Equipment – Part 1-8: General requirements for safety. Collateral standard: General requirements, tests and guidance for alarm systems in medical electrical equipment and medical electrical systems</li> <li>• Class II Special Controls Guidance Document: Cutaneous Carbon Dioxide (tcpCO<sub>2</sub>) and Oxygen (tcpO<sub>2</sub>) Monitors; Guidance for Industry and FDA (December 13, 2002)</li> <li>• Draft Guidance for Industry and Staff: Pulse oximeters premarket notifications submissions [510(k)s]</li> </ul>
pCO <sub>2</sub>	<p><b>Non-linearity and hysteresis:</b> The requirement to non-linearity and hysteresis (±5 mmHg) is fulfilled for gas concentrations between 5 % and 10 % CO<sub>2</sub>.</p> <p><b>Drift:</b> The requirement to drift ≤ 10 % of initial reading over the calibration interval is fulfilled.</p>

*Continued on next page*

## Specifications, *Continued*

### Other approvals and compliance (*continued*)

Item	Description
$p\text{CO}_2$ ( <i>continued</i> )	<p>The following drift per hour has been measured at a sensor temperature of 43 °C:            Max. drift per hour at 5 % <math>\text{CO}_2</math>: 2 %            Max. drift per hour at 10 % <math>\text{CO}_2</math>: 2 %            Calibration interval: 4 hours</p> <p><b>Accuracy:</b>            The device indicates the partial pressure of carbon dioxide (cutaneous <math>p\text{CO}_2</math>) to within 5 mmHg over the measurement range.</p> <p><b>Response time (10 % to 90 % response):</b>            The following max. response times have been measured at a sensor temperature of 43 °C:            E5280/E5260: 26 seconds            E5480: 83 seconds</p>
$p\text{O}_2$	<p><b>Non-linearity and hysteresis:</b>            The requirement to non-linearity and hysteresis (<math>\pm 5</math> mmHg) is fulfilled for gas concentrations between 2 % <math>\text{O}_2</math> and 20.9 % <math>\text{O}_2</math>.</p> <p><b>Drift:</b>            The requirement to drift <math>\leq 5</math> % of initial reading over the calibration interval is fulfilled.            The following drift per hour has been measured at a sensor temperature of 43 °C:            Max. drift per hour at 10 % <math>\text{O}_2</math>: 1 %            Max. drift per hour at 20.9 % <math>\text{O}_2</math>: 1 %            Calibration interval: 4 hours</p> <p><b>Accuracy:</b>            The device indicates the partial pressure of oxygen (cutaneous <math>p\text{O}_2</math>) to within 5 mmHg over the range from 0 % <math>\text{O}_2</math> to 20.9 % <math>\text{O}_2</math>. In the range from 20.9 % <math>\text{O}_2</math> to full scale, the accuracy is better than <math>\pm 10</math> %.</p> <p><b>Response time (10 % to 90 % response):</b>            The following max. response times have been measured at a sensor temperature of 43 °C:            E5280: 18 seconds            E5480: 26 seconds</p>

*Continued on next page*

## Specifications, *Continued*

### Other approvals and compliance (*continued*)

Item	Description		
SpO <sub>2</sub>	<b>Accuracy over 70 % to 100 %:</b>		
	Sensor model:	Weight range:	Accuracy:
	DS100A	> 40 kg	±3 % SpO <sub>2</sub>
	Oxiband A/N (adults)	> 40 kg	±3 % SpO <sub>2</sub>
	Oxiband A/N (neonates)	< 3 kg	±4 % SpO <sub>2</sub>
	Oxiband P/I	3-40 kg	±3 % SpO <sub>2</sub>
	<b>Test considerations and oximeter accuracy – Functional testers and patient simulators</b>		
	Some models of commercially available bench top functional testers and patient simulators can be used to verify the proper functionality of Nellcor pulse oximeter sensors, cables and monitors. See the individual testing device's operator's manual for the procedures specific to the model of tester being used.		
	While such devices may be useful for verifying that the pulse oximeter sensor, cable and monitor are functional, they are incapable of providing the data required to properly evaluate the accuracy of a system's SpO <sub>2</sub> measurements.		
	Fully evaluating the accuracy of the SpO <sub>2</sub> measurements requires, at a minimum, accommodating the wavelength characteristics of the sensor and reproducing the complex optical interaction of the sensor and the patient’s tissue. These capabilities are beyond the scope of known bench top testers, including known devices which claim to measure sensor LED wavelength.		
SpO <sub>2</sub> measurement accuracy can only be evaluated in vivo by comparing pulse oximeter readings with values traceable to SaO <sub>2</sub> measurements obtained from simultaneously sampled arterial blood using a laboratory CO-oximeter.			
Many functional testers and patient simulators have been designed to interface with the pulse oximeter's expected calibration curves and may be suitable for use with Nellcor monitors and/or sensors. Not all such devices, however, are adapted for use with the Nellcor <i>OXIMAX</i> digital calibration system. While this will not affect use of the simulator for verifying system functionality, displayed SpO <sub>2</sub> measurement values may differ from the setting of the test device.			
For a properly functioning monitor, this difference will be reproducible over time and from monitor to monitor within the performance specifications of the test device.			
Pulse	<b>Accuracy:</b> ±3 bpm over 20-250 bpm range		

*Continued on next page*

## Specifications, *Continued*

### Other approvals and compliance (*continued*)

Item	Description
Known sources of interference	<p><b>tcpCO<sub>2</sub>/tcpO<sub>2</sub>:</b></p> <p>After an hour's exposure, halothane (4 % evaporated into carrier gas) will interfere with the tcpO<sub>2</sub> readings of the sensor and thereby also influence the combined effects of non-linearity and hysteresis (beyond ±6 mmHg/±0.8 kPa).</p> <p>Furthermore, if the patient is treated with Halothane, this may cause changes in the skin blood flow and in this way interfere with the tcpO<sub>2</sub> readings.</p> <p>Halothane does <i>not</i> affect the tcpCO<sub>2</sub> readings.</p> <p>The following substances have been tested and shown not to influence the combined effects of non-linearity and hysteresis (within ±6 mmHg/±0.8 kPa):</p> <ul style="list-style-type: none"> <li>• Nitrous oxide (in a gas mixture of 60 % N<sub>2</sub>O, 5 % CO<sub>2</sub>, 20.9 % O<sub>2</sub>, balance N<sub>2</sub>). In the following referred to as carrier gas.</li> <li>• Enflurane: 5 % evaporated into carrier gas.</li> <li>• Isoflurane: 5 % evaporated into carrier gas.</li> <li>• Sevoflurane: 5 % evaporated into carrier gas.</li> <li>• Desflurane: 12 % evaporated into carrier gas.</li> </ul> <p><b>SpO<sub>2</sub>:</b></p> <p>Inaccurate measurements can be caused by:</p> <ul style="list-style-type: none"> <li>• prolonged patient movement</li> <li>• venous pulsations</li> <li>• intravascular dyes, such as indocyanine green or methylene blue</li> <li>• defibrillation</li> <li>• incorrect application of the sensor</li> <li>• placement of the sensor on an extremity with a blood pressure cuff, arterial catheter or intravascular line</li> <li>• ambient light</li> </ul> <p>Loss-of-pulse signal can occur for the following reasons:</p> <ul style="list-style-type: none"> <li>• The sensor is applied too tightly</li> <li>• A blood pressure cuff is inflated on the same extremity as the one with the sensor attached</li> <li>• There is arterial occlusion proximal to the sensor</li> </ul>

*Continued on next page*

## Specifications, *Continued*

### Other approvals and compliance (continued)

Item	Description
Biocompatibility	<p><b>Fixation ring:</b> The suitability of the fixation ring adhesives for use on intact human skin is supported by a series of in vitro and in vivo evaluations.</p> <p>The fixation ring adhesives meet the requirements of ISO 10993-1, "Biological Evaluation of Medical Devices. Part 1: Evaluation and Testing" for surface devices contacting intact human skin. The evaluations include cell cytotoxicity, skin irritation and sensitization potential.</p> <p><b>Nellcor SpO<sub>2</sub> sensors (TCM40 only):</b> The sensors have passed the recommended biocompatibility testing and are therefore in compliance with ISO 10993-1.</p>

### Patents

<p>One or more of the following patents and patent applications may apply:</p> <ul style="list-style-type: none"> <li>• US Patent No.: US7,474,908</li> <li>• US Patent Application No.: US2007/0238943</li> <li>• European Patent Application Nos.: EP1753343, EP2007272</li> <li>• Japanese Patent Application No.: JP2007/537799</li> <li>• Chinese Patent Application No.: CN1988849</li> <li>• International Patent Application No.: WO2007/115568</li> </ul> <p>Other patents pending</p>
---

### Materials and accessories

Item	Description
All materials and accessories	All materials and accessories are latex-free

### E5480 sensor

Item	Description
Dimensions	<p>Diameter:</p> <ul style="list-style-type: none"> <li>• Sensor housing: 15 mm 0.6 in</li> <li>• Silver body: 8.1 mm 0.32 in</li> </ul> <p>Height: 7.81 mm 0.3 in</p> <p>Weight: 2.6 g 0.08 oz</p>
Sensor cable	<p>Length: 2.25 m 88.6 in</p>

*Continued on next page*



## Specifications, *Continued*

### E5280 and E5260 sensors

Item	Description
Dimensions	Diameter: <ul style="list-style-type: none"> <li>• Sensor housing: 15 mm 0.6 in</li> <li>• Silver body: 9.1 mm 0.36 in</li> </ul> Height: 11.3 mm 0.44 in Weight: 2.9 g 0.1 oz
Sensor cable	Length: 2.25 m 88.6 in

### OxiMax sensors

Item	Description
Dissipation	52.5 mW
Wavelength	The wavelength range of the light emitted is near 660 nm and 890 nm.

## Accessories

### TCM4/40 systems accessories

Description	Code no.
E5480 tcpCO <sub>2</sub> /tcpO <sub>2</sub> sensor (tinyTeddy)	945-660
E5280 tcpCO <sub>2</sub> /tcpO <sub>2</sub> sensor	945-377
E5260 tcpCO <sub>2</sub> sensor	945-655
Membraning kit for E5480 sensor, containing: <ul style="list-style-type: none"> <li>• tcpCO<sub>2</sub>/tcpO<sub>2</sub> electrolyte solution</li> <li>• 10 membraning tools</li> <li>• Cleaning paper, bag with 5 pieces</li> </ul>	905-805
Membraning kit for E5280 and E5260 sensors, containing: <ul style="list-style-type: none"> <li>• tcpCO<sub>2</sub>/tcpO<sub>2</sub> electrolyte solution</li> <li>• 12 green membrane units with PP membrane</li> <li>• O-ring remover key</li> <li>• Cleaning paper, bag with 20 pieces</li> </ul>	904-892
Fixation kit for E5480 sensor, containing: <ul style="list-style-type: none"> <li>• 4 × 25 disposable fixation rings</li> <li>• 4 × 20 mL contact liquid</li> </ul>	905-836
Fixation kit for E5280 and E5260 sensors, containing: <ul style="list-style-type: none"> <li>• 4 × 25 disposable fixation rings</li> <li>• 4 × 20 mL contact liquid</li> </ul>	904-891
Large fixation ring for E5480 tcpCO <sub>2</sub> /tcpO <sub>2</sub> sensor (tinyTeddy)	905-901

### SpO<sub>2</sub> starter kits

Description	Code no.
SpO <sub>2</sub> module with Nellcor DS100A sensor	902-876
SpO <sub>2</sub> module with Nellcor Oxiband A/N sensor	902-877
SpO <sub>2</sub> module with Nellcor Oxiband P/I sensor	902-878
SpO <sub>2</sub> module without sensor	902-946

### TCM4/40 documentation

Item	Code no.
TCM4/40 operator's manual, English	990-277

*Continued on next page*

## Accessories, *Continued*

### Line cords

Description	Code no.
Line cord 120 V, USA and JPN	615-407
Line cord 230 V, UK	615-312
Line cord 230 V, ITA	615-313
Line cord 230 V, ISR	615-315
Line cord 230 V, other 230 V countries	615-303
Line cord 230 V, AUS and NZA	615-317
Line cord 230 V, ZAF and IND	615-318

### Additional items

Description	Code no.
Extension cable for tcpCO <sub>2</sub> /tcpO <sub>2</sub> sensor, 3 m (118.1 in)	617-853
Extension cable for tcpCO <sub>2</sub> /tcpO <sub>2</sub> sensor, 6 m (236.2 in)	617-864
TCM4xx ETX Serial Adapter	636-649
TCM4xx ETX Analog Adapter	636-650
TCM4xx ETX VueLink Adapter	636-651
tcpCO <sub>2</sub> /tcpO <sub>2</sub> module	902-778
Calibration chamber gasket for E5480 sensor	837-488
Calibration chamber gasket for E5280 and E5260 sensors	837-159
CAL1 standard calibration gas: 7.5 % CO <sub>2</sub> , 20.9 % O <sub>2</sub> , balance N <sub>2</sub> , 180 mL	962-187
CAL1 standard calibration gas: 7.5 % CO <sub>2</sub> , 20.9 % O <sub>2</sub> , balance N <sub>2</sub> , 180 mL (US and Canada)	962-188
CAL2 standard calibration gas: 10 % CO <sub>2</sub> , balance N <sub>2</sub> , 999 mL <b>NOTICE:</b> Not applicable in Germany.	962-096
CAL2 standard calibration gas: 10 % CO <sub>2</sub> , 1 L <b>NOTICE:</b> Only applicable in Germany.	962-154
Valve key (for CAL2 gas)	922-509
Gas adapter (for CAL2 gas)	847-398
Gasket for adapter (847-398)	837-487
12 V 2AH lead-acid battery	431-018



## 9. Functional description

<b>tcpCO<sub>2</sub>/tcpO<sub>2</sub> measurement.....</b>	<b>9-2</b>
Construction of sensors .....	9-3
pCO <sub>2</sub> measuring principle .....	9-5
pO <sub>2</sub> measuring principle.....	9-6
Calibration of sensor .....	9-7
Solutions and calibration gases .....	9-10
Traceability certificates .....	9-11
<b>Pulse oximetry measurement .....</b>	<b>9-15</b>
Measuring principle .....	9-16
Calibration of sensor .....	9-17

## **tcpCO<sub>2</sub>/tcpO<sub>2</sub> measurement**


Construction of sensors .....	9-3
pCO <sub>2</sub> measuring principle .....	9-5
pO <sub>2</sub> measuring principle.....	9-6
Calibration of sensor .....	9-7
Solutions and calibration gases .....	9-10
Traceability certificates .....	9-11

## Construction of sensors

### Sensor components

The  $\text{tcpCO}_2/\text{tcpO}_2$  sensors (E5480 and E5280) combine a heating element, two temperature sensors, a Clark-type oxygen electrode, and a Severinghaus-type carbon dioxide electrode in a single unit.

### E5480 sensor components




Number	Component
1	NTC resistors – temperature sensors
2	Platinum cathode ( $p\text{O}_2$ part)
3	Electrolyte covering the sensor surface
4	$\text{CO}_2/\text{O}_2$ -permeable membranes
5	Reinforced solid-state glass electrode ( $p\text{CO}_2$ part)
6	Ag/AgCl reference electrode
7	Electrolyte reservoir
8	Heating element

*Continued on next page*

## Construction of sensors, *Continued*

### E5280 sensor components



**NOTICE:** The single tcpCO<sub>2</sub> sensor (E5260) is identical to the tcpCO<sub>2</sub>/tcpO<sub>2</sub> sensor (E5280) except that it does not contain an active platinum cathode (the pO<sub>2</sub> part).

Number	Component
1	NTC resistors – temperature sensors
2	Platinum cathode (pO <sub>2</sub> part)
3	Electrolyte covering the sensor surface
4	CO <sub>2</sub> /O <sub>2</sub> -permeable membranes
5	Reinforced solid-state glass electrode (pCO <sub>2</sub> part)
6	Ag/AgCl reference electrode
7	Electrolyte reservoir
8	O-rings to secure the membranes
9	Heating element

**NTC resistors** The temperature of the sensor is measured by the NTC resistors incorporated in the Ag/AgCl reference electrode. Due to the high thermal conductivity of the silver body, the NTC resistors respond quickly to any changes in temperature. The thermostating system will keep the sensor at the preset temperature.


**Local vasodilation** When the sensor is attached to the skin, the generated heat is transferred from the heating element via the silver body to the skin surface. The heat produces local vasodilation and increases the permeability of the skin to oxygen and carbon dioxide, thus making a measurement on the skin surface possible.



## **$p\text{CO}_2$ measuring principle**

**Definition** The measurement of  $p\text{CO}_2$  is defined as the partial pressure (or tension) of carbon dioxide. It is performed by means of a Stow-Severinghaus electrode based on an electrochemical electrode chain consisting of a pH glass electrode (the sensor electrode) and a silver chloride reference electrode.

**Measuring principle** As  $\text{CO}_2$  is released from the skin, it diffuses through the membrane into the electrolyte, where it reacts with water forming carbonic acid, which immediately dissociates into  $\text{HCO}_3^-$  and  $\text{H}^+$  according to the following equation:



The changes in  $\text{H}^+$  in the electrolyte imply changes in pH.

As the pH in the electrolyte changes, the voltage between the glass electrode and the reference electrode changes. The pH change is converted to a  $p\text{CO}_2$  reading on the basis of the linear relationship between pH and  $\log p\text{CO}_2$ , as expressed by the Henderson-Hasselbalch equation:

$$\text{pH} = \text{pK} + \log \frac{[\text{HCO}_3^-]}{a \times p\text{CO}_2}$$

where

$\text{pK}$  = dissociation constant of carbonic acid

$[\text{HCO}_3^-]$  = concentration of  $\text{HCO}_3^-$

$a$  = solubility coefficient of dissolved  $\text{CO}_2$

$p\text{CO}_2$  = partial pressure of  $\text{CO}_2$

As no charged molecules can penetrate the membrane, the change in pH is strictly due to the carbon dioxide diffusion into the electrolyte.

The potential measured across the combined electrode chain is fed into the  $p\text{CO}_2$  channel, where it is digitized. The digitized signal is then passed on to the microcomputer, where it is converted to display  $p\text{CO}_2$  in mmHg or kPa.

The monitor is based on a microcomputer system in which all data obtained by the sensor are collected, processed and compared with the alarm limits preselected in the monitor's setup programs prior to presentation on the screen.

### **Severinghaus temperature correction**

In most clinical settings, transcutaneous  $p\text{CO}_2$  monitoring is performed using the Severinghaus temperature correction factor.

This means that the  $\text{tcpCO}_2$  readings are corrected to 37 °C (normal body temperature), using the following formula:


$$\text{tcpCO}_2(T) = p\text{CO}_2(37\text{ °C}) \times 10^{-0.019(T-37\text{ °C})}$$

where T is the set sensor temperature (°C).


## **$pO_2$ measuring principle**

**Definition**      The measurement of  $pO_2$  is defined as the partial pressure (or tension) of oxygen. It is performed as a direct polarographic measurement based on an electrochemical electrode chain consisting of the platinum cathode (the sensor electrode) and the silver anode (the reference electrode).

**Measuring principle**      The sensor tip is covered with a thin membrane which stabilizes the conditions of  $pO_2$  diffusion to the sensor. Oxygen diffuses through this membrane to the cathode where a reduction of oxygen occurs as a result of the current-generating process:



At the anode the following reaction takes place:



The reduction of oxygen at the platinum cathode generates a current, which is fed into the  $pO_2$  channel, where it is converted into a voltage and digitized. This digitized signal is then passed on to the microcomputer, where it is reconverted to display  $pO_2$  in mmHg or kPa.

## Calibration of sensor

### Temperature-corrected $p\text{CO}_2$ values

In order for the calibration value to be within the normal physiological range for  $\text{tcpCO}_2$ , Radiometer uses a 7.5 %  $\text{CO}_2$  concentration in the calibration gas for the TCM4/40 monitors.

The 7.5 %  $\text{CO}_2$  calibration gas will give the following  $\text{CO}_2$  calibration value (at a barometric pressure of 760 mmHg):

$$p\text{CO}_2(\text{CAL}) = B \times \frac{\% \text{CO}_2}{100} = 760 \times \frac{7.5}{100} = 57 \text{ mmHg}$$

When the Severinghaus temperature correction factor is activated (Severinghaus corr. "ON"), the above  $p\text{CO}_2$  calibration value will result in the following temperature-corrected  $p\text{CO}_2$  values:

Sensor temp °C	37	38	39	40	41	42	43	44	45
Temperature corr. factor	1.00	1.04	1.09	1.14	1.19	1.24	1.30	1.36	1.42
Temperature corr. value, mmHg*	57	55	52	50	48	46	44	42	40

$$* \text{ Temperature corr. } p\text{CO}_2 \text{ value} = \frac{57}{\text{Temp. corr. factor}}$$

$$\text{Temp. corr. factor} = 10^{-0.019(T-37^\circ\text{C})}$$

### Gas calibration values

The monitor calculates the nominal dry gas  $\text{tcpCO}_2$  calibration values relative to the barometric pressure in accordance with the following formula:

$$p\text{CO}_2(\text{CAL}) = B \times \frac{\% \text{CO}_2}{100}, \text{ where}$$

$B$  = the barometric pressure in mmHg or kPa

$\% \text{CO}_2$  = the percentage of  $\text{CO}_2$  in the calibration gas (i.e. 7.5 %)

The monitor calculates the nominal dry gas  $\text{tcpO}_2$  calibration values relative to the barometric pressure in accordance with the following formula:

$$p\text{O}_2(\text{CAL}) = B \times \frac{\% \text{O}_2}{100}, \text{ where}$$

$B$  = the barometric pressure in mmHg or kPa

$\% \text{O}_2$  = the percentage of  $\text{O}_2$  in the calibration gas (i.e. 20.9 %)


---

*Continued on next page*

## Calibration of sensor, *Continued*

### Calibration process

The calibration process includes the stages described in the table below:



Stage	Description
Heating	The sensor is typically heated to the preset temperature within one minute. If it has not reached this temperature within three minutes, the computer will reject the sensor, and the error message "Calibration error. Sensor temperature error during calibration." will be displayed.
Check of sensor signal range	When the sensor reaches the preset temperature, the monitor checks the sensor $p\text{CO}_2$ and $p\text{O}_2$ signals. If these are not within a specified range within one minute, the monitor rejects the sensor, and the error message "Calibration error. Sensor sensitivity error during calibration." will be displayed. <b>NOTICE:</b> The sensor may also be rejected because of an incorrect calibration value.
Check of stability	After heating and acceptance of the sensor signal range, the monitor checks the sensor stability for maximum 10 minutes or until the change in the sensor signal is below 1 % compared with the signals registered one minute earlier. If, following this, the stability criterion is still not fulfilled, the monitor will reject the sensor, and the error message "Calibration error. Sensor drift during calibration." will be displayed. If the stability is accepted, the monitor will then display "Ready".

*Continued on next page*

## Calibration of sensor, *Continued*

### Patient safety

In order to satisfy the requirements for patient safety, the circuitry that is interconnected with the sensor is galvanically insulated from the rest of the system.

The monitor is equipped with an extensive safety system for controlling and monitoring sensor temperature.

- The microprocessor program gives an audiovisual alert if the sensor temperature deviates  $\pm 0.3$  °C from the preset value.
- The sensor power consumption is permanently supervised. If it exceeds maximum effect for more than two minutes, the heat switches off.
- If the sensor temperature reaches 46 °C or falls below 5 °C, this could indicate a short or an open circuit in the thermostating control loop. In this case, the monitor utilizes one of several hardware and software controls to immediately and permanently switch off the sensor heating.
- If there is a temperature difference of more than 0.6 °C between the two thermistors, heating switches off immediately and permanently.
- If a software error occurs, a "watchdog" circuit in the monitor immediately and permanently switches off the sensor heating.

## Solutions and calibration gases

### Solutions

The following solutions are used with the TCM4/40 systems:

Solution	Description
tcpCO <sub>2</sub> /tcpO <sub>2</sub> electrolyte solution, 10 mL	<p>Use: To be applied on the clean sensor surface.</p> <p>Composition: 1,2-propanediol, propanetriol, potassium chloride, potassium hydrogen carbonate and deionized water.</p> <p>Storage: At room temperature or below.</p> <p>Stability: Expiration date and lot no. are printed on a separate label on the bottle.</p>
Contact liquid, 20 mL	<p>Use: To establish contact between the skin and the sensor during in vivo measurements.</p> <p>Composition: 1,2-propanediol and deionized water.</p> <p>Storage: At room temperature or below.</p> <p>Stability: Expiration date and lot no. are printed on a separate label on the bottle.</p>

**NOTICE:** Keep the bottle caps on when bottles are not in use.

### Calibration gases

- CAL1 standard calibration gas (7.5 % CO<sub>2</sub>, 20.9 % O<sub>2</sub> with N<sub>2</sub> as balance)
- CAL2 standard calibration gas (10 % CO<sub>2</sub> with N<sub>2</sub> as balance)



**WARNING – Risk of explosion**


Calibration gas cylinder: Contents under pressure. Do not puncture. Do not use or store near heat or open flame. Exposure to temperatures above 54 °C (for CAL2) and 50 °C (for CAL1) may cause contents to vent or cause bursting. Never discard container into fire or incinerator as it may cause an explosion.



**WARNING – Risk of explosion**

Before discarding an empty CAL2 gas cylinder, remove the safety valve using the valve key (code no. 922-509). If you do not do this, the calibration gas cylinder may burst if exposed to heat.

## Traceability certificates

Certificate of Traceability			
<b>Product name:</b>	CAL1 standard calibration gas		
<b>Type:</b>	7.5 % CO <sub>2</sub> , 20.9 % O <sub>2</sub> , balance N <sub>2</sub> , 180 mL		
<b>Code:</b>	962-187		
<b>Traceability of parameters:</b>			
Parameter	Unit	Traceable to	Expanded Uncertainty
CO <sub>2</sub>	mol %	Primary, gravimetrically prepared standards. Traceable to NIST traceable weights.	±0.07
O <sub>2</sub>	mol %	Primary, gravimetrically prepared standards. Traceable to NIST traceable weights.	±0.03
<p><b>Certification:</b> Each lot of this product has been tested, and the nominal values, specified on the label of this product, have been established with the above traceability.</p> <p style="text-align: right;">   H.B. Kristensen  Head of Chemical Reference Laboratory </p> <p>The traceability of the above parameters is fully described in booklet AS 117: <i>Traceability to the Primary Reference Standards at Radiometer</i>, available from Radiometer.</p>			

Radiometer™ is a trademark of Radiometer A/S.  
 ©Radiometer A/S, DK-2700 Brønshøj, Denmark, 2004. All Rights Reserved.  
 990-346 • 200409A

**RADIOMETER**  
**COPENHAGEN** 

## Certificate of Traceability

**Product name:** CAL1 standard calibration gas (US and Canada)

**Type:** 7.5 % CO<sub>2</sub>, 20.9 % O<sub>2</sub>, balance N<sub>2</sub>, 180 mL

**Code:** 962-188

**Traceability of parameters:**

Parameter	Unit	Traceable to	Expanded Uncertainty
CO <sub>2</sub>	mol %	Primary, gravimetrically prepared standards. Traceable to NIST traceable weights.	±0.07
O <sub>2</sub>	mol %	Primary, gravimetrically prepared standards. Traceable to NIST traceable weights.	±0.03

**Certification:** Each lot of this product has been tested, and the nominal values, specified on the label of this product, have been established with the above traceability.

*Bjarne Kristensen*

H.B. Kristensen  
Head of Chemical Reference Laboratory

The traceability of the above parameters is fully described in booklet AS 117: *Traceability to the Primary Reference Standards at Radiometer*, available from Radiometer.

Radiometer™ is a trademark of Radiometer A/S.  
©Radiometer A/S, DK-2700 Brønshøj, Denmark, 2004. All Rights Reserved.  
990-347 • 200409A

**RADIOMETER**  
**COPENHAGEN** 



**NOTICE:** The below-mentioned CAL2 calibration gas is not applicable to Germany. For information on the CAL2 calibration gas that applies to Germany, see next page.

## Certificate of Traceability

**Product name:** CAL2 standard calibration gas

**Type:** 10 % CO<sub>2</sub>, balance N<sub>2</sub>, 999 mL

**Code:** 962-096

**Traceability of parameters:**

Parameter	Unit	Traceable to	Expanded Uncertainty
CO <sub>2</sub>	mol %	Primary, gravimetrically prepared standards. Traceable to NIST traceable weights.	±0.02
O <sub>2</sub>	mol %	Primary, gravimetrically prepared standards. Traceable to NIST traceable weights.	±0.03

**Certification:** Each lot of this product has been tested, and the nominal values, specified on the label of this product, have been established with the above traceability.

*Bjarne Kristensen*

H.B. Kristensen  
Head of Chemical Reference Laboratory

The traceability of the above parameters is fully described in booklet AS 117: *Traceability to the Primary Reference Standards at Radiometer*, available from Radiometer.

Radiometer™ is a trademark of Radiometer A/S.  
©Radiometer A/S, DK-2700 Brønshøj, Denmark, 2004. All Rights Reserved.  
990-348 • 200409A

**RADIOMETER**  
**COPENHAGEN** 

**NOTICE:** The below-mentioned CAL2 calibration gas applies to Germany only.

## Certificate of Traceability

**Product name:** Cal 2 Gas, 10 % CO<sub>2</sub>

**Type:** Gas mixture, 1 L

**Code:** 962-154

**Traceability of parameters:**

Parameter	Unit	Traceable to	Expanded Uncertainty
CO <sub>2</sub>	mol %	Primary, gravimetrically prepared standards. Traceable to NIST traceable weights.	0.02
O <sub>2</sub>	mol %	Primary, gravimetrically prepared standards. Traceable to NIST traceable weights.	0.02

**Certification:** Each lot of this product has been tested, and the nominal values, specified on the label of this product, have been established with the above traceability.

*Bjarne Kristensen*

H.B. Kristensen  
Head of Chemical Reference Laboratory

The traceability of the above parameters is fully described in booklet AS 117: *Traceability to the Primary Reference Standards at Radiometer*, available from Radiometer.

Radiometer™ is a trademark of Radiometer A/S.  
©Radiometer A/S, DK-2700 Brønshøj, Denmark, 2003. All Rights Reserved.  
989-853 • 0303A

**RADIOMETER**  
**COPENHAGEN** 

## Pulse oximetry measurement

Measuring principle .....	9-16
Calibration of sensor .....	9-17

## Measuring principle

### Introduction

The monitor uses pulse oximetry to measure functional oxygen saturation in the blood. Pulse oximetry works by applying a sensor to a pulsating arteriolar vascular bed, such as a finger or toe. The sensor contains a dual light source and a photodetector.

Bone, tissue, pigmentation and venous vessels normally absorb a constant amount of light over time. The arteriolar bed normally pulsates and absorbs variable amounts of light during the pulsations. The ratio of light absorbed is translated into a measurement of functional oxygen saturation (SpO<sub>2</sub>).

Because a measurement of SpO<sub>2</sub> is dependent upon light from the sensor, excessive ambient light can interfere with this measurement.

Specific information about ambient conditions, sensor application and patient conditions is contained throughout this manual.

### Measuring principles

Pulse oximetry is based on two principles: that oxyhemoglobin and deoxyhemoglobin differ in their absorption of red and infrared light (spectrophotometry), and that the volume of arterial blood in tissue (and hence, light absorption by that blood) changes during the pulse (plethysmography). A pulse oximeter determines SpO<sub>2</sub> by passing red and infrared light into an arteriolar bed and measuring changes in light absorption during the pulsatile cycle. Red and infrared low-voltage light-emitting diodes (LED) in the oximetry sensor serve as light sources; a photodiode serves as the photodetector.

**NOTICE:** Information about the range of wavelength can be useful to clinicians in particular. For information on the range of wavelength, see *OxiMax sensors* in the section *Specifications* in chapter 8.

Because oxyhemoglobin and deoxyhemoglobin differ in light absorption, the amount of red and infrared light absorbed by blood is related to hemoglobin oxygen saturation. To identify the oxygen saturation of arterial hemoglobin, the monitor uses the pulsatile nature of arterial flow. During systole, a new pulse of arterial blood enters the vascular bed, and blood volume and light absorption increase.

During diastole, blood volume and light absorption reach their lowest point. The monitor bases its SpO<sub>2</sub> measurements on the difference between maximum and minimum absorption (measurements at systole and diastole). By doing so, it focuses on light absorption by pulsatile arterial blood, eliminating the effects of non-pulsatile absorbers such as tissue, bone and venous blood.

### Oxygen saturation vs oxyhemoglobin fraction

The monitor measures oxygen saturation, also called functional saturation (oxygenated hemoglobin expressed as a percentage of the hemoglobin that can transport oxygen). It does not detect significant amounts of dysfunctional hemoglobin, such as carboxyhemoglobin or methemoglobin. To compare oxygen saturation measurements with those from an instrument that measures oxyhemoglobin fraction, also erroneously called fractional saturation (oxygenated hemoglobin expressed as a percentage of all measured hemoglobin, including measured dysfunctional hemoglobins), oxyhemoglobin fraction must be converted as follows:

$$\text{oxygen saturation} = \frac{\text{oxyhemoglobin fraction}}{100 - (\% \text{ carboxyhemoglobin} + \% \text{ methemoglobin})} \times 100$$

## Calibration of sensor

**Description** Because light absorption by hemoglobin is wavelength dependent and because the mean wavelength of LEDs varies, an oximeter must know the mean wavelength of the sensor's red LED to accurately measure SpO<sub>2</sub>.

During monitoring, the monitor's software selects coefficients that are appropriate for the wavelength of that individual sensor's red LED; these coefficients are then used to determine SpO<sub>2</sub>.

Additionally, to compensate for differences in tissue thickness, the light intensity of the sensor's LEDs is adjusted automatically.



# Index

## A

Accessories  
 additional items ..... 8-15  
 fixation kit ..... 8-14  
 line cords ..... 8-15  
 membranizing kit ..... 8-14  
 SpO<sub>2</sub> starter kits ..... 8-14  
 TCM4/40 documentation ..... 8-14  
 TCM4/40 systems ..... 8-14

Accuracy  
 pCO<sub>2</sub> ..... 8-9  
 pO<sub>2</sub> ..... 8-9  
 pulse ..... 8-10  
 SpO<sub>2</sub> ..... 8-10

Alarm  
 pCO<sub>2</sub> ..... 3-10  
 pO<sub>2</sub> ..... 3-11  
 Pulse ..... 3-12  
 reset ..... 6-35  
 SatSeconds ..... 6-36  
 silence ..... 6-35  
 SpO<sub>2</sub> ..... 3-12  
 tone ..... 6-35  
 VueLink ..... 6-27

Alarm mode  
 latching ..... 3-20  
 non-latching ..... 3-20

Alarm sound level 3-12

Alarms 6-34

Alert  
 tone ..... 6-36  
 VueLink ..... 6-27

Alerts 6-34, 6-35

All rights reserved 1-2

Ambient light sources 6-5

Ambient temperature 4-2

Analog output 2-4, 6-21

Analog range  
 pCO<sub>2</sub> ..... 3-21  
 pO<sub>2</sub> ..... 3-21

Anemia 6-4

Application of sensor  
 DS100A ..... 6-9  
 E5260 ..... 6-7  
 E5280 ..... 6-7  
 E5480 ..... 6-6  
 Oxiband A/N and P/I ..... 6-9

## B

Battery 2-4, 4-7  
 change ..... 4-7

indicator ..... 2-3  
 level ..... 4-7  
 recharging ..... 4-7

Biocompatibility 8-12

Blood gas comparison 6-44

Blood gas values  
 calibrating transcutaneous  
 measurements against ... 3-13, 6-45  
 comparing transcutaneous  
 measurements with ..... 3-13, 6-44

## C

CAL1 9-10

CAL2 9-10

Calibration  
 automatic ..... 5-2  
 frequency ..... 5-2  
 in vivo ..... 6-45  
 manual ..... 5-3  
 materials ..... 5-2  
 patient safety ..... 9-9  
 process ..... 9-8  
 recommendation ..... 5-2  
 SpO<sub>2</sub> sensor ..... 9-17  
 with calibration gas ..... 5-3

Calibration chamber 2-3

Calibration gas 9-10

Calibration setup 3-14

Calibration status 3-16

Calibration values 9-7

Certificates of traceability 9-11

Changing  
 battery ..... 4-7  
 gas bottle ..... 4-8  
 gasket in calibration chamber... 4-8

Cleaning  
 cable ..... 4-15  
 E5480 sensor head ..... 4-11  
 exterior ..... 4-6  
 sensor head ..... 4-15  
 SpO<sub>2</sub> sensors ..... 4-16  
 soak method ..... 4-16  
 wipe method ..... 4-16  
 touch screen ..... 4-6

Combi module 2-3

Common name 1-2

Compliance 8-8

Connecting to  
 external chart recorder ..... 6-21  
 external equipment ..... 4-3  
 external PC ..... 6-22, 6-30  
 polysomnograph ..... 6-21

PPMS monitor..... 6-24  
 printer ..... 6-38  
 Construction of sensors 9-3  
 Contact liquid 6-7, 9-10  
 Continuous data output 3-20, 6-22  
   MonLink..... 6-29  
   standard ..... 6-22  
   VueLink ..... 6-24  
 CSA 8-7  
 Cursor 6-19  
 Curve report 6-41, 6-43

**D**

Data dump 6-30, 6-33  
 Data export 6-30  
   serial ..... 6-30  
   USB ..... 6-33  
 Data Management System 6-13  
 Date/time settings 3-22  
 Default values 3-23  
 Dialog box messages 7-6  
 Disinfection  
   outer surfaces ..... 4-6  
   sensor ..... 4-15  
 Disinfection solutions 4-15  
 Display brightness 3-21  
 Documentation 8-14  
 Drift  
    $p\text{CO}_2$  ..... 8-8  
    $p\text{O}_2$  ..... 8-9  
 Dysfunctional hemoglobins 6-4

**E**

E5260 sensor components 9-4  
 E5280 sensor components 9-4  
 E5480 sensor  
   membraning ..... 4-12  
 E5480 sensor components 9-3  
 E5480 sensor head  
   cleaning ..... 4-11  
 Electrocautery 6-3  
 Electrolyte solutions 9-10  
 EMC approvals and compliance 8-4  
 Enriched oxygen atmosphere 4-2,  
   6-3  
 Environment of use 1-2  
 Environmental requirements 4-2  
 Error messages 7-3  
    $\text{SpO}_2$  ..... 7-6  
   tc ..... 7-3  
 Event 6-11  
 Excel 6-23, 6-30  
 Exporting data to  
   memory stick ..... 6-33

PC ..... 6-30  
 External chart recorder 6-21

**F**

Factory defaults 3-23  
 Fixation kit 6-6  
 Fixation ring 6-7  
 Flammable anesthetics 4-2  
 Front view 2-3  
 Functional oxygen saturation 9-16

**G**

Gas adapter 8-15  
 Gas bottle  
   maintenance ..... 4-8  
 Gas calibration values 9-7  
 Gas level 3-16, 5-3  
 Gas status 5-3  
 Gasket  
   change of ..... 4-8

**H**

Handle 2-3  
 Henderson-Hasselbalch equation  
   9-5  
 High-frequency electrical signals  
   6-3  
 Hyperbaric chamber 4-2, 6-3  
 HyperTerminal 6-22, 6-30  
 Hyperthermia 6-3

**I**

Importing data files into Microsoft  
   Excel 6-31  
 In vivo calibration 3-20, 6-45  
 In vivo monitoring 6-11  
 Inaccurate measurements 6-4,  
   6-5, 8-11  
 Installation 4-3  
 Intended use 1-2  
 Interference  
   from ambient light ..... 6-5  
   known sources ..... 8-11

**K**

Keyboard 6-15

**L**

Latching alarm 6-34  
 Light absorption 9-16  
 Line fuse 2-4  
 Linking session 6-12  
 Local vasodilation 9-4  
 Loss-of-pulse signal 6-4, 6-5, 8-  
   11



**M**

- Main screen
  - when not activated for
    - 30 seconds ..... 3-2, 3-5
- Markings
  - normal view ..... 6-17
  - trend curve ..... 6-20
  - trend table ..... 6-18
- Measurement
  - $p\text{CO}_2$  ..... 9-5
  - $p\text{O}_2$  ..... 9-6
  - $\text{SpO}_2$  ..... 9-16
- Measuring principle
  - $p\text{CO}_2$  ..... 9-5
  - $p\text{O}_2$  ..... 9-6
  - $\text{SpO}_2$  ..... 9-16
- Measuring sites 6-2
- Membraning 4-10
  - E5260 sensor ..... 4-13
  - E5280 sensor ..... 4-13
  - E5480 sensor ..... 4-12
- Menu screen
  - configuration ..... 2-7
  - touch keys ..... 2-7
- Menu structure 3-2
- Meta. corr. factor 3-20
- Module release 2-4
- Modules
  - $\text{SpO}_2$  ..... 2-3
  - $\text{tcpCO}_2/\text{tcpO}_2$  ..... 2-3
- MonLink 6-29

**N**

- Name 1-2
- Nominal dry gas calibration values
  - 9-7
- Non-latching alarm 6-34
- Non-linearity and hysteresis
  - $p\text{CO}_2$  ..... 8-8
  - $p\text{O}_2$  ..... 8-9
- Normal view 3-5, 6-16
- NTC resistors 9-4

**O**

- On-screen keyboard 6-15
- Operating requirements 4-2
- Operator profile 1-2
- O-ring 4-13
- O-ring remover 4-13
- Oxygen saturation
  - fractional ..... 9-16
  - functional ..... 9-16
- Oxygenated hemoglobin 9-16
- Oxyhemoglobin fraction 9-16

**P**

- Parameter setup 3-9
- Parameters
  - $p\text{CO}_2$  ..... 3-10
  - $p\text{O}_2$  ..... 3-11
  - Pulse ..... 3-12
  - $\text{SpO}_2$  ..... 3-12
- Parts and functions
  - rear ..... 2-4
  - top and front ..... 2-3
- Patents 8-12
- Patient DMS 6-12, 6-13
- Patient ID 6-12, 6-13
- Patient monitoring 6-11
- Patient movement 6-5
- Patient safety 8-7
- Patient stabilization time 6-11
- PCL-compatible PC printer 6-38
- $p\text{CO}_2$ 
  - measurement ..... 9-5
  - measuring principle ..... 9-5
  - sensitivity ..... 5-4
- $p\text{CO}_2$  3-10
- $p\text{CO}_2$  alarm
  - high ..... 3-10
  - low ..... 3-10
- Performance considerations
  - $\text{SpO}_2$  sensors ..... 6-5
  - TCM40 monitor ..... 6-5
- Philips Patient Monitoring System
  - 6-24
- Photodetector 9-16
- Photodiode 9-16
- Physiological alarms 6-34
- Plethysmography 9-16
- $p\text{O}_2$  3-11
  - measurement ..... 9-6
  - measuring principle ..... 9-6
  - zero current ..... 5-4
- $p\text{O}_2$  alarm
  - high ..... 3-11
  - low ..... 3-11
- $p\text{O}_2$  parameter only 3-20
- Polysomnograph 6-21
- Power socket 2-4
- Power switch 2-4
- PPMS monitor 6-24
  - configuration ..... 6-27
- Print reports 6-43
- Printer setup 3-17
- Printer type 3-18
- Printing 6-38
  - USB ports ..... 6-38
- Printouts 3-18

- Proprietary name 1-2
- Pulse alarm
  - high..... 3-12
  - low..... 3-12
- Pulse oximetry
  - measuring principles ..... 9-16
- Pulse rates 6-4
- Pure text messages 7-8
- R**
- Rear view 2-4
- Recharging battery 4-7
- Relative humidity 4-2
- Remembraning 4-10
- Removal of wraps 6-10
- Report
  - curve..... 6-41
  - table..... 6-39
- Requirements
  - ambient temperature..... 4-2
  - environmental..... 4-2
  - operating..... 4-2
  - relative humidity ..... 4-2
  - ventilation..... 4-2
- Response time
  - $p\text{CO}_2$ ..... 8-9
  - $p\text{O}_2$  ..... 8-9
- Results
  - Normal view..... 6-16
  - Trend curve view..... 6-19
  - Trend table view..... 6-18
- S**
- SatSeconds 3-12, 6-36
  - safety net ..... 6-37
- Saturation 6-4, 9-16
- Saturation monitoring 6-8
- Screen
  - general elements..... 2-5
  - menu..... 2-5
  - types ..... 2-5
  - view ..... 2-5
  - when not activated for
    - 30 seconds..... 3-2, 3-5
- Screen saver 3-6, 3-21
- Sensor
  - E5480 ..... 6-6
- Sensor
  - DS100A..... 6-9
  - E5260 ..... 6-7
  - E5280 ..... 6-7
  - Oxiband A/N ..... 6-9
  - Oxiband P/I ..... 6-9
  - storage ..... 4-15
  - temperature..... 6-2
- Sensor application procedure
  - DS100A..... 6-9
  - E5260 ..... 6-7
  - E5280 ..... 6-7
  - E5480 ..... 6-6
  - Oxiband A/N and P/I..... 6-9
- Sensor components
  - E5260 ..... 9-4
  - E5280 ..... 9-4
- Sensor components
  - E5480 ..... 9-3
- Sensor socket
  - $\text{SpO}_2$  ..... 2-4
  - $\text{tcpCO}_2/\text{tcpO}_2$ ..... 2-4
- Sensor temperature 3-10, 3-11
- Serial port (RS232) 2-4
- Session 6-12, 6-13
- Session number 6-13
- Setting up
  - analog output..... 6-21
- Setup programs 3-3
  - Printer..... 3-17
- Severinghaus corr. 3-20
- Severinghaus temperature correction
  - 9-5
- Shutting down 4-5
- Site time 3-10, 3-11, 6-11
- Site time heat 3-10, 3-11
- SmartCal 3-15, 5-2
- SmartHeat 3-10
- Solutions
  - composition..... 9-10
  - stability..... 9-10
  - storage ..... 9-10
  - use ..... 9-10
- Specifications
  - E5260 sensor ..... 8-13
  - E5280 sensor ..... 8-13
  - E5480 sensor ..... 8-12
  - monitor ..... 8-2
  - OxiMax sensors..... 8-13
- Spectrophotometry 9-16
- $\text{SpO}_2$ 
  - measurement ..... 9-16
  - measuring principles ..... 9-16
  - module..... 9-16
  - sensor performance
    - considerations..... 6-5
- $\text{SpO}_2$  alarm
  - high ..... 3-12
  - low..... 3-12
- $\text{SpO}_2$  module 2-3
- $\text{SpO}_2$  monitoring
  - required materials..... 6-8

- Stabilization time of patient 6-11
- Standard calibration gas 5-4
- Standard data output 6-22
- Storage of sensor 4-15
- Symbols
  - used in the manual..... 1-6
  - used on the monitor..... 1-4
- T**
- Table report 6-39, 6-40
- TCM4 Series monitoring system 1-2
- TCM4/40 monitoring systems 1-2
- TCM40
  - performance considerations ..... 6-4
- TCM40 SpO<sub>2</sub> module 9-16
- tcpO<sub>2</sub>/tcpCO<sub>2</sub> electrolyte solution 9-10
- Technical alarms 6-34
- Technical settings 3-20
- Technical setup 3-19
- Temperature, sensor 6-2
- Time settings 3-22
- Touch key glossary 2-9
- Touch screen 2-3
- Traceability certificates 9-11
- Trend curve view 3-8, 6-19
- Trend table view 3-7, 6-18
- Troubleshooting
  - cause of error..... 7-2
  - dialog box messages..... 7-6
  - error messages..... 7-3
  - errors ..... 7-2
  - pure text messages ..... 7-8
  - recommended operator actions 7-2
  - symptoms ..... 7-2
- Tutorials 2-6, 2-8
- U**
- USB 6-33
- USB ports 2-4, 6-38
- V**
- Valve key 8-15
- Ventilation requirements 4-2
- View screen
  - configuration ..... 2-5
  - touch keys ..... 2-6
- View setup 3-4
- Views
  - Normal ..... 3-5
  - Trend curve ..... 3-8
  - Trend table ..... 3-7
- VueLink 6-24
- W**
- Waste of Electrical and Electronic Equipment 1-4
- WEEE 1-4
- What is what 2-1
- Z**
- Zero current 5-4



## Date of issue

---

**Radiometer representative:**

**Manufacturer:**

**RADIOMETER** 

If you have any questions or need assistance, please contact your local Radiometer representative.



Radiometer Medical ApS  
Åkandevvej 21  
2700 Brønshøj  
Denmark  
[www.radiometer.com](http://www.radiometer.com)

**TCM4/40 monitoring systems  
operator's manual  
from software version 3.00**

Publication 201101

Edition C

Code number 994-941



0459

